
medi+physics™

Announces An Ideal Radioisotope For The Study of Pulmonary Ventilation

- A half-life of 13 seconds and decay by Isomeric Transition means low radiation exposure to patients and staff.
- The monoenergetic gamma emission of 191 keV is well suited for the gamma camera.
- No special radioactive gas collection or disposal system required.
- Completely portable system allows studies in ICU, CCU, and Post-Surgical departments with portable camera.



- Studies can be conducted on comatose, uncooperative, or mechanically vented patients.
- Distribution of radioactive gas is mainly to the lungs.
- Elaborate delivery system is not required.
- The only radioisotope that can be administered ON and OFF as needed.
- Easy to license when compared to Xenon Xe 133 gas.

MPI Krypton Kr 81m Gas Generator
Krypton Kr 81m

The Pulmonary Profile

THE CONCEPT

The pulmonary profile is a series of matched perfusion and ventilation studies done consecutively on a patient using the MPI Krypton Kr 81m Gas Generator and Technetium Tc 99m Albumin Aggregated. Following administration of the two products you are able to switch the energy window on the gamma camera and scan the patient in the same position for each of the isotopes before you move the patient to the next view. Thus, a complete series of matching views may be accumulated for any number of patient positions.

THE PURPOSE

To increase the diagnostic sensitivity and specificity of lung imaging procedures by providing an easy means of obtaining matched perfusion-ventilation images in one patient visit.

THE RESULT

A new patient study which combines ventilation and perfusion imaging procedures into one study called the *Pulmonary Profile Study*.

For information regarding the MPI Krypton Kr 81m Gas Generator Krypton Kr 81m please call Medi-Physics at (415) 658-2184, Outside California (800) 227-0492 or Inside California at (800) 772-2477.

medi+physics™

For complete prescribing information please consult package insert, a summary of which follows:

MPI KRYPTON Kr 81m GAS GENERATOR KRYPTON Kr 81m

DESCRIPTION: The Krypton Kr 81m Gas Generator consists of Rubidium Rb 81 fixed to a solid support from which the Krypton Kr 81m is eluted by passage of humidified oxygen or air through the generator. Other rubidium radio-isotopes which do not decay to radioactive Krypton Kr 81m in their decay are present in the generator (Rubidium Rb 82m, for example, is present at a concentration of 30-40%).

INDICATIONS AND USAGE: The Krypton Kr 81m Gas Generator is indicated for use in the study of pulmonary ventilation.

CONTRAINDICATIONS: None known.

WARNINGS: None known.

PRECAUTIONS:

General

The Krypton Kr 81m Gas Generator as well as other radioactive drugs, must be handled with care to minimize radiation exposure to clinical personnel. Also care should be taken to minimize radiation exposure to the patient consistent with proper patient management.

Carcinogenesis, Mutagenesis, Impairment of Fertility

No long-term animal studies have been performed to evaluate carcinogenic potential, mutagenicity potential, or whether Krypton Kr 81m gas affects fertility in males or females.

Pregnancy-Category C

Animal reproduction studies have not been conducted with Krypton Kr 81m gas. It is also not known whether Krypton Kr 81m gas can cause fetal harm when administered to a pregnant woman, or can affect reproduction capacity. Krypton Kr 81m gas should be given to a pregnant woman only if clearly needed.

Nursing Mothers

It is not known whether this drug is excreted in human milk. Because many drugs are excreted in human milk, caution should be exercised when Krypton Kr 81m gas is administered to a nursing woman.

Ideally, examinations using radiopharmaceuticals, especially those elective in nature, in women of childbearing capability, should be performed during the first few (approximately ten) days following the onset of menses.

Pediatric Use

Safety and effectiveness in children have not been established.


Radiopharmaceuticals should be used only by physicians who are qualified by training and experience in the safe use and handling of radionuclides and whose experience and training have been approved by the appropriate government agency authorized to license the use of radionuclides.

ADVERSE REACTIONS: None known.

DOSAGE AND ADMINISTRATION: The recommended dose range for Krypton Kr 81m is 1-10 millicuries and should be administered by continuous inhalation for a sufficient time to provide desired diagnostic information. The multiplication product of the radioactivity and the time of continuous inhalation of Krypton Kr 81m generally should not exceed 100 millicurie-minutes.

HOW SUPPLIED: The Krypton 81m Gas Generator is supplied in the form of Rubidium Rb 81, bound to a solid support, with an activity of 2-10 millicuries at calibration time. The generator is enclosed in a lead shielded filter assembly surrounded by a capped plastic canister to which a handle is affixed. The generator should be stored at room temperature. The generator expires 12 hours after date and time of calibration.

Technical Support



Technical questions answered by experts. Just give us a toll-free call.

Mallinckrodt gives you an open line to our comprehensive technical support system... a system designed to meet all your technical needs.

If you have a question about test procedures, applications or need technical assistance on any of our radiopharmaceuticals, we're here to help. If you encounter any problems, no matter how minor, we want to make things right... right away. Need information on the U.S. Nuclear Regulatory Commission or other State or Federal agencies? Call us for that too.

Whatever technical support you need, Mallinckrodt is just a toll-free call away.

800-325-8181, Your Technical Support Number
(In Missouri, 314-895-2405 collect)

For ordering information, call: 800-325-3688
(In Missouri, 314-344-3880 collect)

THE MALLINCKRODT COMMITMENT

to Nuclear Medicine

MALLINCKRODT, INC., Post Office Box 5840, St. Louis, MO 63134



Yesterday



1970

Today



1980

Tomorrow



**Brought to you
in part
by NEN**

For the past decade, nuclear medicine has enjoyed a continuing stream of new radiopharmaceuticals, new isotopes, new diagnostic procedures — and new patients. Many of these new diagnostic procedures resulted directly or indirectly from the investments in product research and development, testing, production, and promotion by a single company: New England Nuclear.

We supported investigators with grants to develop their ideas into agents suitable for animal and human testing... we invested in the production facilities to manufacture sufficient quantities of radiopharmaceuticals and isotopes to perform the studies necessary to bring new products to you.

And then, we underwrote an effort unique in nuclear medicine — we began spending hundreds of thousands of dollars each year to inform primary-care physicians and specialists why they should send their patients to nuclear departments for these new studies.

Such investments in new product development and physician education are common among traditional pharmaceutical companies producing proprietary products that can be patented. However, all NEN's investments were made on products for which no exclusivity of patent protection was available. Some of NEN's investments were not successful. A few were, however — and they profoundly changed nuclear medicine.

Of course, NEN could have waited for other companies to develop new

procedures and products... to carry the risk and investment of pioneering trial and error. We could have waited until someone else created a demand for new isotopes, and then capitalized on their efforts.

Instead, we built *four* of our own cyclotrons, and are currently building a multimillion-dollar linear accelerator — further evidence of NEN's unique commitment to research and development innovation in isotope and radiopharmaceutical production.

If NEN had not been so committed to advancing nuclear diagnostics, perhaps bone scans might still be done with strontium... and techniques such as tumor, abscess, and myocardial perfusion imaging might still be subjects for academic — not clinical — consideration.

NEN has maintained a high level of customer acceptance of its isotopes and radiopharmaceuticals, thanks to physicians and technologists who understand that when they trust their business to NEN they are sharing our investment in future nuclear diagnostics... in the profession's future ability to diagnose diseases for which medicine has no agents today... and in the effort to communicate the benefits of nuclear diagnostics to the medical community.

We're committed to you. We're **NEN** New England Nuclear®

THE STABLE SOLUTION TO YOUR BONE IMAGING NEEDS

"Tc-99m-MDP with ascorbate had a higher affinity for hydroxylapatite than did MDP without ascorbate."

"Based on these studies and previous shelf life studies, the authors conclude that Tc-99m-MDP with ascorbate is different and is the agent of choice."

From: DIFFERENCES BETWEEN Tc-99m-MDP STABILIZED AND UNSTABILIZED. K.T. Study, K.A. Reed, and D.L. Laven. University of New Mexico, Albuquerque, NM; reprinted with author's permission from J.N.M.T., Vol. 8, No. 2, p. 127. Study based on biodistribution in Swiss-Webster mice.



- One Year Shelf Life
- No Refrigeration Required
- Full 6 Hour Use After Preparation
- Contains Ascorbic Acid as an Antioxidant

For ordering, customer service, and technical information, call toll-free 800-431-1146 (in NYS call 800-942-1986).

CintiChem[®]
MDP^{KIT}

Technetium Tc 99m Medronate Kit

BRIEF SUMMARY OF PRESCRIBING INFORMATION

indications and usage

Technetium Tc 99m Medronate may be used as a bone imaging agent to delineate areas of altered osteogenesis.

contraindications

None known.

warnings

This class of compound is known to complex cations such as calcium. Particular caution should be used with patients who have, or who may be predisposed to, hypocalcemia (i.e., alkalosis).

This radiopharmaceutical drug product should not be administered to children, to pregnant women, or to nursing mothers, unless the expected benefit to be gained outweighs the potential risk.

Ideally, examinations using radiopharmaceuticals, especially those elective in nature, of a woman of childbearing capability should be performed during the first few (approximately 10) days following the onset of menses.

precautions

general

Technetium Tc 99m Medronate as well as other radioactive drugs, must be handled with care and appropriate safety measures should be used to minimize radiation exposure to clinical personnel. Also, care should be taken to minimize radiation exposure to the patients consistent with proper patient management.

To minimize radiation dose to the bladder, the patient should be encouraged to void when the examination is completed and as often thereafter as possible for the next 4-6 hours.

This preparation contains no bacteriostatic preservative.

Technetium Tc 99m Medronate should be formulated within six (6) hours prior to clinical use.

pregnancy category C

Adequate reproductive studies have not been performed in animals to determine whether this drug affects fer-

tility in males or females, has teratogenic potential, or has other adverse effects on the fetus. Technetium Tc 99m Medronate should be used in pregnant women only when clearly needed.

nursing mothers

It is not known whether this drug is excreted in human milk. As a general rule nursing should not be undertaken while a patient is on the drug since many drugs are excreted in human milk.

pediatric use

Safety and effectiveness in children have not been established.

adverse reactions

No adverse reactions specifically attributable to the use of Technetium Tc 99m Medronate have been reported.

how supplied

Union Carbide's Technetium Tc 99m Medronate Kit is supplied as a sterile, pyrogen-free kit containing 5 vials.

Each 10 ml vial contains 10 mg medronic acid, 0.17 mg (minimum) stannous chloride (maximum stannous and stannic chloride 0.29 mg), and 2 mg ascorbic acid. The pH has been adjusted to 4-8 with either HCl or NaOH prior to lyophilization. Following lyophilization, the vials are sealed under a nitrogen atmosphere.

Product #17500502 Multidose vial shield with cap and retainer ring available separately.




FROM ATOM TO IMAGE

Manufactured For:

Union Carbide Corporation • Medical Products Division •
Nuclear Products • P.O. Box 324 • Tuxedo, New York 10987

CintiChem is a registered trademark of Union Carbide Corporation.



**700
TIMES
MORE
SENSITIVE**

¹²⁵I METHOTREXATE RADIOIMMUNOASSAY KIT

Our ¹²⁵I Methotrexate Radioimmunoassay Kit provides a rapid, simple method with an unexcelled level of sensitivity and specificity.

Here is a comparison chart that speaks for itself.

Select the proven DBI ¹²⁵I MTX-RIA kit to monitor the circulating methotrexate levels in serum, plasma, cerebral spinal fluid or urine.

Also available:

- 125I Doxorubicin-RIA Kit
- 125I Digoxin-Stat-RIA Kit
- 125I Folate Kit
- 125I T₄-One Step-RIA Kit
- 125I T₃-Uptake Kit


Call or write for our low priced introductory kit.

	DBI RADIOIMMUNOASSAY	IMMUNOENZYME ASSAY
STAT INCUBATION:	15 minutes at 37°C	1 minute
SENSITIVITY:	0.0004 μM (700 times more sensitive)	0.3 μM
EXOGENOUS INTERFERENCE:	None	Lypemic Icterus Hemolysis
STANDARDS SUPPLIED:	7	6
PRICE:	*57½ cents per tube	\$1.86 per tube

*In units of 200

**Diagnostic
Biochemistry
Inc.**

10457-H Roselle Street, San Diego, CA 92121
Tel. (714) 452-0950



MPI

Thallous Chloride Tl 201

Thallous Chloride Tl 201

Medi-Physics announces the addition of MPI Thallous Chloride Tl 201 to its product line. Medi-Physics can offer our customers a complete line of cyclotron and reactor produced radiopharmaceuticals.

medi+physics™

5801 Christie Ave., Emeryville, Calif. 94608, (415) 658-2184

multicrystal count rates with single-crystal image quality

Another unique feature of the Apex Line




In the practice of modern Nuclear Medicine, physicians have learned that a camera's major diagnostic advantage is often negated by a parallel disadvantage.

High count rate is one such case.

Until Apex, high count rates were achievable only with multicrystal cameras—at the expense of image quality. Only Elscint's Apex Line provides count rates as high as 500,000 CPS and resolutions as fine as 1.8mm bars (Apex 215M).

Elscint Inc.
138-160 Johnson Avenue,
Hackensack, N.J. 07602, U.S.A.
Call Toll Free: 800-631-1694




High Count Rates—The Clinical Need

As Nuclear Medicine techniques become more sophisticated, they require higher count rates. Cardiac first-pass studies, for example, can only be effectively accomplished with count rates exceeding the limitations of most present day gamma cameras. Apex systems, however, do perform these studies—with remarkable image clarity.

Some Impressive Apex Qualities

The remarkable *Count Rate* performance of the Apex Line is supported by a high *Dynamic Frame Rate* of 64 FPS for 64×64 pixels, and a *Multigated Frame Rate* of 64 frames per heart cycle for 64² matrix.



Pick a n

Any of these four ADAC nuclear medicine systems will deliver useful clinical information faster and more accurately than ever before possible.

And one has exactly the capacity your clinic or hospital needs today.

That's the one to buy.

Later, as you expand, your ADAC expands with you.

ADAC
System

I.



It processes and acquires.

It has a Diagnostic Acquisition/Processing Console, a Computer Section, and a Winchester disc drive.

It can easily be expanded to a System II.

ADAC
System

II.



It processes one study while acquiring or processing another.

It's a System I — plus a Remote Acquisition/Processing Terminal, a second Computer Section, and a second Winchester.

It can easily be expanded to a System III.

umber.

ADAC
System

III.



It processes and acquires in two places at once.
It looks just like a System II.
But there's more capability inside the Computer Section.
It can easily be expanded to a System IV.

ADAC
System

IV.



It has a three-location option.
With two Consoles and an expanded Computer Section, you
can process and acquire in two places at once.
Add an optional Remote Terminal and you can process and
acquire studies at three locations.

For more information on ADAC Systems I, II, III, & IV, write
or call collect. ADAC Laboratories, 255 San Geronimo Way,
Sunnyvale, California 94086. (408) 736-1101.

ADAC
Nuclear Medicine Computers

FIRST AGAIN.



Being first can be habit forming. We were first with the 2-mm resolution mobile gamma camera. We were first with 12-inch (30 mm) FOV resolution good enough for cardiology. We were first with instant on-line ejection fraction. We were first with energy calibration. Question: Who's first with small field of view resolution in large field of view cameras? Picker. Again.


It took 61 photomultiplier tubes and an ingenious optimization of crystal design to improve resolution 25%. INTRODUCING THE DYNACAMERA 4/15/61 CLINICAL SYSTEM —

- Intrinsic resolution and linearity so good, myocardial tomography is now a practical reality.
- Improved clarity in small-lesion detection.
- Better reading of low-contrast liver, lung, and renal studies.
- Over 30,000 clinical imaging hours in nuclear medicine departments.

DynaCamera 4/15/61: 61-tube resolution in a compact, newly designed chassis. It's not only first, but it's good enough to remain first. There's more you should know. Write: Picker Corporation, 12 Clintonville Rd., Northford, CT 06472, or Picker International, 595 Miner Rd., Highland Hts., OH 44143.

**THE
IMAGE
OF
VALUE.**

PICKER®



A² Image Processing Systems

Total Capability and a Built-In Future

Every A² System offers total capability. Whether you select a single terminal system, mobile or stationary, or a multiple terminal system, you can perform all existing acquisition processing and display protocols with the same ease of operation and superior image quality.

Of course you can expand A² Systems to operate at one, two or three locations. Or, place an addi-


tional terminal at a remote viewing and processing station. All of this from one central console without duplication of central processors or storage media.

A² Systems offer more than a "built-in future". They have your future built-in today. Please write or call for our new A² brochure—or ask for a demonstration of "Total Capability".

m d s


Medical Data Systems

division of Medtronic, Inc.
2311 Green Road
Ann Arbor, MI 48105
313 769 9353
Telex 235794



MDS products, hardware and software, are tools for diagnosis and research which do not come in contact with, and cannot cause direct injury to the patient. Refer to the operation manual and instructions accompanying the gamma camera and injectable imaging agent for further information on their use. To ensure proper clinical results, an MDA product must be used under the direction of, and using procedures verified by a qualified physician.

Elution Flexibility



4 vial sizes for one-step elution & dilution.

With a Mallinckrodt generator, you can elute 5, 10, 20 or 30 ml volumes with no time-consuming dilutions. And you don't have to load saline vials, so there's less exposure time. When generator activity is high, simply use a larger vial. By use of the proper size vial, the internal saline reservoir automatically provides the proper volume for the concentration you need. Saves preparation time, minimizes handling and reduces contamination risk. And eluting flexibility is part of every package, because all the vial sizes you need are available at no extra charge. Just let us know the sizes you need.

Find out more about the generator that gives you a choice of elution volumes with no extra steps. Call your Mallinckrodt representative or this toll-free number:

800-325-3688 (In Missouri, 314-344-3880 collect)

For technical assistance it's 800-325-8181
(In Missouri, 314-895-2405 collect)

THE MALLINCKRODT COMMITMENT



to Nuclear Medicine
MALLINCKRODT, INC., Post Office Box 5840, St. Louis, MO 63134

**When a
SQUIBB
product
comes
into your lab,
SERVICE
comes
with it.**

Ever since the first Squibb radiopharmaceutical was delivered over two decades ago, we have known that good products accompanied by good service result in customer satisfaction. We are proud of our record of professional service to nuclear medicine technologists and physicians and of our efforts to help fulfill their unique needs. Our current programs continue this tradition.


Squibb Health Physics Service

Specially trained Technical Associates aid laboratory personnel with a pre-inspection of the laboratory to help ensure compliance with radiation safety regulations.



Customtec®

Exclusive computer service that custom-tailors generator size and delivery schedule to a laboratory's daily technetium Tc 99m requirements.



Technical Associates

They have had extensive training and experience in nuclear medicine, radiopharmaceuticals, RIA and instrumentation. They can aid in license procurement and renewal, laboratory design and expansion, technologist training, trouble-shooting and instrument calibration.



Squibb National Nuclear Medicine Management Seminars

Two-and-a-half-day program in management training and development. Designed for chief technologists, supervisors and persons who hold positions of equivalent or greater responsibility.



Squibb National Nuclear Medicine Seminars

Educational course for technologists. Two days on *in vivo* procedures and two days on *in vitro* procedures. Approved for continuing education credits.



Technologist Education Plan

A lead recovery program in which Squibb puts funds into a customer's education account each time spent Minitec® (Technetium Tc 99m) Generators are returned.



Your Squibb Representative will be happy to give you more information about these programs and other Squibb services. Technical assistance can be obtained from Technical Customer Service, 609-921-4100.



Back to Basics!

The Assayer 1 by Radx

The never ending struggle for product popularity often leads a manufacturer to add gadgets. It's called "one-upmanship." We sometimes lose sight of what YOU, the user, wants.

By customer demand, Radx has gone "Back to Basics" and developed the Assayer 1, a simple dosecalibrator, a reliable dose-calibrator, an economical dosecalibrator.

The return to basics does not require a

return to the 1960's technology. The Assayer 1 is microprocessor controlled, totally solid state, with a method of isotope selection way ahead of its time (an optical scanner) which is so precise, reproducible, and reliable that it will soon be copied.


It is not a gadget, it calibrates doses accurately, with precision and unprecedented reliability. It's the Assayer 1—\$2950.

Call today for the last dosecalibrator you'll ever own.

RADX

P.O. Box 19164 • Houston, Texas 77024 • (713) 468-9628

Sterile Technique



**Quick-change elution needle
reduces risk of contamination
for greater patient safety.**

A fresh needle every time you elute helps insure sterility. Mallinckrodt generators have a unique, quick-change system that makes it easy. Each generator comes with 6 sterile-packed, pyrogen-free needles. After eluting, just remove the old needle and insert the new, leaving the protective cap in place. That way, your generator is always ready to supply sterile eluate. And if a needle is damaged, you can change it in seconds... no waiting for a whole new needle assembly to be delivered and installed, no disruption of your schedule.

For more information, call your Mallinckrodt representative or this toll-free number:

800-325-3688 (In Missouri, 314-344-3880 collect)

For technical assistance it's 800-325-8181
(In Missouri, 314-895-2405 collect)

THE MALLINCKRODT COMMITMENT

to Nuclear Medicine

MALLINCKRODT, INC., Post Office Box 5840, St. Louis, MO 63134



Thallium imaging in acute myocardial infarction

Lewis C. Becker, MD

Associate Professor of Medicine
Director, Nuclear Cardiology
The Johns Hopkins Medical Institutions
Baltimore, Maryland



One of the most significant findings to come from our clinical research over the past several years has been the observation that thallium-201 imaging, performed early after onset of symptoms, can reliably distinguish high-risk and low-risk groups of hemodynamically stable patients with acute myocardial infarction. The value of such a prognostic indicator in the management of acute MI is evident. Patients determined to be at low risk could be ambulated earlier and perhaps discharged sooner than in current practice; in the future, such patients might be placed early in a progressive-care-type unit rather than be maintained in the more expensive coronary care unit.

Patients at higher risk might be found to require more intensive monitoring for even longer periods than today. And following discharge, these patients could justifiably be subjected to much closer and long-lasting followup. Most important, reliable identification of patients at high risk would permit earlier initiation of aggressive treatment directed at limiting the extent of infarction.

Predicting mortality

Our recently reported study¹ covered 42 consecutive patients determined by conventional means (history, ECG, serum creatine kinase) to have suffered an acute MI. These were Killip class I or II patients—the largest group of MI patients, and

those normally considered to be at relatively low risk. All 42 patients were admitted within 12 hours of onset of chest pain, and underwent thallium imaging within 15 hours of onset.


The thallium images—in the anterior, 40° LAO, and 60° LAO views—were interpreted both subjectively and by a computer-assisted quantitative technique.² For each interpretive approach, scores for all views were summed to give a total “defect score”—the lower the score, the smaller the area of thallium defect, with a total defect score of 7 corresponding to reduction in thallium uptake involving approximately 40% of the left ventricle in at least two views. The total defect scores were then correlated with the patients’ subsequent clinical course and with other clinical indices believed to have prognostic value—previous history of MI, anterior location of MI, alveolar infiltrates on admission, peak CK greater than 1,000 IU/liter, age, and sex.

Of the 42 patients, 35 survived the initial hospitalization. These survivors were followed for 6 to 20 months after discharge.

What were our results? Nonsurvivors had significantly larger thallium defects than survivors. The mean score for nonsurvivors was 14.3 vs 2.3 for survivors. In the 13 patients with a score greater than 7—ie, 40% or more involvement—the in-hospital mortality was 46%; at 6 months it was 62%; and at last followup (mean 9 months) it was 92%. In the group of patients with a total defect score less than 7, the in-hospital mortality was 3%; at 6 months and at last followup, it was, respectively, 7% and 7%.

These data conclusively showed that the thallium study performed within hours of admission could identify apparently stable MI patients at high-risk for mortality. In addition, when we compared the predictiveness of the thallium score with the other clinical indices—history, MI location, enzymes, etc—singly *and* in combination, the thallium study was significantly better.

We were, of course, very excited by our results. But, because this was a retrospective study, we felt it important to validate the findings prospectively. Over a 6-month period, we studied more than 90 consecutive patients admitted to the CCU with documented or strongly suspected MI. We applied the same scoring system and same dividing line (score 7)—and confirmed our ability to use thallium imaging to distinguish between high-risk and low-risk groups. The mortality rates of the two groups were almost identical to those established in the earlier retrospective study.



Admission studies in a 66-year-old male with known previous acute MI. The patient was hemodynamically stable. Chest X-ray showed slightly increased heart size, no pulmonary congestion. ECG showed ST elevation in II, III, and aV_F, indicative of acute inferior MI, as well as changes suggestive of old anterior wall damage. Thallium imaging disclosed markedly diminished uptake involving the inferior wall, apex, distal anterior wall, and septum. The total defect score was 40.7. The patient's condition gradually deteriorated, despite aggressive vasopressor administration. He died on the sixth postadmission day.

Irreversible damage and reversible ischemia

We believe the thallium study accurately predicts prognosis in MI patients because the size of the defect reflects the total hypoperfused mass of the left ventricular myocardium—both infarcted and ischemic areas. We know from observations of other investigators that the thallium defect tends to diminish with time after an acute MI. Thus the image recorded immediately after admission will show a larger defect than those recorded on serial followup over subsequent days. Our own pathologic studies have demonstrated that large thallium defects seen on post-MI images may be associated with small areas of infarction on postmortem examination.

Together, these findings strongly support the concept that areas of reduced or absent thallium on the initial post-MI images represent both ischemia and infarction, and that the "filling-in" seen on followup imaging represents resolution of ischemia due to resolution of coronary artery spasm or enlargement of coronary collaterals. Thus, the post-MI study identifies myocardium irreversibly damaged by the acute event, myocardium damaged

by previous infarction, and surrounding areas of severe ischemia that are at risk for necrosis either immediately or at some future time.

Clinical implications

In our patients, the highest percentage of in-hospital deaths was due to sudden pump failure—possibly due to the large total volume of compromised myocardium. Postdischarge deaths were generally related to a new ischemic event. In both of these high-risk groups, the thallium study might have helped in patient management decisions. For those patients who died while in the hospital, more aggressive support might have been indicated; those whose deaths occurred posthospitalization might have been identified as candidates for coronary artery bypass.

References

1. Silverman KJ, Becker LC, Bulkley BH, *et al*: Value of early thallium-201 scintigraphy for predicting mortality in patients with acute myocardial infarction. *Circulation* 61, 996–1003, 1980.
2. Burrow RD, Pond M, Schafer AW, Becker L: "Circumferential profiles:" A new method for computer analysis of thallium-201 myocardial perfusion images. *J Nucl Med* 20, 771–777, 1979.

Please see following page for brief summary of prescribing information.

Thallous Chloride TI 201

INDICATIONS AND USAGE: Thallous Chloride TI 201 may be useful in myocardial perfusion imaging for the diagnosis and localization of myocardial infarction.

It may also be useful in conjunction with exercise stress testing as an adjunct in the diagnosis of ischemic heart disease (atherosclerotic coronary artery disease).

CONTRAINDICATIONS: None known.

WARNINGS: In studying patients in whom myocardial infarction or ischemia is known or suspected, care should be taken to assure continuous clinical monitoring and treatment in accordance with safe, accepted procedure. Exercise stress testing should be performed only under the supervision of a qualified physician and in a laboratory equipped with appropriate resuscitation and support apparatus.

Ideally, examinations using radiopharmaceutical drug products — especially those elective in nature — of women of childbearing capability should be performed during the first ten days following the onset of menses.

PRECAUTIONS: Data are not available concerning the effect of marked alterations in blood glucose, insulin, or pH (such as is found in diabetes mellitus) on the quality of thallium TI 201 scans. Attention is directed to the fact that thallium is a potassium analog, and since the transport of potassium is affected by these factors, the possibility exists that the thallium may likewise be affected.

Thallous Chloride TI 201, as all radioactive materials, must be handled with care and used with appropriate safety measures to minimize external radiation exposure to clinical personnel. Care should also be taken to minimize radiation exposure to patients in a manner consistent with proper patient management.

No long-term animal studies have been performed to evaluate carcinogenic potential.

Adequate reproduction studies have not been performed in animals to determine whether this drug affects fertility in males or females, has teratogenic potential, or has other adverse effects on the fetus. Thallous Chloride TI 201 should be used in pregnant women only when clearly needed.

It is not known whether this drug is excreted in human milk. As a general rule nursing should not be undertaken when a patient is administered radioactive material.

Safety and effectiveness in children have not been established.

The expiration date for Thallous Chloride TI 201 is six days postcalibration.

ADVERSE REACTIONS: Adverse reactions related to use of this agent have not been reported to date.

DOSAGE AND ADMINISTRATION: The recommended adult (70kg) dose of Thallous Chloride TI 201 is 1-1.5mCi. Thallous Chloride TI 201 is intended for intravenous administration only.

For patients undergoing resting thallium studies, imaging is optimally begun within 10-20 minutes after injection. Several investigators have reported improved myocardial-to-background ratios when patients are injected in the fasting state, in an upright posture, or after briefly ambulating.

Best results with thallium imaging performed in conjunction with exercise stress testing appear to be obtained if the thallium is administered when the patient reaches maximum stress and when the stress is continued for 30 seconds to one minute after injection. Imaging should begin within ten minutes post-injection since target-to-background ratio is optimum by that time. Several investigators have reported significant decreases in the target-to-background ratios of lesions attributable to transient ischemia by two hours after the completion of stress testing.

The patient dose should be measured by a suitable radioactivity calibration system immediately prior to administration.

Radiopharmaceuticals should be used by persons with specific training in the safe use and handling of radionuclides produced by nuclear reactor or particle accelerator and whose experience and training have been approved by the appropriate government agencies authorized to license the use of radionuclides.

HOW SUPPLIED: Thallous Chloride TI 201 for intravenous administration is supplied as a sterile, non-pyrogenic solution containing at calibration time, 1mCi/ml of Thallous TI 201, 9mg/ml sodium chloride, and 9mg/ml of benzyl alcohol. The pH is adjusted to between 4.5-6.5 with hydrochloric acid and/or sodium hydroxide solution. Vials are available in the following quantities of radioactivity: 1.5, 3.0, 4.5, 6.0 and 9.0 millicuries of Thallous TI 201.

The contents of the vial are radioactive. Adequate shielding and handling precautions must be maintained.

Catalog Number NRP-427

May 1980



**New England Nuclear
Medical Diagnostics Division**

601 Treble Cove Rd., North Billerica, MA 01862

Call Toll-Free: 800-225-1572 / Telex: 94-0996
(In Mass. and International: 617-482-9595)

Canada: NEN Canada, 2453 46th Avenue, Lachine, Que. H8T 3C9
Tel 514-636-4971

Europe: NEN Chemicals GmbH, D-6072 Dreieich, W. Germany, Postfach 401240
Tel (06103) 85034 Order Entry (06103) 81011

24A

NEW

AUDIOVISUALS IN NUCLEAR CARDIOLOGY

- SI-18 **Basic Concepts in Cardiac Anatomy and Physiology**
by Glen W. Hamilton, M.D.
- SI-19 **The Measurement of Ejection Fraction**
by William Ashburn, M.D.
- SI-20 **Intracardiac Shunts and Cardiac Output**
by William Ashburn, M.D.
- SI-21 **Perfusion Studies of the Ischemic Heart**
by Glen W. Hamilton, M.D.
- SI-22 **Detection of Acute Myocardial Infarction**
by B. Leonard Holman, M.D.
- SI-23 **Instrumentation for Nuclear Cardiology**
by Trevor D. Craddock, Ph.D.

Each audiovisual kit comes complete with expert narration and carefully selected supporting visual materials. Consisting of 35 mm color slides and standard audio cassette, each kit forms a complete self-teaching package. Suitable for individual or group instruction, these units offer active learner participation to reinforce the most important concepts. Each kit has been prepared by an authority in the field, making expert instruction available to you in your home, office or hospital.

SNM Audiovisuals cost \$55.00 each for members of the Society of Nuclear Medicine, \$75.00 each for nonmembers. **There is a 10% discount if all six nuclear cardiology units are ordered at once.** A complete list of SNM Audiovisuals is available on request.

MAIL TO: Audiovisual Department, Society of Nuclear Medicine, 475 Park Avenue South, New York, NY 10016.

Please send the following Audiovisual units. (Check units desired.)

_____ SI-18 _____ SI-20 _____ SI-22
_____ SI-19 _____ SI-21 _____ SI-23

\$55.00 each for members; \$75.00 each for nonmembers.

Total _____ Audiovisual units @ _____ each.

- Please send the complete list of SNM Audiovisuals.
- I plan to use the Audiovisual units on a machine that **automatically** advances the slides. Send one side only audio tapes.

TOTAL \$ _____
DEDUCT 10% IF ORDERING ALL SIX UNITS \$ _____
FOREIGN ORDERS ADD \$7.50 \$ _____
TOTAL ENCLOSED \$ _____


SEND TO:

NAME _____
ADDRESS _____

ZIP _____

Check or purchase order must accompany all orders.
Make checks payable to the Society of Nuclear Medicine.
Inc. **U.S. funds only, please.** **JNM 10/80**

UNIQUE



XenoGard™ Xenon-133 Room Air & Trap Monitor*

The ONLY instrument that integrates ^{133}Xe concentrations (MPC • HRS) for a full week, as required by NRC and Agreement States.¹

Chart recorder option permits continuous hard-copy recording of xenon levels.

- Audible and visual alarms alert you BEFORE a hazardous xenon concentration is reached.
- Monitors and displays ^{133}Xe concentration exiting from gas trap to indicate when filter needs replacement:

⁽¹⁾The Code of Federal Regulations† clearly limits the permissible ^{133}Xe exposure to 1 MPC for 40 hours per week for 13 weeks. The data is continuously updated and displayed by the "XenoGard."

†10 CFR, Part 20, Sec. 20.103 and Appendix B, Table 1.

80% OF USERS FOUND XENON IN ROOM AIR

80% of a sample of XenoGard™ owners reported finding xenon gas releases of which they were previously unaware. Discovery of varying xenon concentrations during "routine" ventilation studies had been virtually impossible to detect prior to using the "XenoGard" Monitor.

VICTOREEN NUCLEAR ASSOCIATES



100 Voice Road
Carle Place, N.Y. 11514
(516) 741-6360

Demonstrations available.

Send for full details. Ask for Bulletin 2660-B

Quality, Comfort, and Protection...

accept nothing less for your eyes

Nuclear Pacific's Wrap-Around Shielding Eyeglasses provide the best eye protection available. The process of grinding, polishing, tempering, and anti-reflective coating the higher density glass is more than just a science. It is a skill. The optics and prescription you receive in these lenses are guaranteed for accuracy and quality. Adjustable rotors are installed in the frames for a personalized and comfortable fit.

For impact resistance standards, the FDA requires that every lens survive the steel ball drop test after tempering. Because of the high lead crystalline content of the glass, three out of four lenses are lost. For this reason, many other manufacturers will not attempt to meet industrial impact standards. Every lens produced for your purchase

has exceeded all federal eyeglass specifications.

Nuclear Pacific has set the standards in radiation shielding glass for 30 years. The optics grinder has 25 years of experience in lens process and systems design. The result is warranted quality from Nuclear Pacific, Inc.

Standard style and clip-ons: \$155. Wrap-around \$197.




"See us at
RSNA in Dallas,
November 16-21"
"Booth #2313"

Nuclear Pacific, Inc.

6701 Sixth Avenue South, Seattle, WA 98108 (206) 763-2170 Telex Number 32-8891

Diagnostic Exclusivity



**TechneScan® PYP® Technetium Tc99m
Pyrophosphate Kit...the only
pyrophosphate kit indicated for
gated cardiac blood pool imaging.**

Mallinckrodt's TechneScan PYP is the only pyrophosphate kit available that gives you the additional diagnostic capability of an advanced method for the dynamic assessment of cardiac function.

Other indications in nuclear cardiology include use as an adjunct in the diagnosis of acute myocardial infarction.

For more information about TechneScan PYP — and all the other organ-imaging kits available from Mallinckrodt — just call your Mallinckrodt representative or this toll-free number:

800-325-3688 (In Missouri, 314-344-3880 collect)

For technical assistance it's 800-325-8181
(In Missouri, 314-895-2405 collect)

See brief summary on following page.

THE MALLINCKRODT COMMITMENT

to Nuclear Medicine

MALLINCKRODT, INC., Post Office Box 5840, St. Louis, MO 63134



TechneScan[®] PYP[®]

Technetium Tc-99m Pyrophosphate Kit

BRIEF SUMMARY

CLINICAL PHARMACOLOGY

When injected intravenously **TechneScan PYP Tc 99m** has a specific affinity for areas of altered osteogenesis. It is also concentrated in the injured myocardium, primarily in areas of irreversibly damaged myocardial cells.

One to two hours after intravenous injection of **TechneScan PYP Tc 99m**, an estimated 40 to 50 percent of the injected dose has been taken up by the skeleton and approximately 0.01 to 0.02 percent per gram of acutely infarcted myocardium. Within a period of one hour, 10 to 11 percent remains in the vascular system, declining to approximately 2 to 3 percent twenty-four hours post injection. The average urinary excretion was observed to be about 40 percent of the administered dose after 24 hours.

TechneScan PYP also has an affinity for red blood cells. When administered 30 minutes prior to the intravenous administration of sodium pertechnetate Tc-99m approximately 76 percent of the injected activity remains in the blood pool providing excellent images of the cardiac chambers.

INDICATIONS AND USAGE

TechneScan PYP Tc 99m is a skeletal imaging agent used to demonstrate areas of altered osteogenesis, and a cardiac imaging agent used as an adjunct in the diagnosis of acute myocardial infarction.

As an adjunct in the diagnosis of confirmed myocardial infarction (ECG and serum enzymes positive), the incidence of false negative images has been found to be 6 percent. False negative images can also occur if made too early in the evolutionary phase of the infarct or too late in the resolution phase. In a limited study involving 22 patients in whom the ECG was positive and serum enzymes questionable or negative, but in whom the final diagnosis of acute myocardial infarction was made, the incidence of false negative images was 23 percent. The incidence of false positive images has been found to be 7 to 9 percent. False positive images have also been reported following coronary by-pass graft surgery, in unstable angina pectoris, old myocardial infarcts and in cardiac contusions.

TechneScan PYP is a blood pool imaging agent which may be used for gated cardiac blood pool imaging. When administered intravenously 30 minutes prior to the intravenous administration of sodium pertechnetate Tc-99m approximately 76 percent of the injected activity remains in the blood pool.

CONTRAINDICATIONS

None.

WARNINGS

This radiopharmaceutical should not be administered to patients who are pregnant or lactating unless the information to be gained outweighs the potential hazards.

Ideally, examinations using radiopharmaceuticals, especially those elective in nature, of a woman of childbearing capability should be performed during the first few (approximately 10) days following the onset of menses.

Warning: Preliminary reports indicate impairment of brain scans using pertechnetate Tc-99m which have been preceded by bone scan. The impairment may result in false positives or false negatives. It is recommended, where feasible, that brain scans precede bone imaging procedures.

Radiopharmaceuticals should be used only by physicians who are qualified by specific training in the safe use and handling of radionuclides produced by nuclear reactor or particle accelerator and whose experience and training have been approved by the appropriate government agency authorized to license the use of radionuclides.

The **TechneScan PYP Kit** must be maintained at refrigerator temperature until use.

The contents of the **TechneScan PYP** reaction vial are intended for use in the preparation of Technetium Tc 99m Stannous Pyrophosphate. **TechneScan PYP** may also be reconstituted with sterile, pyrogen-free normal saline containing no preservatives and injected intravenously prior to the administration of sodium pertechnetate Tc-99m.

Sodium pertechnetate Tc-99m solutions containing an oxidizing agent are not suitable for use with the **TechneScan PYP Kit**.

The contents of the kit are not radioactive. However, after the sodium pertechnetate Tc-99m is added, adequate shielding of the final preparation must be maintained.

TechneScan PYP Tc 99m should not be used more than six hours after preparation.

PRECAUTIONS

As in the use of any other radioactive material, care should be taken to insure minimum radiation exposure to the patient, consistent with proper patient management, and to insure minimum radiation exposure to occupational workers.

Bone Imaging

Both prior to and following **TechneScan PYP Tc 99m** administration, patients should be encouraged to drink fluids. Patients should void as often as possible after the **TechneScan PYP Tc 99m** injection to minimize background interference from accumulation in the bladder and unnecessary exposure to radiation.

Cardiac Imaging

Patient's cardiac condition should be stable before beginning the cardiac imaging procedure.

If not contraindicated by the cardiac status, patients should be encouraged to ingest fluids and to void frequently in order to reduce unnecessary radiation exposure.

Interference from chest wall lesions such as breast tumors and healing rib fractures can be minimized by employing the three recommended projections.

Blood Pool Imaging

TechneScan PYP should be injected by direct venipuncture. Heparinized catheter systems should be avoided.

ADVERSE REACTIONS

None.

HOW SUPPLIED

Catalog Number—094 **TechneScan PYP Technetium Tc-99m Pyrophosphate Kit.**


Kit Contains:

5—Stannous Pyrophosphate Reaction Vials (lyophilized) for the preparation of Technetium Tc-99m Stannous Pyrophosphate.

Reaction Vial Contains:


12.0 mg sodium pyrophosphate and 3.4 mg stannous chloride (lyophilized). Hydrochloric acid is added for pH adjustment prior to lyophilization.

5—Radioassay Information String Tags.



Mallinckrodt, Inc.
P.O. Box 5840, St. Louis, Missouri 63134

Stressing your needs... and your patients!



"Connect-able" Unit

Hi-Lo Cardiac Stress Imaging System

A sensitive design for a sensitive test. The 8407 Hi-Lo stress system's "drop-leaf" ergometer is essential for radionuclide imaging with large field-of-view cameras. Drop the leaf for that needed margin of safety between patient and camera head—even with standard or mobile cameras.

Just one example of EDC's response to clinical problems—design solutions.

With your patient in a supine position, you may image continuously during myocardial stress perfusion scintigraphy or radionuclide exercise ventriculography. Add patient participation and eliminate incessant directions via a color-coded meter that allows self-monitoring of pedal speed.

A more relaxed, better-imaged patient. It adds up to better medicine.

The 8407 control console displays actual workload, pedal speed and elapsed time. Optional automatic heart-rate-control-of-workload enables the clinician to select the specific, individual heart rate at which cardiac stress imaging will be conducted.

The versatile table can be used for other nuclear imaging studies, radiology and ultrasound. Image from above or below table. Vary patient position—supine or sitting.

"Connect-able" Unit

To stay current and within your budget, attach the 8414 modular stress imaging unit to your present table for use in nuclear imaging or cardiac catheterization. The "Connect-able" includes multi-

position ergometer and electronic control console with all options. Specially-designed casters provide easy mobility as well as rigid locking power.

EDC stress imaging systems ... for the progressive nuclear medicine department ... for the cardiologist who demands the diagnostic sensitivity of stress imaging.



**Engineering
Dynamics Corporation**
120 Stedman Street
Lowell, Massachusetts 01851
Tel. (617) 458-1456

COLLEGE OF AMERICAN PATHOLOGISTS

presents



ENDORSED BY
AMERICAN COLLEGE
OF RADIOLOGY

NUCLEAR MEDICINE PROGRAMS

in 1981



ENDORSED BY
SOCIETY OF
NUCLEAR MEDICINE

TRANSMISSION IMAGING SIMULATOR SERIES TIS PRICE \$370.00

Endorsed by: AMERICAN COLLEGE OF RADIOLOGY
and the SOCIETY OF NUCLEAR MEDICINE

SHIPPING DATES: TIS•A/May 13 TIS•B/October 20

Each of two shipments will contain an imaging
device permitting:

- Interlaboratory comparison of technique and instrumentation.
- Initiation and maintenance of a continuing quality control program of imaging procedures.

The purpose of the program is to interrelate the demonstration of abnormal scintiphoto findings with various techniques and different imaging devices. On receipt of the subscriber's completed questionnaire, the CAP will send a preliminary report of the specifications of the imaging device. This information will permit the subscriber to immediately use his device with knowledge of what defects are present for daily quality control of the instruments. The imaging device does not contain radioactivity; therefore, no specific license or safety equipment is needed. A summary report will be sent later detailing participant results.

BASIC LIGAND ASSAY SERIES K SHIPPING DATES: PRICE \$208.00

K•A/April 21 K•C/October 6
K•B/July 21 K•D/January 19 ('82)

Each of four shipments contains three samples (10 ml) of different concentrations to be tested for the following constituents. This will provide 168 challenges per year.

The CAP cannot guarantee that all shipments in a series will be available for orders received after December 31, 1980. There is an administrative charge of \$20 for orders receiving partial shipments.

To order the CAP Nuclear Medicine Survey program or to receive a brochure describing all the CAP Survey programs, please write or call:

COLLEGE OF AMERICAN PATHOLOGISTS
7400 NORTH SKOKIE BLVD./SKOKIE, ILLINOIS 60077
PHONE: (312) 677-3500

BASIC LIGAND REFERENCE MATERIAL SERIES KRM PRICE \$132.00

Participants subscribing to the Basic Ligand Assay Survey may purchase additional series of the material used in this program.

BASIC LIGAND REFERENCE MATERIAL (KRM) MUST BE ORDERED IN ADVANCE AND IS AVAILABLE ONLY TO THOSE LABORATORIES ENROLLED IN THE BASIC LIGAND ASSAY PROGRAM.

KRM (four mailings for Basic Ligand Assay Survey participants, a total of twelve 10 ml vials) will be mailed out three weeks after the Basic Ligand Assay Survey kits.

WHY IS KRM UNIQUE?

- No other quality control material available in this form has had more extensive analysis. The concentration of each measured analyte in KRM has been established after critical analysis of the data from the many participating laboratories in the Survey program.
- Virtually every method in use on any system has contributed concentration values to the data base. This allows specific comparison with your own methods to a degree of reliability not found in any other quality control material.
- The data base provided you is unique, and considering the limited numbers of standard materials available for ligand assay, this product provides you the most up to date approach toward insuring a high quality analytical effort.


ADVANCED LIGAND ASSAY SERIES Y PRICE \$276.00

SHIPPING DATES:

Y•A/May 12 Y•C/November 9
Y•B/August 11 Y•D/February 9 ('82)

Each of four shipments contains two specimens (two 10 ml vials/specimen) to be tested for the listed constituents, a total of 112 tests.

Aldosterone	Digitoxin	hFSH	Progesterone
Alpha feto protein	Estradiol	hGH	Prostatic acid
CEA	Gastrin	hLH	Phosphatase*
	hCG	Prolactin	Henin
			Testosterone




Improve your image with a different slant.

With the EDC Rotatable 30° Slant Hole Collimator, . . . achieve optimal clinical evaluation of left ventricular performance by imaging in the MLAO view.

During radionuclide ventriculography of the cardiac blood pool in the modified left anterior oblique (MLAO) view, EDC's slant hole collimator permits placement of the detector flat against the patient's chest because of the holes' 30° caudal angulation. The improved outcome, demonstrated by the above ungated cardiac images of a patient with aortic stenosis and left ventricular hypertrophy, includes:

- Complete resolution and separation of left ventricle and atrium.
- Viewing of ventricular septum normal to its longitudinal axis, with no foreshortening.
- Effective separation between the cardiac apex and base, as well as between distributions of the left anterior descending and left circumflex arteries.

A study by Dr. J. Anthony Parker,* Harvard Medical School, found that radionuclide ventriculography with the EDC 30° slant hole collimator provides "an accurate measure of



ejection fraction at equilibrium and a qualitative assessment of regional changes in ventricular volume."

Refined resolution of the cardiac apex is obtained in the RAO view.

Other applications of the slant hole collimator include imaging of the spleen and posterior cranial fossa.

The EDC slant hole can be mounted on any commercial Anger scintillation camera. Rotatability of the slant hole inserts facilitates correct positioning. Computed tomography is made possible via indexing of the collimator with detents at up to 24 angles.

Other collimators available from EDC: Seven Pinhole, Bifocal Diverging, Div/Con, Parallel Hole.

EDC collimators. . . . To improve your image and your patients'. Write or call EDC for further information:

edc

Engineering Dynamics Corporation

120 Stedman Street
Lowell, Massachusetts 01851
Tel. (617) 458-1456

*PARKER, J.A. et al: Radionuclide left ventriculography with the slant hole collimator. J Nucl Med 18:848-851, 1977.

AccuSync

The finest R-wave Triggering device available for computerized gated cardiac studies.

FEATURES

- Exclusive Double Discrimination provides precise definition of R-wave.
- ECG Strip Chart Recorder
- Four digit LED Display
- Trigger Pulse LED
- Unlimited Heart Rate Capability
- ONE YEAR WARRANTY

BENEFITS

- Computer is gated only on the R-wave. High amplitude T-waves are ignored.
- Provides permanent record of patient ECG. Insures proper lead placement.
- Indicates R-R Interval or Heart Rate during stress studies.
- Monitors presence of output signals to the computer.
- Both Heart Rate display and R-trigger pulses have unlimited tracking capability during stress studies.
- ONE YEAR WARRANTY

MODEL

FEATURES

AccuSync-I

R-Trigger pulse output, ECG output, Heart Rate/R-R int. display, Strip Chart Recorder and Isolation Amplifier for patient safety.



AccuSync-II

All above features incorporated into a Module designed to fit into certain Mobile cameras.



AccuSync-III

All AccuSync-I features with the exception of the Strip Chart Recorder.




AccuSync-IV

All AccuSync-III features with the exception of the Heart Rate R-R int. display.



Advanced Medical Research Corp./P.O. Box 3094 PPS/301 Brewster Road
Milford, CT 06460/Telephone : (203) 877-1610

Personal Service



**Knowledgeable response
to your needs... right now!**

Mallinckrodt has sales representatives who are strategically located throughout the country; there's one near you. You have a representative who's responsive to you... to answer questions, to iron out problems, to place orders.

Our representatives are easy to reach through their own individual 800 numbers and they'll respond to your needs quickly — even when they're away from their phones. They check their answering services at least twice a day... you'll usually hear from your representative within two hours, if not sooner.

For knowledgeable service and a quick response to your needs, your Mallinckrodt representative is as close as your phone. Call your representative or this toll-free number:

800-325-3688 (In Missouri, 314-344-3880 collect)
For technical assistance it's 800-325-8181
(In Missouri, 314-895-2405 collect)


THE MALLINCKRODT COMMITMENT

to Nuclear Medicine

MALLINCKRODT, INC., Post Office Box 5840, St. Louis, MO 63134



PRECISION ACCURACY PERFORMANCE



For the
quantitative
determination of
Prostatic Acid
Phosphatase levels
in serum.

Introducing GammaDab® [¹²⁵I] PAP RIA Kit

- Affinity purified PAP in human serum-based standards
- Superior maximum binding
- Greatest sensitivity in clinically significant range
- Ready-to-use reagents—no reconstitution
- Convenience—choice of same-day or overnight procedures, 60 or 125-tube kits

Innovating for Life



CLINICAL ASSAYS
DIVISION OF TRAVENOL LABORATORIES, INC.


620 Memorial Drive, Cambridge, MA 02139
(617) 492-2526 TLX: 921461 CLASS CAM
Toll free outside Mass: (800) 225-1241

For other worldwide locations, please contact your local Clinical Assays/Travenol representative or the International Sales Department, Clinical Assays, Cambridge, MA 02139 U.S.A.

Complete directions are provided with each product. These directions should be read and understood before use. Particular attention should be paid to all warnings and precautions. Additional performance data are available. Should you have any questions, contact your Clinical Assays/Travenol representative.

3 Big Favorites

OF OVER 500 NUCLEAR PRODUCTS
available in our free catalog




**a cardiac stress system
that does more and costs less.**

Designed for exercise imaging
SPECIFY MODEL 056-180 IN YOUR INFORMATION REQUEST



radiochromatogram scanner

Computerized analysis of
radiopharmaceutical purity
SPECIFY "ATOMASTER" (149-200) IN YOUR INQUIRY



pulmanex xenon system

A single unit
with an integrated gas trap
SPECIFY "PULMONEX" (130-500) IN YOUR INQUIRY

FOR COMPLETE INFORMATION WRITE OR CALL—

Atomic Products Corporation

ATOMLAB DIVISION • ESTABLISHED 1949
P.O. BOX 657 CENTER MORICHES, NEW YORK 11934 USA
(516) 878-1074
TWX #510-228-0449

“In-depth experience keeps Harshaw on the frontiers of scintillation technology. Our warranty is your ticket to join us there.”

“If you’re developing applications dependent on state of the art scintillation detection, chances are you’ll reach your objectives faster by talking with Harshaw.

“The proprietary knowledge we’ve accumulated during four decades of research and production is a primary resource that’s shared with every Harshaw detector customer.

“As participants in scientific programs ranging from space shuttle development to the Los Alamos Scientific Laboratory laser fusion program, Harshaw physicists and engineers have specialized expertise in nuclear detection, dosimetry and electronics.

“And these unique human resources — plus the patented techniques and processes, the innovative facilities, and the precise synthesis and purification controls behind them — mean Harshaw has the flexibility to provide the best detector for your needs.

“We offer a full selection of standard NaI detectors built for premium resolution. The sizes, shapes and configurations in highest demand are stocked for one week delivery from inventory.

We also stock designs utilizing a dozen other scintillators, including high-Z, non-hygroscopic and fast counting crystals. And we have the experience and ability to custom-design and build exotic crystal scintillators to agreed-upon specifications, for delivery on or before the promised date.

“Behind each Harshaw detector are warranties that are second to none. We don’t just tell you our detectors will meet the high performance requirements of nuclear medicine and endure the shock and vibration tests encountered in non-laboratory

applications — we guarantee it in writing.

“First, your detector is guaranteed to meet or exceed specified resolution requirements. Second, it carries a two-year warranty against malfunction due to faulty construction or failure of the hermetic seal. All photo-multiplier tubes coming with detectors we build carry the manufacturer’s warranty extending from the date we ship to you.

“You’ll probably never have to exercise those warranties. But a Harshaw field representative will be in touch after delivery to make certain your detector is performing up to expectations, and to remind you that we’ll be there when and if you need service or consultation help.

“Pioneering experimenters in radiation detection come to us for their scintillators. Shouldn’t you?”



Elmer C. Stewart
Vice President

HARSHAW

THE HARSHAW CHEMICAL COMPANY
Crystal & Electronic Products
6801 Cochran Road
Solon, Ohio 44139
(216) 248-7400



Medi-Ray announces . . .

SURVEY METER

CALIBRATION and REPAIR SERVICE

The Medi-Ray Survey Meter Calibration and Repair Service is designed to provide reliable, competent calibration and repair for the areas of Nuclear Medicine, Radiology, Research and Industry. Our service incorporates the latest techniques and facilities, as well as a staff of highly qualified personnel functioning in the latest and most modern of environments. The result is the highest quality service at a reasonable cost to the customer.

Types of Meters:


- Ionization Chamber
- Geiger - Mueller
- Scintillation

Features:

- New York State Licensed Laboratory
- Three calibration points on each range
- Accuracy \pm 10% of indicated reading
- Low cost - \$50.00 meter calibration
\$50.00 repair service (excluding GM tube replacement)
- Rapid turnaround

For information, write or call collect:

Medi-Ray, Inc. / 150 Marbledale Rd. / Tuckahoe, N.Y. 10707
(914) 961-8484




Medi-Ray, Inc.

Medi-Ray, Inc. / 150 Marbledale Rd. / Tuckahoe, N.Y. 10707
Please send information on calibration service.

Name	_____	Title	_____
Hospital	_____	Dept.	_____
Address	_____	City	_____
State	_____	Zip	_____
Phone	_____		_____

ALL FILM BADGE SERVICES ARE NOT CREATED EQUAL



On the surface, most radiation monitoring services seem pretty much alike. They read the exposed film accurately and reliably. They satisfy your legal record-keeping requirements. And most of the services are reasonably priced.

Look deeper! The Siemens Nuclibadge II Radiation Monitoring Service really takes off where everyone else quits.

Take the speed of reporting, for example. Comprehensive Siemens reports are in the mail only a few days after receipt of the badges for processing.

Our color coding system lets you know at a glance that each person is wearing the right badge.

New badges or replacement film packets are sent to you, on time.

Even if the exposed film has been damaged in the mail, chances are that our

special program will succeed in getting an accurate reading.

Most important of all, our customer service is without equal.

- You get immediate attention to your problems via our HOTLINE—toll free WATS phone service—that puts you in direct contact with helpful experts. Call for any reason at no cost.

- Additional badges are shipped within 24 hours of receipt of your order. Again, there's no extra charge.

- If we can't answer your question immediately, you will receive a call back within 24 hours with an answer.

- You are automatically notified by telephone of any unusually high readings without charge.

- Your reports can be custom designed to fit your procedures.

- Accounts may be subdivided by department. You take advantage of quantity pricing and receive a single monthly bill with charges itemized by department. Yet, for your convenience badges may be sent directly to the individual departments.

LET US SHOW YOU THAT WE MEAN SERVICE

SIEMENS GAMMASONICS, INC.
Health Physics Services
2000 Nuclear Drive
Des Plaines, Illinois 60018

**call toll-free
800/323-6015**

(In Illinois, call collect, 312/635-3387)

Freight Savings



Kits and generators shipped together to save shipping cost and paper work.

One order assures dependable delivery of your nuclear imaging needs... shipped together for freight savings. You also get the convenience of one-source responsibility for shipping, billing and follow-through support from the most responsive service organization in the business.

Mallinckrodt offers a full range of generator sizes and the organ-imaging kits you use most often. Find out how Mallinckrodt's efficient shipping can save time and money in your department. Call your Mallinckrodt representative or this toll-free number:

800-325-3688 (In Missouri, 314-344-3880 collect)
For technical assistance it's 800-325-8181
(In Missouri, 314-895-2405 collect)


THE MALLINCKRODT COMMITMENT

to Nuclear Medicine


MALLINCKRODT, INC., Post Office Box 5840, St. Louis, MO 63134



THE LOWEST-PRICED ERGOMETER SYSTEM ON THE MARKET ... AND ALSO THE BEST!




- Retractable Casters provide mobility and absolute stability (shown at right)
- Large Table can also be used for conventional imaging
- Left Lateral Offset eliminates table/camera head interference
- Table mounted patient RPM meter
- Adjustable Incline for semi-upright imaging (optional)



O'NEILL

O'NEILL ENTERPRISES 221 FELCH ST. ANN ARBOR, MICHIGAN 48103 (313) 973-2335

Bone



Diagnosis: hypertrophic
pulmonary osteoarthropathy

Imaging information: Instrument: GE MaxiCamera™ 535 Dose: 20 mCi OSTEOLITE
Scan time: 2.5-3.0 hours postinjection Acquisition time: 6 minutes/view

OSTEOLITE™
Technetium Tc 99m Medronate Sodium Kit (MDP)

NEN New England Nuclear®

Please see following page for brief prescribing information.

OSTEOLITE™

Tchnetium Tc 99m Medronate Sodium Kit (MDP)

INDICATIONS AND USAGE: Tchnetium Tc 99m OSTEOLITE may be used as a bone imaging agent to delineate areas of altered osteogenesis.

CONTRAINDICATIONS: None known.

WARNINGS: The contents of the OSTEOLITE vial are intended only for use in the preparation of Tchnetium Tc 99m medronate sodium and are NOT to be directly administered to the patient.

Ideally, examinations using radiopharmaceuticals — especially those elective in nature — of women of childbearing capability should be performed during the first ten days following the onset of menses.

PRECAUTIONS: A thorough knowledge of the normal distribution of intravenously administered Tchnetium Tc 99m medronate sodium is essential in order to accurately interpret pathologic studies.

Tchnetium Tc 99m medronate sodium, as well as any radioactive agent, must be handled with care. Once sodium pertechnetate Tc 99m is added to the kit, appropriate safety measures should be used to minimize external radiation exposure to clinical personnel. Care should also be taken to minimize radiation exposure to patients in a manner consistent with proper patient management.

Since 50–75% of the administered dose is renally excreted, good patient hydration and frequent voiding for 4–6 hours post-injection will significantly reduce the bladder wall dose.

The Tchnetium Tc 99m labeling reaction involved in preparing Tchnetium Tc 99m medronate sodium depends on the maintenance of tin in the divalent state. Any oxidant present in the sodium pertechnetate Tc 99m employed may adversely affect the quality of the prepared agent. Thus, sodium pertechnetate Tc 99m containing oxidants should not be used without first demonstrating that it is without adverse effect on the properties of the resulting agent.

The use of bacteriostatic sodium chloride as a diluent for sodium pertechnetate Tc 99m may adversely affect the biologic distribution of the prepared agent, and its use is not recommended.

Adequate reproduction studies have not been performed in animals to determine whether this drug affects fertility in males or females, has teratogenic potential, or has other adverse effects on the fetus. Tchnetium Tc 99m medro-

nate sodium should be used in pregnant women only when clearly needed.

It is not known whether this drug is excreted in human milk. As a general rule nursing should not be undertaken when a patient is administered radioactive material.

Safety and effectiveness in children have not been established.

ADVERSE REACTIONS: None reported.

DOSAGE AND ADMINISTRATION: The recommended dose for the average 70kg adult patient is 15mCi with a range of 10-20mCi. The patient dose should be measured by a suitable radioactivity calibration system immediately prior to administration.

Optimal imaging results are obtained within one to four hours after administration.

OSTEOLITE should be used within six hours after aseptic reconstitution with sodium pertechnetate Tc 99m. For optimum results this time should be minimized.

The vial contains no bacteriostat.

Radiopharmaceuticals should be used by persons who are qualified by specific training in the safe use and handling of radionuclides produced by nuclear reactor or particle accelerator and whose experience and training have been approved by the appropriate governmental agencies authorized to license the use of radionuclides.

HOW SUPPLIED: NEN's OSTEOLITE™ Tchnetium Tc 99m Medronate Sodium Kit is supplied as a set of five or thirty vials, sterile and non-pyrogenic. Each nitrogen-flushed vial contains in lyophilized form:

Medronate Disodium — 10mg

Stannous Chloride Dihydrate — 0.85mg

The pH is adjusted to between 7.0-7.5 with hydrochloric acid and/or sodium hydroxide solution. The contents of the vial were lyophilized under nitrogen.

Store at room temperature (15°-30° C). Included in each five (5) vial kit is one (1) package insert and six (6) radiation labels. Included in each thirty (30) vial kit is one (1) package insert and thirty-six (36) radiation labels.

The contents of the kit vials are not radioactive; however, after reconstitution with sodium pertechnetate Tc 99m the contents are radioactive and adequate shielding and handling precautions must be maintained.

Do not use if there is a vacuum in the immediate drug container or if air is injected into the container when the dose is withdrawn.

Catalog Number NRP-420 (5 vial kit)

Catalog Number NRP-420C (30 vial kit)

April 1978

Gallium Citrate Ga67

INDICATIONS AND USAGES: Gallium Citrate Ga-67 may be useful in demonstrating the presence and extent of the following malignancies: Hodgkins disease, lymphomas and bronchogenic carcinoma. Positive Ga-67 uptake in the absence of prior symptoms warrants follow-up as an indication of a potential disease state.

Gallium Citrate Ga 67 may be useful as an aid in detecting some acute inflammatory lesions.

CONTRAINDICATIONS: None known.

WARNINGS: Gallium Citrate Ga 67 should not be administered to children or to patients who are pregnant or to nursing mothers unless the information to be gained outweighs the potential hazards. Ideally, examinations using radiopharmaceutical drug products, especially those elective in nature of a woman of childbearing capability should be performed during the first few (approximately ten) days following the onset of menses.

PRECAUTIONS: A thorough knowledge of the normal distribution of intravenously administered Gallium Citrate Ga 67 is essential in order to accurately interpret pathologic studies.

The findings of an abnormal gallium concentration usually implies the existence of underlying pathology, but further diagnostic studies should be done to distinguish benign from malignant lesions. Gallium Citrate Ga 67 is intended for use as an adjunct in the diagnosis of certain neoplasms. Certain pathologic conditions may yield up to 40% false negative gallium studies. Therefore a negative study cannot be definitively interpreted as ruling out the presence of disease.

Lymphocytic lymphoma frequently does not accumulate Gallium Ga 67 sufficiently for unequivocal imaging; and the use of gallium with this histologic type of lymphoma is not recommended at this time.

Gallium Citrate Ga 67, as well as other radioactive drugs, must be handled with care and appropriate safety measures should be used to minimize external radiation exposure to clinical personnel. Also, care should be taken to minimize radiation exposure to patients consistent with proper patient management.

No long term animal studies have been performed to evaluate carcinogenic potential.

Adequate reproduction studies have not been performed in animals to determine whether this drug affects fertility in males or females, has teratogenic potential, or has other adverse effects on the fetus. Gallium Citrate Ga 67

should be used in pregnant women only when clearly needed.

Gallium Citrate Ga 67 has been found to accumulate in breast milk and should not be used in nursing mothers.

Safety and effectiveness in children have not been established.

Gallium Ga 67 localization cannot differentiate between tumor and acute inflammation; and other diagnostic studies must be added to define the underlying pathology.

The expiration date of the drug is seven days after the date of calibration.

ADVERSE REACTIONS: Severe itching, erythema and rash were observed in one patient of 300 studied.

DOSAGE AND ADMINISTRATION: The recommended adult (70kg) dose of Gallium Citrate Ga 67 is 2-5mCi. Gallium Citrate Ga 67 is intended for intravenous administration only.

Approximately 10% of the administered dose is excreted in the feces during the first week after injection. Daily laxatives and/or enemas are recommended from the day of injection until the final images are obtained in order to cleanse the bowel of radioactive material and minimize the possibility of false positive studies.

Studies indicate the optimal tumor to background concentration of ratios are often obtained about 48 hours post-injection. However, considerable biological variability may occur in individuals, and acceptable images may be obtained as early as 6 hours and as late as 120 hours after injection.

The patient dose should be measured by a suitable radioactivity calibration system immediately prior to administration.

Radiopharmaceuticals should be used by persons who are qualified by specific training in the safe use and handling of radionuclides produced by nuclear reactor or particle accelerator and whose experience and training have been approved by the appropriate governmental agencies authorized to license the use of radionuclides.

HOW SUPPLIED: Gallium Citrate Ga 67 is supplied sterile and non-pyrogenic for intravenous use. Each ml contains 2mCi of Gallium Ga 67 on the calibration date, as a complex formed from 9ng gallium chloride Ga 67, 2mg of sodium citrate, 6.8mg sodium chloride, and 0.9% benzyl alcohol w/v as preservative. The pH is adjusted to between 4.5-7.5 with hydrochloric acid and/or sodium hydroxide solution.

Vials are available from 3mCi to 18mCi in increments of 3mCi on calibration date.

The contents of the vial are radioactive and adequate shielding and handling precautions must be maintained.

Catalog Number NRP-121

December 1979

NEN New England Nuclear®

601 Treble Cove Rd., North Billerica, MA 01862


Call Toll-Free: 800-225-1572 Telex 94-0996

(In Mass. and International 617-482-9595)

Canada: NEN Canada, 2453 46th Avenue, Lachine, Que. H8T 3C9 Tel: 514-636-4971

Europe: NEN Chemicals GmbH, D-6072 Dreieich, W. Germany, Postfach 401240. Tel: (06103) 85034 Order Entry: (06103) 81011

Tumor



Diagnosis: plasmacytoma

Imaging information: Instrument: Cleon 760 Whole Body Imager
Scan time: 48 hours postinjection Speed: 5 cm/min

Dose: 5 mCi Gallium Citrate Ga 67

Gallium Citrate Ga67

 **New England Nuclear**

Please see preceding page for brief prescribing information.

SPECIAL OFFER
(see box below)

Announcing... BOOKS FROM SNM

Just Published!

NUCLEAR MEDICINE REVIEW SYLLABUS Peter T. Kirchner, M.D., Editor

Designed to help physicians bring themselves up to date in all areas of clinical practice in nuclear medicine, this brand new, 619 page book provides a thorough update on methodological advances that have occurred in nuclear medicine since the early 1970's.

The *Nuclear Medicine Review Syllabus* has chapters titled: Radiopharmacology; Instrumentation; Radiation Effects and Radiation Protection; Cardiovascular; Central Nervous System; Endocrinology; Gastroenterology; Genito-Urinary System; Hematology-Oncology; Pulmonary; Radioassay; and Skeletal System.

The clear prose of each of the book's 12 chapters describes advances and outlines current practice, with a detailed bibliography at the end of each chapter serving as a guide to additional information. A 32-page index makes the *Nuclear Medicine Review Syllabus*'s wealth of information instantly accessible. Individuals seeking a vehicle for final review prior to taking a certification (or recertification) examination will find the *Nuclear Medicine Review Syllabus* particularly valuable.

Soft cover, 619 pages, \$30.00 plus \$2.50 postage and handling.

RADIOPHARMACEUTICALS II: Proceedings of the Second International Symposium on Radiopharmaceuticals.

This 867 page, copiously illustrated, large format volume has chapters titled: Quality Control; Organic Radiopharmaceuticals; Inorganic Radiopharmaceuticals; Functional Imaging; Radioimmunoassay; Oncology; Hematology; Pharmacokinetics; Renal; Cardiopulmonary System; RES Biliary; Skeletal; Thyroid; Pancreas, Prostate, and Adrenals; and Radionuclide Production. For each of these chapters, *Radiopharmaceuticals II* has an introductory paper summarizing the state of the science in the field. The introductory papers are supplemented by papers describing current research. Also included in the book are papers from a panel discussion entitled "International Regulatory Affairs Relating to Pharmaceuticals," and excerpts from the Keynote Address given by former AEC Chairman and now Governor of the State of Washington, Dixy Lee Ray.

Soft cover, 867 pages,
\$40.00 plus \$2.50 for postage and handling.

RADIOPHARMACEUTICALS, Gopal Subramanian, Ph.D. et al, Editors.

Hardcover, 555 pages,
\$30.00 plus \$2.50 postage and handling.

THE HERITAGE OF NUCLEAR MEDICINE

Soft cover, 185 pages,
\$14.50 plus \$2.50 postage and handling.

NUCLEAR MEDICINE SCIENCE SYLLABUS

Loose-leaf plus binder, 169 pages,
\$30.50 plus \$2.50 postage and handling.

NUCLEAR CARDIOLOGY: Selected Computer Aspects

Soft cover, 213 pages,
\$12.50 plus \$2.50 postage and handling.

SPECIAL OFFER! Purchasers of *Radiopharmaceuticals II* may order *Radiopharmaceuticals* for only \$10.00 more. A \$20.00 savings! (Total cost: \$50.00 plus \$5.00 postage and handling.) Just check off "Radiopharmaceuticals Special Offer" on the coupon below.

MAIL TO: Book Order Dept., Society of Nuclear Medicine, 475 Park Avenue South, New York, NY 10016.

____ Nuclear Medicine Review Syllabus (\$30.00)
____ Radiopharmaceuticals II (\$40.00)
____ Radiopharmaceuticals (\$30.00)
____ Radiopharmaceuticals SPECIAL OFFER
(2 books \$50.00)

____ The Heritage of Nuclear Medicine (\$14.50)
____ Nuclear Medicine Science Syllabus (\$30.50)
____ Nuclear Cardiology: Selected Computer Aspects (\$12.50)

____ FOREIGN ORDERS ADD \$7.50
____ POSTAGE AND HANDLING (\$2.50 per book)
____ TOTAL ENCLOSED.

SEND TO:

NAME _____

ADDRESS _____

CITY _____ STATE _____ ZIP _____

Check or purchase order must accompany all orders. Please make checks payable to Society of Nuclear Medicine, Inc.
U.S. funds only, please.

Generator



Technetium Tc 99m Generator

NEN New England Nuclear®

ANOTHER FIRST FROM RADX...



ISOTRON

INVENTORY CONTROL COMPUTER

This small desk top microprocessor computer provides complete inventory control and NRC record keeping functions for the nuclear medicine department.

It is user programmable — you program it to fit your requirements even down to the half-life of the radionuclide so the Isotron never becomes obsolete in the rapidly changing field of nuclear medicine.

The Isotron can keep track of up to 20 different radiopharmaceuticals simultaneously by both radionuclide and chemical form! Updates the quantity of radioactivity every minute to reflect radioactivity decay.

The Isotron performs patient dose/volume calculations.

RADX gave you the first calculating dose calibrator, the first printing dose calibrator, and now the first desk top inventory control computer, the ISOTRON.



The Isotron subtracts the administered dose from the decayed activity and provides a running total of remaining activity.

The Isotron performs future time calculations. If it is 8:00 A.M. and you want to draw up a dose for 1:00 P.M. the calculation is simply and rapidly performed.

An optional hard copy data printer is available with the Isotron, known as the Isocord, which provides three copies of all pertinent data for your record keeping.

The Isotron may be used with any manufacturer's dose calibrator.


The Cost? Very reasonable. When combined with the Isocord and our Assayer 1 Dose calibrator the total price is less than competitive systems with 50% of the capabilities.

For more information or to arrange a demonstration call our toll free number 800-231-1747 (Texas customers call 713-468-9628.)

RADX

P. O. Box 19164 Houston, TX 77024

Quick Response



Call in your radiopharmaceutical orders by 2:30 PM today... get what you need by 8:00 AM tomorrow.

Through the years, Mallinckrodt's customer service system has been refined until it's second to none for prompt, dependable deliveries.

For example, our schedule is designed for your convenience... and your time zone. Just call by 2:30 PM your time, specify early delivery, and you'll have the radiopharmaceuticals you need no later than 8:00 your time the next morning.

Find out how our schedule can help keep your department on schedule. Call your Mallinckrodt representative or this toll-free number:

800-325-3688 (In Missouri, 314-344-3880 collect)

For technical assistance it's **800-325-8181**
(In Missouri, 314-895-2405 collect)

THE MALLINCKRODT COMMITMENT

to Nuclear Medicine

MALLINCKRODT, INC., Post Office Box 5840, St. Louis, MO 63134



INTRODUCING NEW MULTI-PURPOSE COLLIMATOR




PINHOLE ATTACHMENT

Announcing a new Multi-Purpose collimator developed by NSI that converts from a pinhole to an uptake collimator with a simple change of attachments.

NSI

P.O. BOX 454
BRANFORD, CONN. 06405
(212) 895-4076
(516) 752-9270



UPTAKE ATTACHMENT

NOW FULLY UTILIZE YOUR GAMMA CAMERA

- Replace your obsolete rectilinear scanner used for thyroids with our dual purpose pinhole collimator and do your thyroid uptakes and imaging on your gamma camera
- Change easily from pinhole to uptake system without removing main assembly from gamma camera
- Eliminates the need for a separate uptake system
- Pinhole attachment has three interchangeable apertures available
- Available for most gamma cameras and compatible with both high and low energy isotopes

FINALLY ...

A chair for your Gamma Camera!


Now rapid, convenient positioning can be done on ambulatory patients for brain, lung or liver scans.


Fits all CAMERAS, requires no electrical connections, firmly locks in all positions, Patient securely held with seat belt.

Enhance your current Camera investment by reducing the time required for these predominant exams.

HUMANETICS, INC.

214-242-2164 Box 185 CARROLLTON, TEXAS 75006





Preserve your copies of *The Journal of NUCLEAR MEDICINE* for years of reference in a durable, custom-designed Library Case or Binder. These storage units will hold an entire 12-issue volume. The case supplied is an attractive blue with a gold-embossed spine. Each unit also includes a gold transfer so that the volume and year can be recorded.

CASE: Holds 12 issues/\$4.95 each
three for \$14.00; six for \$24.00

BINDER: Holds 12 issue/\$6.50 each
four for \$25.00



TO: Jesse Jones Box Corp.
P.O. Box 5120 Dept. JNM
Philadelphia, PA 19141

I enclose my check or money order for \$_____ (Orders outside the U.S. add \$1.00 per file for postage and handling)

Please send me _____ *JOURNAL OF NUCLEAR MEDICINE*

_____ Files _____ Binders

Name _____

Address _____


City _____ State _____ Zip _____

Note: Satisfaction guaranteed or money refunded. Allow 5 weeks for delivery.



CMS SCINTISLICE™ TOMOGRAPHY

Multiple, simultaneous imaging.
with Software Rotation



CMS PROVIDES

Seven Pinhole and Rotating Slant Scintislice™ Tomography

- SOFTWARE FOR BOTH SYSTEMS
SPEEDS FROM 3 SECONDS / SLICE
3 DIMENSIONAL ROTATION
- COLLIMATORS
QUAD MODE PANORAMIC 7 PINHOLE
ROTATING BILATERAL & SINGLE SLANT
- UPDATE FROM EITHER SYSTEM
- CALIBRATION EQUIPMENT
- ON-SITE INSTALLATION and TRAINING

BILATERAL TOMOGRAPHY

- TWO ROTATIONS for SIX VIEWS
- IDEAL FOR LARGE or DEEP ORGANS

CMS

WRITE FOR LITERATURE

CARDIAC MEDICAL SYSTEMS CORPORATION
Main Office: 3710 Commercial Ave., Northbrook, IL 60062 USA
UK Branch: 60 Palmerston Pl., Edinburgh EH12 5AY, Scotland
Telephone: USA 312-564-4644, UK 031-225-3097

GENERIC PURCHASING ISN'T JUST FOR ASPIRIN


Nova Medical now brings the savings of generic purchasing to Nuclear Medicine.

Nova can save you as much as forty percent on the price of collimators and other Nuclear Medicine accessories, and give you a money-back guarantee.

How can we do this? Buying directly from the accessory manufacturers themselves, we can supply you with collimators and other products identical to those you now buy from scintillation camera companies. The accessory products we sell are not made by scintillation camera manufacturers. Scintillation camera manufacturers purchase these products from job shops, adding only their names and their higher overhead. If you are buying from a third party, why not choose a supplier with low overhead?

For example, there are three manufacturers presently producing collimators for scintillation camera companies. We represent each of these companies, and can provide comparable (and, in some cases, superior) collimators for forty percent less. Write or call us for a complete user's list, detailing the savings you will enjoy.

Let Nova save you money on your next purchase of nuclear medicine equipment. Our salesmen will be happy to show you the wide range of products and services, including computer products, available through Nova's generic purchasing program.



NOVA MEDICAL

Route 171 and Hughes Avenue
Lockport, Illinois 60441
815/838-7050


- OEM and specialty collimators
- Nuclear cardiology accessories
- Nuclear Tomography
- Computer products
- Nuclear Medicine accessories

Does Sulfur Colloid Preparation get you "Boiling Mad"?

If your present method of compounding Sulfur Colloid takes more than 5 minutes to prepare, you'll welcome the Diagnostic Isotopes "5-minute boil." Yes, with the DI kit for preparing Technetium Tc 99m Sulfur Colloid injection, you boil only once, and just for 5 minutes. The compounding is truly simple with the DI kit. For instance, cooling time may be reduced by placing the solution in a shielded ice bath. Everything needed for 5 multidose preparations is packaged in a ready-to-use syringe system. The DI kit comes in a convenient work station arrangement and has a one-year shelf-life.

Our discount policy has been liberalized to permit discounts on smaller quantities. Discounts may be earned by ordering as few as five of one kit or any combination of kits listed below.

Diagnostic Kits Available: Sulfur Colloid, MDP, DTPA, MAA, HSA, Diphosphonate, Polyphosphate.



Diagnostic Isotopes

SULFUR COLLOID

for imaging liver, spleen and bone marrow

di

Diagnostic Isotopes, Incorporated

225 Belleville Avenue, Bloomfield, NJ 07003

Toll free (800)631-1260, (201)429-7590

Telex 133393 Diagnostic BLFD

Our Quality Helps Your Image.

See Opposite Page For Summary
Of Prescribing Information

TECHNETIUM 99m

TSC

Kit for the preparation of
TECHNETIUM Tc 99m
SULFUR COLLOID INJECTION
DIAGNOSTIC—FOR INTRAVENOUS USE

Brief summary of prescribing information.
See package insert for full disclosure.

description

Each kit contains sufficient material to prepare (5) five formulations. Each formulation consists of a reaction vial containing 0.5 ml 1.0 N hydrochloric acid, and two syringes, one containing a 1.1 ml aqueous solution of 1.9 mg sodium thiosulfate anhydrous and the other containing 5.3 mg gelatin in 2.1 of an aqueous buffer solution containing 177 mg sodium acetate. All components are sterile and pyrogen-free. When a solution of sterile and pyrogen-free Sodium Pertechnetate Tc 99m in isotonic saline is mixed with these components, following the instructions provided with the kit, Technetium Tc 99m Sulfur Colloid Injection is formed. The product so derived is intended for intravenous injection. The precise structure of Technetium Tc 99m Sulfur Colloid Injection is not known at this time.

clinical pharmacology

Following intravenous administration, Technetium Tc 99m Sulfur Colloid Injection is rapidly cleared by the reticuloendothelial system from the blood with a nominal clearance half-life of approximately 2½ minutes. Uptake of the radioactive colloid by organs of the reticuloendothelial system is dependent upon both their relative blood flow rates and the functional capacity of the phagocytic cells. In the average patient 80 to 90% of the injected colloidal particles are phagocytized by the Kupffer cells of the liver, 5 to 10% by the spleen and the balance by the bone marrow.

indications and usage

Technetium Tc 99m Sulfur Colloid Injection is used as an agent for imaging areas of functioning reticuloendothelial cells in the liver, spleen, and bone marrow.

contraindications

None known.

warnings

The contents of the two syringes, one syringe containing the sodium thiosulfate solution and the second syringe containing the appropriate buffer solution, are intended *only* for use in the preparation of the Technetium Tc 99m Sulfur Colloid Injection and *are not to be directly administered to the patient.*

adverse reactions

Hypersensitivity reactions, including anaphylaxis, have been reported in patients receiving sulfur colloid preparations. One death and several cases of lung and soft tissue uptake

other than RES have been reported in the association with the use of Technetium Tc 99m Sulfur Colloid Injection

dosage and administration

The suggested intravenous dose range used in the average patient (70 kg) is 1 to 8 millicuries of Technetium Tc 99m Sulfur Colloid Injection.

When orally administered, the Technetium Tc 99m Sulfur Colloid Injection is not absorbed from the G.I. tract.

The patient dose should be measured by a suitable radioactivity calibration system immediately prior to administration.

Radiopharmaceuticals should be used only by physicians who are qualified by specific training in the safe use and handling of radionuclides produced by nuclear reactor or particle accelerator and whose experience and training have been approved by the appropriate government agency authorized to license the use of radionuclides.

radiation dosimetry

The estimated absorbed radiation doses ⁽²⁾ to an average patient (70 kg), or to patients with diffuse parenchymal liver disease from an intravenous injection of a maximum dose of eight millicuries of Technetium Tc 99m Sulfur Colloid Injection are shown in Table IV.


table IV. radiation doses

tissue	absorbed radiation doses (rads) 8 millicuries Tc 99m		
	normal liver	diffuse parenchymal disease early- intermediate	intermediate advanced
Liver	2.7	1.7	1.3
Spleen	1.7	2.2	3.4
Bone			
Marrow	0.22	0.36	0.63
Testes	0.0088	0.017	0.026
Ovaries	0.045	0.065	0.096
Whole-body	0.15	0.15	0.14

(2) Modified from: Summary of Current Radiation Dose Estimates to Humans with Various Liver Conditions from ^{99m}Tc-Sulfur Colloid. MIRD Dose Estimate Report No. 3, *J. Nucl. Med.*, 16, No. 1, 108A-B (1975).

manufactured for:
diagnostic isotopes incorporated
225 belleville ave./bloomfield, n.j. 07003

Always Available



**On call 24 hours a day
to take your order.**

You can place an order with Mallinckrodt any time of the day or night, 7 days a week. Our order department is open a full ten hours every business day. But you can call even after hours or on weekends and your order will be recorded. Since our people are scheduled to come in early, order processing starts before normal business hours.

Ours is a time-proven system which processes and ships literally hundreds of orders a day. And isn't it nice to know you can always get through to Mallinckrodt should you need to order a technetium Tc99m generator, or radiopharmaceuticals?

For more information, call your Mallinckrodt representative or this toll-free number:

800-325-3688 (In Missouri, 314-344-3880 collect)

For technical assistance it's 800-325-8181
(In Missouri, 314-895-2405 collect)

THE MALLINCKRODT COMMITMENT



to Nuclear Medicine

MALLINCKRODT, INC., Post Office Box 5840, St. Louis, MO 63134

GENERATOR TECHNETIUM Tc 99m

GENERATOR FOR THE PRODUCTION OF SODIUM PERTECHNETATE Tc 99m

description—The Union Carbide TECHNETIUM Tc 99m Generator provides a means of obtaining a sterile, pyrogen-free solution of Sodium Pertechnetate Tc 99m in isotonic saline from elution of the generator containing Molybdenum Mo 99. Hydrochloric acid and/or sodium hydroxide may have been used for pH adjustment. The carrier-free solution may be used as is, or with proper dilution to prepare the studies described herein. Over the life of the generator, an elution will contain a yield of 80% to 100% of the theoretical amount of Technetium Tc 99m available from the Molybdenum Mo 99 on the generator column.

physical characteristics

Technetium Tc 99m decays by isomeric transition with a physical half-life of 6.02 hours. Photons that are useful for imaging studies and the principle radiations contributing to the internal dose rate are listed in Table I.

Table I. principle radiation emission data

radiation	mean %/disintegration	mean energy (keV)
Gamma-2	88.96	140.5

*Martin, M.J., ed., Nuclear Decay Data for Selected Radionuclides, ORNL-5114, p. 24, March 1976.

external radiation

The specific gamma ray constant for Technetium Tc 99m is 0.8 R/millicurie-hour at 1 cm. The first half-value thickness of lead (Pb) for Technetium Tc 99m is 0.2 mm. A range of values for the relative attenuation of the radiation emitted by this radionuclide that results from interposition of various thicknesses of Pb is shown in Table II. For example, the use of 2.5 mm of Pb will decrease the external radiation exposure by a factor of 1,000.

Table II. radiation attenuation by lead (Pb) shielding

shield thickness (Pb) mm	coefficient of attenuation
0.2	0.5
0.8	10 ⁻¹
1.6	10 ⁻²
2.5	10 ⁻³
3.3	10 ⁻⁴

Molybdenum Mo 99 decays to Technetium Tc 99m with a Molybdenum Mo 99 half-life of 2.75 days. The physical decay characteristics of Molybdenum Mo 99 are such that only 86.8% of the decaying Molybdenum Mo 99 atoms form Technetium Tc 99m. Generator elutions may be made at any time, but the amount of Technetium Tc 99m available will depend on the interval from the last elution. Approximately 47% of maximum Technetium Tc 99m is reached after 6 hours and 95% after 24 hours. To correct for physical decay of each radionuclide, the fractions that remain at selected intervals of time are shown in Table III.

Table III. physical decay chart

Molybdenum Mo 99 half-life 2.75 days		Technetium Tc 99m half-life 6.02 hours	
days	fraction remaining	hours	fraction remaining
0	1.000	0	1.000
1	.777	1	.891
2	.604	2	.784
3	.469	3	.708
4	.365	4	.631
5	.284	5	.562
6	.220	6	.501
7	.171	7	.447
8	.133	8	.398
9	.103	9	.355
10	.080	10	.316
11	.063	11	.282
12	.049	12	.251
13	.038		

*Calibration time.

clinical pharmacology—Following intravenous administration, the pertechnetate ion distributes in the body similarly to the iodide ion, but it is not organified when trapped in the thyroid gland. Sodium Pertechnetate Tc 99m tends to accumulate in intracranial lesions with excessive neovascularity or an altered blood-brain barrier. It also concentrates in the thyroid gland, stomach and choroid plexus.

After intravenous administration, it remains in the circulatory system for sufficient time to permit blood pool, organ perfusion, and major vessel studies. It gradually equilibrates with the extracellular space. A fraction is promptly excreted via the kidneys.

indications and usage—Sodium Pertechnetate Tc 99m is used as an agent for brain imaging, thyroid imaging, salivary gland imaging, placenta localization and blood pool scans.

contraindications—None known.

warnings—This radiopharmaceutical preparation should not be administered to children or to patients who are pregnant or to nursing mothers unless the expected benefits to be gained outweigh the potential hazards.

Ideally, examinations using radiopharmaceuticals, especially those elective in nature, of a woman of childbearing capability should be performed during the first few (approximately 10) days following the onset of menses.

precautions—Sodium Pertechnetate Tc 99m, as well as other radioactive drugs, must be handled with care and appropriate safety measures should be used to minimize external radiation exposure to clinical personnel. Also, care should be taken to minimize radiation exposure to patients consistent with proper patient management.

Adequate reproduction studies have not been performed in animals to determine whether this drug affects fertility in males or females, has teratogenic potential, or has other adverse effects on the fetus. Sodium Pertechnetate Tc 99m should be used in pregnant women only when clearly needed.

It is not known whether this drug is excreted in human milk. As a general rule, nursing should not be undertaken while a pa-

tient is on a drug since many drugs are excreted in human milk.

Safety and effectiveness in children have not been established.

The generator should not be used after 16 days from the date and time of calibration.

At time of administration, the solution should be crystal clear.

adverse reactions—No adverse reactions have been reported with the use of this radiopharmaceutical.

dosage and administration—Sodium Pertechnetate Tc 99m is usually administered by intravenous injection, but can be given orally. The dosage employed varies with each diagnostic procedure. The suggested intravenous dose range employed in the average adult (70 kg) in millicuries of Sodium Pertechnetate Tc 99m for various diagnostic indications is as follows:

Brain Scan	10 to 20 millicuries
Thyroid Gland Scan	1 to 10 millicuries
Salivary Gland Scan	1 to 5 millicuries
Placenta Localization	1 to 3 millicuries
Blood Pool Scan	10 to 20 millicuries

NOTE: Up to 1 gram of reagent grade potassium perchlorate in a suitable base or capsule may be given orally prior to administration of Sodium Pertechnetate Tc 99m injection for brain scan, placenta localization and blood pool scan for the purpose of blocking uptake of Sodium Pertechnetate Tc 99m by the choroid plexus.

The patient dose should be measured by a suitable radioactivity calibration system immediately prior to administration.

Radiopharmaceuticals should be used only by physicians who are qualified by training and experience in the safe use and handling of radionuclides and whose experience and training have been approved by the appropriate government agency authorized to license the use of radionuclides.

radiation dosimetry

The estimated absorbed radiation doses¹ to an average patient (70kg) from an intravenous injection of a maximum dose of 20 millicuries of Sodium Pertechnetate Tc 99m distributed uniformly in the total body of subjects not pretreated with blocking agents such as NaClO₄, KClO₄, or iodide are shown in Table IV. For placental localization studies when a maximum dose of 3 millicuries is used it is assumed to be uniformly equilibrated between maternal and fetal tissues.

Table IV. radiation doses

tissue	absorbed radiation dose (rads/20 millicuries)		(rads/3 millicuries)
	Resting Population	Active Population	
Bladder wall	1.06	1.70	
Gastrointestinal tract			
Stomach wall	5.00	1.02	
Upper large intestine wall	1.36	2.40	
Lower large intestine wall	1.22	2.20	
Red marrow	0.38	0.34	
Testes	0.18	0.18	
Ovaries	0.44	0.60	
Thyroid	2.60	2.60	
Whole-body	0.28	0.22	
* Brain	0.28	0.24	
* Placenta			0.05
* Fetus			0.05

*Method of Calculation: A Schema for Absorbed Dose Calculations for Biologically Distributed Radionuclides, Supplement No. 1, MIRD Pamphlet No. 1, J. Nucl. Med., p. 7 (1968).

*Summary of Current Radiation Dose Estimates to Normal Humans From 99mTc as Sodium Pertechnetate, MIRD Dose Estimate Report No. 8, J. Nucl. Med., 17:1, 1976.

Table V. Generator dosimetry readings

Generators up to 4140 mCi internal lead shield			Generators 4970 to 16600 mCi internal depleted uranium shield		
days from calibration	mR/hr 2" 12"	mCi 99Mo	days from calibration	mR/hr 2" 12"	mCi 99Mo
0*	425 57	4410	0*	174 33	16800
1	330 44	3430	1	135 26	13100
2	256 34	2660	2	105 20	10200
3	199 27	2070	3	81 16	7900
4	155 21	1610	4	63 12	6100
5	120 16	1250	5	49 9	4800
6	94 12	970	6	38 7	3200
7	73 10	750	7	30 6	2900

*Day of calibration at 12:00 hrs E.T. is the day of shipment from Tuxedo, N.Y.

Table VI. elution vial radiation dosimetry

vial distance from probe	dosimetry		dosimetry shielded vial*
	contact	bare vial	
30.5 cm	47200mR/hr	13000mR/hr	4 mR/hr 0.8mR/hr

*Union Carbide Elution Vial Shield Cat. No. 17500500, Shield 6.35mm Lead.

how supplied—Sodium Pertechnetate Tc 99m is supplied as a Molybdenum Mo 99/Technetium Tc 99m generator in sizes from 830 millicuries up to 16,600 millicuries (in approximately 830 millicurie increments) of Molybdenum Mo 99 as of noon of the day of calibration. The TECHNETIUM Tc 99m GENERATOR consists of:

- 1) sterile generator, 2) Sodium Chloride Injection source, 3) 10 cc sterile evacuated vials, 4) sterile needles, 5) elution vial shield* 6) finished drug labels. Elution vials in 5 cc and 20 cc sizes are available upon request.

*initial order only.

preparation


The following instructions must be carefully followed for optimum preparation of Sodium Pertechnetate Tc 99m.

Union Carbide Generators are sterile and pyrogen-free at the time of shipment. Aseptic technique must be observed during the use of the generator to maintain a sterile and pyrogen-free system. Gloves should be worn during all elution procedures.

The sealed column and fluid path MUST NOT be removed from the shielding system.
*CAUTION: It is recommended that elution vial shields be used when eluting the generators, shielded syringes be used when preparing formulations, and appropriate vial shields be used for the formulations.

First Elution

1. Remove generator system and accessories from carton.
 2. Lift hinged cover exposing dispenser end. Remove protective cap from dispenser end and attach a sterile needle—REMOVE PLASTIC NEEDLE COVER (Figure 2). Return cover to closed position.
 3. Place an elution vial in the elution shield (Figure 1) and clean septum of elution vial with an antiseptic swab. Position elution shield on dispensing platform (Figure 3).
 4. Rotate fluid path shut off valve several full turns counterclockwise until loose. Valve is located on left side of generator.
 5. Slide elution shield to far left position (Figure 4). The dispensing needle will pierce the septum of the evacuated elution vial. The elution will begin immediately.
 6. Step away to reduce your radiation exposure. Allow 3 to 5 minutes for complete elution.
- NOTE: If vacuum in elution vial is lost, i.e., no eluate present in vial, discard vial and use a new elution vial.
7. When elution is complete, slide elution shield to far right position. Remove elution shield, containing vial with Sodium Pertechnetate Tc 99m eluate, from dispensing platform.
 8. Replace dispensing needle with sterile needle with plastic cover in place. DO NOT REMOVE COVER FROM NEEDLE until next elution.
 9. Affix the pressure-sensitive label to the dose vial shield. Sodium Pertechnetate Tc 99m is ready for use. Maintain adequate shielding of the radioactive preparation.



storage

Store generator at room temperature (18-25 °C).

Caution: Avoid Freezing.

subsequent elutions

1. Lift hinged cover exposing dispenser needle. Remove plastic needle cover from dispensing needle and discard. Return cover to closed position.
2. Repeat steps 3, 5, 6, 7, 8 and 9.

20 ml elutions—To use the larger size elution vial, remove the spacer in the elution shield and replace with the spacer designed for 20 cc vials.

The radioactivity concentration of the final Sodium Pertechnetate Tc 99m preparation may be calculated by using the following formula:

C = AV where C equals radioactivity concentration of the Sodium Pertechnetate Tc 99m preparation (millicuries/ml),
A = Technetium Tc 99m activity added to the reaction mixture vessel (millicuries),
V = Total volume in the final mixture (ml).

Technetium Tc 99m assay procedure

1. Determine the equivalent Technetium Tc 99m value for a Cobalt Co 57 standard by multiplying the number of millicuries of Cobalt Co 57 standard by the appropriate equivalent factor. This equivalent value of Cobalt Co 57 for the standard need only be decayed daily for use as a secondary standard.
2. Place the standard in the chamber and record μ amp reading.
3. Transfer the Technetium Tc 99m sample from the shield to the chamber. Record the μ amp reading.
4. Calculate activity:
 μ amps of Tc 99m Sample / μ amps of ⁶⁰Co std. = x millicuries Cobalt Co 57 std. = millicuries Technetium Tc 99m
where millicuries Cobalt Co 57 std. = the equivalent millicurie value for Cobalt Co 57 from 1. above, corrected for decay.

direct readout procedure—A direct readout dose calibrator is used.

1. Determine the equivalent millicurie Technetium Tc 99m value for a Cobalt Co 57 std. using method 1. above. Correct millicurie value for decay.
2. Place Cobalt Co 57 standard in chamber and adjust the calibrator to the proper reading according to the manufacturer's instructions.
3. Transfer sample vial to chamber and read directly millicuries Technetium Tc 99m.

Molybdenum Mo 99 breakthrough test

1. Determine the amount of Technetium Tc 99m eluted (millicuries).
 2. Place the Technetium Tc 99m elution in a lead container. Place lid on container and put the entire container in the chamber.
 3. Record the amount of Molybdenum Mo 99 (microcuries) on the most sensitive scale.
 4. Divide the microcuries Molybdenum Mo 99 by the millicuries Technetium Tc 99m. Correct for decay and shielding effect, if necessary.
- The acceptable limit is 1.0 microcurie Molybdenum Mo 99/millicurie Technetium Tc 99m, not to exceed 5 microcuries per human dose at the time of injection.

disposal

The TECHNETIUM Tc 99m GENERATOR should not be discarded in ordinary trash within 70 days of the calibration date. Vials and needles used for eluting may be discarded after two (2) days. It is suggested that all identification labels be destroyed before discarding the generator or vials.

TECHNETIUM Tc 99m GENERATORS OF \leq 4140 millicuries may be returned to the manufacturer; while those of 4970 to 16,600 millicuries must be returned to the manufacturer. Please refer to the instructions included with each shipment.


This generator is approved for use by persons licensed by the U.S. Nuclear Regulatory Commission pursuant to Sec. 35.14 and Sec. 35.100 Group III of 10 CFR Part 35 or under equivalent licenses of Agreement States.

New High Activity Union Carbide CintiChem®

Technetium 99m Generators

(Technetium Tc 99m Generators for the Production of Sodium Pertechnetate Tc 99m)

AVAILABLE EXCLUSIVELY FROM UNION CARBIDE



INCORPORATE THE FOLLOWING FEATURES:

COST EFFECTIVENESS

The use of a single high activity CINTICHEM® Technetium 99m Generator can significantly reduce your operating expenses by replacing from two to five lower activity Technetium 99m Generators.

CALIBRATIONS IN Mo 99

8,280 mCi to 16,600 mCi calibration sizes, in increments of approximately 830 mCi.

Tuesday and Thursday calibration days.

SAFETY

The depleted uranium internal shielding of the column possesses greater density and therefore superior shielding properties than the lead shielding used in lower activity generators.

The CINTICHEM® 8,280 mCi Technetium 99m Generator provides an approximate dosimetry of 9.7 mr/hr on the day of calibration, following elution, 18 inches from the surface, *without* a secondary external shield.¹

The weight of a 16,600 mCi Generator is 45.5 pounds. All high activity generators are shipped in the D.O.T. Transport Index II classification.

QUALITY

A special glass column design incorporated in the high activity CINTICHEM® Technetium 99m Generators provides for high yields with as low as a 5 cc elution volume. Furthermore, the specially designed column reduces the potential Al + + + and Mo 99 content in the eluate.

VERSATILITY

High activity CINTICHEM® Technetium 99m Generators can reduce the time required to perform quality control because in each single elution levels of activity are provided that would require eluting several lower activity Generators, and quality controlling each eluate.

CONVENIENCE

High activity CINTICHEM® Technetium 99m Generators can dramatically reduce shielded shelf space requirements.

High activity CINTICHEM® Technetium 99m Generators can eliminate long term decay storage.

¹Data on file at Union Carbide Corporation, Tuxedo, New York, and with the State of New York, Division of Safety and Health.

For full prescribing information, refer to preceding page.



UNION CARBIDE...INVOLVED IN NUCLEAR MEDICINE FOR OVER 19 YEARS

From Atom to Image


UNION CARBIDE CORPORATION • P.O. BOX 324 • TUXEDO, NEW YORK 10987
FOR PRODUCT INFORMATION CALL TOLL FREE (800) 431-1146. IN N.Y.S. CALL (800) 942-1986

CintiChem®

Technetium 99m Generators

(Technetium Tc 99m
Generators for the
Production of Sodium
Pertechnetate Tc 99m)

ARE DESIGNED TO MAXIMIZE
RADIATION PROTECTION AND
FOR EASY ELUTION



Rapid Elution
Vial-Needle
Engagement
Reduces the
Radiation Exposure
Time Factor

INCORPORATE THE FOLLOWING ADVANTAGES:

- Only UNION CARBIDE CINTICHEM® Technetium 99m Generators are made at one domestic production site that possesses its own Nuclear Reactor for the production of Fission Product Mo 99, manufactures and purifies by a patented process high specific activity Mo 99, loads it onto columns, assembles the Generators, performs quality control, and ships them directly to the user. This provides you with a reliable product supply and a uniformly high-quality product.
- The UNION CARBIDE Fission Product Mo 99 used in CINTICHEM® Technetium 99m Generators provides Sodium Pertechnetate Tc 99m activity concentrations sufficient for bolus injections.
- CINTICHEM® Technetium 99m Generators come in 32 activity and day of calibration combinations, which can satisfy a wide range of activity needs.
- A new sterile needle is utilized for each elution, reducing the chances of a septic or pyrogenic situation occurring in routine clinical usage. This method offers an advantage compared to competitive dry column systems where the same needle assembly is used for the life of the product.
- Rigid Quality Control testing, which includes an elution check on each Generator, assures that your UNION CARBIDE CINTICHEM® Technetium 99m Generator meets our high internal specifications. Our experience obtained in over 19 years of involvement in Nuclear Medicine assures you of high quality products.

For full prescribing information, refer to page opposite preceding page.



UNION CARBIDE...INVOLVED IN NUCLEAR MEDICINE FOR OVER 19 YEARS

From Atom to Image

UNION CARBIDE CORPORATION • P.O. BOX 324 • TUXEDO, NEW YORK 10987
FOR PRODUCT INFORMATION CALL TOLL FREE (800) 431-1146. IN N.Y.S. CALL (800) 942-1986

OUR CRC-30 RADIOISOTOPE CALIBRATOR. IT'S PART EINSTEIN, PART FREUD, AND PART GUTENBURG.




The CRC-30 calibrates and computes, analyzes radiochemical purity, and puts it all in print.

Computes radiopharmaceutical dose to assure that activity is exactly as prescribed.

Analyzes imaging preparations to assure radiochemical purity.

Prints permanent records in triplicate. Gives molybdenum assay printout.



Simplifies compliance with regulatory and hospital accreditation standards.

Capintec, Inc., 136 Summit Avenue,
Montvale, NJ 07645. (201) 391-3930.
Toll Free (800) 631-2557. TELEX 642375
CAPINTEC MTLE.


CI CAPINTEC
THE MEASURE OF EXCELLENCE

Gamma Camera Upgrade

1/10 Inch or Better Resolution at a fraction of new system cost.

BEFORE


NEN Thallium 201 phantom at 2" distance from collimator.
500K



Picker 2C with ultrafine collimator.

AFTER

NEN Thallium 201 phantom at 2" distance from collimator.
500K



Picker 2C with ultrafine collimator and SX-11 detector head.

The picture on your left does not provide adequate resolution for cardiac work. The picture on the right is more than adequate!

Enjoy new camera performance without a major investment. Nuclear Service Inc. can upgrade your existing gamma camera system to provide you with 1/10" or better intrinsic resolution.

With NSI upgrade you not only receive State-of-the-Art resolution, but in most cases your

converted system will provide you with **IMPROVED FIELD UNIFORMITY.**

Learn more about this efficient and economical method available from one of the country's largest independent nuclear medicine service organizations. Call or write NSI for complete infor-

mation on gamma camera upgrade.

Up to 75% Better Resolution.

Picker	Improved Resolution
2	75%
2C	50%
3C, 4-12	40%
1/8"	20%

*Leasing plans and reconditioned upgraded systems also available.

Nuclear Services Inc.



Nuclear Services Inc.
P.O. Box 454
Branford, CT 06405

(516) 752-9270
(212) 895-4076

PLACEMENT

POSITIONS OPEN

ACADEMIC POSITION AT THE ASSOCIATE or Assistant Professor level available in the Nuclear Radiology Division of the Department of Radiology at the University of Texas Medical School at Houston. Certification in Radiology and Nuclear Medicine, or in Radiology with Special Competence in Nuclear Radiology is required. Applicant should have a sincere interest and a performance record in relevant clinical or basic nuclear research. Please send curriculum vitae to Robert W. McConnell, M.D., Director, Division of Nuclear Radiology, Department of Radiology, The University of Texas Medical School at Houston, 6431 Fannin Street, Houston Texas 77030.

NUCLEAR MEDICINE TECHNOLOGIST. Challenging position for registered technologist in progressive dept. Must have a working knowledge of nuclear imaging and radioimmunoassays. Good salary and superior benefits. Send resume or contact Director, Employee Relations, University Community Hospital, 3100 E. Fletcher Avenue. Tampa, Florida 33612.

NUCLEAR MEDICINE RESIDENCY Available July 1981. Two-year accredited affiliated program including 700 bed VA General Hospital, 400 bed County Hospital and 1,000 bed Air Force Medical Center; an equal opportunity employer. Comprehensive training in basic sciences, laboratory sciences, computer technology, patient care services, and research. Contact Martin L. Nusynowitz, M.D., Division of Nuclear Medicine, University of Texas Health Science Center, San Antonio, Texas 78284. (512) 691-7265.

NUCLEAR MEDICINE TECHNOLOGIST Full time position for registered or registry eligible Nuclear Medicine Technologists to join Nuclear Medicine Department of progressive 286 bed hospital in sunny Arizona. Full range of in-vivo procedures, 2 gamma cameras with computer. Good salary and employee benefits. Send resume to: Personnel Dept., Mesa Lutheran Hospital, 525 West Brown Rd., Mesa, Arizona 85201.

CALIFORNIA-EXPERIENCED NUCLEAR Medicine Technologists are invited to call us collect regarding immediate career opportunities at our 540-bed, acute care, university affiliated teaching hospital. In addition to excellent salary and benefits, we provide the finest in equipment and a very professional work environment. Call or submit resume to: Ms. Lynne Luboviski, Employment Manager, St. Mary Medical Center, 1050 Linden Avenue, Long Beach, CA 90801. (213) 435-4441, Ext. 420. An equal opportunity employer-male/female/handicapped.

NUCLEAR MEDICINE RESIDENT. Two-year residency program in Nuclear Medicine at The New York Hospital-Cornell Medical Center. Positions available July 1, 1981. Contact Jerome G. Jacobstein, M.D., Division of Nuclear Medicine, The New York Hospital-Cornell Medical Center, 525 East 68th St., New York, NY 10021.

RADIOLOGIST, BOARD CERTIFIED IN Nuclear Medicine, to join large multi-specialty prepaid medical group. Opportunity to expand department and plan department for new hospital in 1984. Salary negotiable. Liberal fringe benefits. Contact: Hawaii Permanente Medical Group, Inc., 1697 Ala Moana Boulevard, Honolulu, HI 96815. an equal opportunity employer.

NUCLEAR MED TECH-ARE YOU LOOKING for a challenging position in the field of nuclear medicine? If you are registered or registry eligible and need a full time position, then this may be the opportunity you've been looking for. This is a beautiful university community with several lakes and parks in close proximity. Bloomington Hospital is a 314 bed Hospital that services Bloomington and the surrounding areas. If interested, please call or write: Bloomington Hospital, 619 W. First St., Bloomington, IN 47401. (812) 336-9535. An Equal Opportunity Employer M/F.

DIRECTOR OF SECTION OF NUCLEAR Medicine in 500 bed University affiliated hospital in Metropolitan Chicago area. Newly remodeled, fully equipped lab with latest in computer and other equipment. Several years of experience required with strong background in Nuclear Cardiology and certification American Board of Nuclear Medicine. C.V. to Box 1002 Society of Nuclear Medicine, 475 Park Ave. South, New York, N.Y., 10016.

NUC. MED. TECHNOLOGIST POSITIONS open immediately at the University of Texas Medical Branch in Galveston, Texas. UTMB is a 1200 bed teaching hospital system that performs the latest of procedures with modern equipment. We are a rapidly growing organization that offers opportunities for professional growth and advancement. Applicants should have formal training and be registry eligible. Salary commensurate with background, training and experience. Benefits include 15 days of paid vacation, paid sick leave, state pays most of Social Security Tax, no state income tax plus many others. Casual living in a semitropical climate. Send resume to William K. Otte, Jr., Assistant to the Chairman, Department of Radiology, P.O. Box 146, University of Texas Medical Branch, Galveston, Texas 77550. For more direct information call collect (713) 765-2009. UTMB is an affirmative action/equal opportunity employer.

NUCLEAR MEDICINE TECHNOLOGIST Immediate opening in our growing Nuclear Medicine Department for a registered Nuclear Medicine Technologist or registry eligible. Experience with computers and nuclear cardiology desirable. At our modern Midwest hospital we offer a competitive salary and comprehensive benefit package. For prompt consideration please contact Teresa Jenkins, Personnel, St. Mary's Medical Center, 3801 Spring St., Racine, WI 53405. (414) 636-4495.

NUCLEAR MEDICINE TECHNICIAN A large 575 bed teaching hospital, 10 miles south of Pittsburgh has an immediate opening for a Nuclear Medicine Technician (registered or eligible). Excellent salary and fringe benefit program. For interview contact: McKeesport Hospital, 1500 5th Avenue, McKeesport, PA 15132. Phone (412) 664-2415. An equal opportunity/affirmative action Employer.

NUCLEAR MEDICINE TECHNOLOGIST Immediate full-time opening for experienced technologist. Excellent opportunity with progressive mobile service company. Send resume to: Universal Nuclear, Inc., 6900 Knott Avenue, Suite F, Buena Park, CA 90621. Telephones: (714) 523-5431 or (213) 921-7931.

FELLOWSHIP PROGRAMS IN DIAGNOSTIC Imaging. One (ABR) or two (ABNM) year approved programs offer pediatric and adult experience in all aspects of imaging including Ultrasound, CT, and Cardiovascular Nuclear Medicine, and training in computers, radioassay and basic sciences. Facilities include state of the art instrumentation, centralized nuclear pharmacy and well equipped in vitro and animal laboratories for excellent research opportunities. Available July 1, 1981. Contact: K.M. Spicer, M.D., Ph.D., Director, Division of Diagnostic Imaging, Medical University of South Carolina, 171 Ashley Avenue, Charleston, S.C. 29403.

NUCLEAR PHYSICIAN-JULY 1981. Diagnostic Imaging ABNM, ABR program seeking recently trained nuclear medicine physician to assist one full-time and two part-time staff with educational and clinical responsibilities in Medical University setting. Excellent research opportunities and facilities available. Contact: K.M. Spicer, M.D., Ph.D., Director, Diagnostic Imaging, Medical University of South Carolina, 171 Ashley Avenue, Charleston, S.C. 29403.

FLORIDA-NUCLEAR MEDICINE TECH- nologist. Immediate opening presently exists in our 714 bed non-profit, community hospital, located between Miami and Ft. Lauderdale. Individual must be registered or registry eligible. Hospital experience preferred. Excellent benefit program. Salary commensurate with experience and education. Send resume or call COLLECT (305) 987-2000 Ext. 5470 for Ann Cunniff, Employment Manager, Memorial Hospital, 3501 Johnson Street, Hollywood, FL. 33021. Equal Opportunity Employer.

NUCLEAR RADIOLOGIST-NEEDED for expanding University affiliated Radiologic group. Experience and/or board certification desired. Opportunity to share other imaging responsibilities. Excellent opportunity for a young aggressive person with primary interest in Nuclear Radiology. Northeast Ohio Region. Reply Box 1003, Society of Nuclear Medicine, 475 Park Ave. South, New York, NY 10016.

NUCLEAR MEDICINE TECHNOLOGIST. Full-time opening for a Registered Nuclear Medicine Technologist in a 715 bed critical-care hospital. Strong background in nuclear cardiology and computers preferred. For further information, contact Ann Dungan, Personnel Dept., St. Anthony Hospital Systems, 4231 W. 16th Ave., Denver CO 80204; Phone (303) 629-3700.

ASSISTANT PROFESSOR OF NUCLEAR Pharmacy/Health Physics at Purdue University. Responsibilities include research and teaching in both areas. Ph.D. in nuclear pharmacy or health physics required. Preference will be given to individuals with appropriate board certification in health physics or licensure in nuclear pharmacy. Application deadline November 1, 1980. Contact Paul L. Ziemer, Associate Head, Department of Bionucleonics, Purdue University, West Lafayette, IN 47909. An Equal Opportunity/Affirmative Action Employer.

RADIOLOGICAL PHYSICIST: A PHYSICIST at the M.S. or B.S. level wishing to practice all aspects of applied radiological physics required for a professional consulting firm. Experience in radiation therapy treatment planning desired. Our activities include service in (1) Radiation Therapy Planning, (2) Nuclear Medicine, (3) Diagnostic Radiology, and (4) Radiation Protection. Excellent salary and benefits. Submit resume to Box 1001, Society of Nuclear Medicine, 475 Park Avenue South, New York, NY 10016.

NUCLEAR MEDICINE RESIDENCY: THE Washington, D.C. Veterans Administration Medical Center is offering a 2 year A.M.A. approved program affiliated with George Washington University beginning July 1981. The center is a 700 bed general medical and surgical hospital. The program includes training in Radionuclide in vivo and invitro procedures computer application as well as diagnostic ultrasound. At least 2 years of prior training in radiology, internal medicine or pathology is required. Equal Opportunity Employer. Contact: B.J. Sauerbrunn, M.D., Chief, Nuclear Medicine Services, VA Medical Center, 50 Irving St., N.W., Washington, D.C. 20422.

POSITIONS WANTED

ABNM CERTIFIED INTERNIST WITH one year training in general diagnostic radiology seeks directorship of nuclear medicine section (private practice). Eight years experience as director of an active department of nuclear medicine and diagnostic ultrasound. Experienced in both clinical and research aspects of nuclear medicine imaging, including cardiovascular and computer applications, RIA, and B-mode ultrasound. All positions will be considered. For further information and C.V. reply Box 1000, Society of Nuclear Medicine, 475 Park Ave., So., New York, NY 10016.

FOR SALE

SEARLE LEM MOBILE CAMERA WITH new heavy-duty power unit and 1977 Ford 1 ton lift gate truck for transport. Camera has latest performance specs. Will sell separately. Make offer. K. Meier (415) 664-7400.

R.I.A. PROFICIENCY TESTING PROGRAM

RAS-1
(\$100/yr.)

Single vial providing 5 ml. when reconstituted
Constituents:

Cortisol, Digoxin, Triiodothyronine (T3), T3 uptake (developmental), Thyroxine (T4), Free T4, Compensated T4 (developmental), Thyroid Stimulating Hormone (TSH), Thyroxine Binding Globulin (TBG), Insulin, Human Growth Hormone (HGH), Follicle Stimulating Hormone (FSH), Luteinizing Hormone (LH), Folic Acid, Vitamin B-12, Gastrin, Ferritin, Tobramycin.

RAS-2
(\$140/yr.)

Two identical vials, each providing 5 ml. when reconstituted. Constituents identical to RAS-1.

Shipped Quarterly
Enrollment Accepted Now

For information call (202) 857-1135 or write:

**AMERICAN COLLEGE
OF NUCLEAR PHYSICIANS**
1101 Connecticut Avenue N.W. Box 831
Suite 700
Washington, D.C. 20036

Pharmatopes, Inc.
NUCLEAR PHARMACY SERVICES


NUCLEAR PHARMACISTS

We are one of the leaders in this exciting field and now have positions available nationwide for managers and staff nuclear pharmacists. We are especially interested in R. Ph.'s for the northeast & midwest area.

We offer the best incentive programs for managers, excellent fringe benefits and salaries commensurate with experience.

SEND RESUME TO: PHARMATOPES, INC.
25721 COOLIDGE
OAK PARK, MI 48237
OR CALL 313-543-8400 EOE.

**NUCLEAR MEDICINE
TECHNOLOGIST
MIAMI, FLORIDA**

Cedars of Lebanon, a 700-bed, private, non-profit health care center, is seeking an ARRT or NMTCB registered Senior Staff Technologist with a minimum of two years experience. Modern equipment including complete cardiovascular phase and portable camera.

Our specialty referral facility is located in the heart of Miami's medical professional center within easy access to residential, beach and cultural areas.

We offer excellent starting salary, comprehensive benefits and relocation assistance program. For immediate consideration, please send resume, including salary history, in strict confidence to:



PERSONNEL DEPARTMENT
**cedars of lebanon
health care center**

1295 Northwest 14th St. Miami, Florida 33136
An equal opportunity employer M/F/H

AN OK OPPORTUNITY

**NUCLEAR MEDICINE
TECHNOLOGIST**

Two full-time openings for Registered or Registry Eligible Nuclear Medicine Technologists. Both positions are in the Nuclear Medicine Department to do large varieties of studies. Baptist Medical Center is a 565 bed modern hospital located in the pleasant suburbs of Oklahoma City. Excellent salary and employee benefits await qualified applicants.

For information call (collect) Gene Stanford at (405) 949-3202 or send resume to:

Baptist Medical Center



Employment Office
3300 NW Expressway
Oklahoma City, OK 73112

**NUCLEAR PHYSICIAN
CERTIFIED ABNM OR ABR**

Faculty position available July 1981 in well equipped diagnostic imaging division with educational and clinical Nuclear Medicine responsibilities. Excellent research opportunities and facilities are available. Some Ultrasound and computer experience desired. Medical University complex is located in historic Charleston on the Atlantic Coast bordered by three rivers and barrier islands. Contact: **K.M. Spicer, M.D., Ph.D.,**
Director, Division of Diagnostic Imaging
Medical University of South Carolina
171 Ashley Avenue
Charleston, SC 29403

**FELLOWSHIP PROGRAMS IN
DIAGNOSTIC IMAGING**

One (AB Radiology) or two (AB Nuclear Medicine) year approved programs offer pediatric and adult experience in all aspects of imaging including Ultrasound, CT, and Cardiovascular Nuclear Medicine, and training in computers, Radioassay and basic sciences. Facilities include state of the art instrumentation, centralized Nuclear Pharmacy and well equipped In-Vitro and animal laboratories for excellent research opportunities.
Available July 1, 1981
Contact: K.M. Spicer, M.D., Ph.D.,
Director, Division of Diagnostic Imaging,
Medical University of South Carolina
171 Ashley Avenue, Charleston, SC 29403

INSTRUCTOR

State University of New York at Upstate Medical Center is seeking candidates for a position as Instructor of Nuclear Medicine Technology. This is a full time position which carries a faculty appointment in the College of Health Related Professions.

Qualifications include graduation from AMA approved Nuclear Medicine Technology Program, certification by AART or NMTCB and some teaching experience. Baccalaureate Degree desirable, but not required. Some computer experience useful. Excellent starting salary plus an outstanding benefit program.

To apply please forward resume to:
Charles E. Duxbury, R.T.
Assistant Administrator
Division of Nuclear Medicine
750 East Adams Street
Syracuse, New York 13210

RESIDENCY IN NUCLEAR MEDICINE

Two-year approved program offering broad clinical experience including tertiary care and community hospitals; oncology and pediatrics; ultrasound and CT; strong basic science teaching; radiation safety; central radiopharmacy and RIA; opportunity for research; an integrated program at State University of New York at Buffalo School of Medicine; available July 1, 1981.

Contact: **M.A. Bender, M.D.,** Program Director, Dept. of Nuclear Medicine, Roswell Park Memorial Institute, 666 Elm Street, Buffalo, NY 14263; or **M. Blau, Ph.D.,** Chairman, Dept. of Nuclear Medicine, SUNY/Buffalo, 3495 Bailey Avenue, Buffalo, NY 14215.

**CHIEF
Nuclear Medicine Technologist**



JOHNSTON-WILLIS HOSPITAL
2908 Kensington Avenue
RICHMOND, VIRGINIA 23221

Johnston-Willis Hospital located in historic Richmond, Virginia is seeking a Chief Nuclear Medicine Technologist to work in our new 292-bed acute care facility.

Required background includes experience in cardiac procedures, In-vivo procedures, quality control program, Shillings Testing, Blood Volume testing and supervision.

Equipment includes Siemens (Searle) and Ohio Nuclear LFOV Cameras and Siemens Scientiview and MDS Computer Systems.

Salary is commensurate with experience and excellent company benefits are offered.

If interested, please contact **Bonnie Shoaf,** Employment Coordinator, Johnston-Willis Hospital, 1401 Johnston-Willis Drive, Richmond, Virginia 23235.

An Affiliate of Hospital Corporation of America

NUCLEAR MEDICINE TECHNOLOGIST

Full-time day position available for a graduate of an A.M.A. approved school of radiologic technology who is registered with the AART. Graduate of 2-year approved school in nuclear medicine technology approved by the American Society of Radiologic Technologists. Registered by the American Registry of Nuclear Technologists. Minimum of one year experience as registered staff nuclear medicine technologists in performing radionuclide procedures on patients including organ imaging, function studies, dynamic studies, and computer data storage and techniques requiring computer application in nuclear medicine.

For an opportunity to join an expanding department and to receive the excellent benefits and salary which accompanies this position, please send resumes to:

Raleigh Clemons III
Sr. Personnel Asst.
Kaiser Permanente Medical Program
2475 E. Blvd.
Cleveland, Ohio 44120

V.A. Medical Center, Long Beach, CA
 Affiliated with
University of California at Irvine

NUCLEAR MEDICINE RESIDENCIES

Position available July 1981 for first year resident in A.M.A. approved program. Professional staff includes radiopharmacist and physicist offering broad opportunity for clinical and research experience. Equal opportunity employer. English language proficiency required (PL 95-201). Contact:

Kenneth P. Lyons, M.D.
Chief, Nuclear Medicine Services
 or
Norah Milne, M.D.
Director of Training
V.A. Medical Center
Long Beach, CA 90822
(213) 498-6237

JNM CLASSIFIED PLACEMENT SERVICE SECTION


This section in the *Journal of Nuclear Medicine* contains "Positions Open", "Positions Wanted", and "For Sale" listing. Nondisplay "Positions Wanted" ads by members of the Society are billed at 60¢ per word for each insertion with no minimum rate. Nondisplay "Positions Wanted" ads by nonmembers and all nondisplay "Positions Open" and "For Sale" ads by members and nonmembers are charged at 85¢ per word. Display advertisements are accepted at \$125 for 1/8 page, \$185 for 1/4 page, \$295 for 1/2 page, and \$510 for a full page. Closing date for each issue is the 1st of the month preceding publication. Agency commissions and cash discounts are allowed on display ads only. Box numbers are available for those who wish them.

All classified ads must be prepaid or accompanied by a purchase order. Send orders to:

Journal of Nuclear Medicine
475 Park Ave. South
New York, NY 10016

You saw it in Detroit

AT THE SNM ANNUAL MEETING



Pharmatopes, Inc.

Mail to: T-SHIRT
 PHARMATOPES, INC
 25721 COOLIDGE
 OAK PARK, MI 48237


Sizes & Quantities S M L XL

Name

Shipping Address City State Zip

Representing Hospital City State Zip

Each order must be accompanied by a check or money order payable to Pharmatopes, Inc. Price is \$6.00 per shirt. .50 postage per shirt.



**The
Superior Nuclide
For
Thyroid Studies**

**Sodium
Iodide I 123**

I 123... A Superior Thyroid Agent

Sodium Iodide I 123 is superior to I 131 because of its low radiation dose to the patient, its short half-life of 13.2 hours, and its imaging energy of 159 KeV.

Sodium Iodide I 123 is superior to Tc99m because it is trapped and organified by the thyroid gland and, therefore, will image the "cold," non-functioning nodule that may appear "hot" or "cold" with Tc99m.^{1,2}

For a Consistent Quality Image.....Sodium Iodide I 123

medi+physics™

5801 Christie Avenue, Emeryville, CA 94608

For More Information, Please Call (415) 658-2184

Inside California Toll Free (800) 772-2446 • Outside California Toll Free (800) 227-0483

¹Steinbach, HL, Kundy, D, Moss M, et al: A comparison of three agents in thyroid uptake and scintigraphy. Scientific Exhibit, Society of Nuclear Medicine, Philadelphia, June 16-20, 1975.

²"Information for Physicians—Irradiation-Related Thyroid Cancer" prepared by the Division of Cancer Control and Rehabilitation National Cancer Institute, DHEW Publication No. (NIH) 77-1120, p 13

For complete prescribing information consult package insert, a summary of which follows:

SODIUM IODIDE I 123 CAPSULES AND SOLUTION FOR ORAL ADMINISTRATION

DESCRIPTION: Sodium iodide I 123 for diagnostic use is supplied as capsules and in vials as an aqueous solution for oral administration. At calibration time, each capsule has an activity of 100 microcuries and each vial contains solution with a total specific concentration of two millicuries per ml.

INDICATIONS: Sodium iodide I 123 is indicated for use in the diagnosis of thyroid function and imaging.

CONTRAINDICATIONS: None known.

WARNINGS: This radiopharmaceutical should not be administered to children or to patients who are pregnant or to nursing mothers unless the information to be gained outweighs the potential hazards. Ideally, examinations using radiopharmaceuticals, especially those elective in nature, in women of child-bearing capability should be performed during the first few (approximately 10) days following the onset of menses. However, when studies of thyroid function are clinically indicated for members of these special population groups, use of I 123 would be preferable to the use of I 131 in order to minimize radiation dosage.

Adequate reproduction studies have not been performed in animals to determine whether this drug affects fertility in males or females, has teratogenic potential, or has other adverse effects on the fetus. Sodium iodide I 123 should be used in pregnant women only when clearly needed.

PRECAUTIONS: Sodium iodide I 123, as well as other radioactive drugs, must be handled with care. Appropriate safety measures should be used to minimize radiation exposure to clinical personnel. Care should also be taken to minimize radiation exposure to the patient consistent with proper patient management. The prescribed sodium iodide I 123 dose should be admin-

istered as soon as practicable in order to minimize the fraction of radiation exposure due to relative increase of radionuclidic contaminants with time. The uptake of I 123 may be decreased by recent administration of iodinated contrast materials, by intake of stable iodine in any form, or by thyroid, anti-thyroid, and certain other drugs. Accordingly, the patient should be questioned carefully regarding diet, previous medication, and procedures involving radiographic contrast media.

ADVERSE REACTIONS: There were nine adverse reactions reported in a series of 1,393 administrations. None of these were attributed to I 123. Five adverse reactions, consisting of gastric upset and vomiting, were attributed to a filler in the capsule. Two cases of headache and one case of nausea and weakness were attributed to the fasting state. One case of garlic odor on the breath was presumed to be attributable to the presence of tellurium.

DOSAGE AND ADMINISTRATION: The recommended oral dose range for diagnostic studies of thyroid function in the average adult patient (70 kg) is from 100 to 400 microcuries. The patient dose should be measured by a suitable radioactivity calibration system immediately prior to administration. Concentration of I 123 in the thyroid gland should be measured in accordance with standardized procedures.

SPECIAL CONSIDERATION: Radiopharmaceuticals should be used only by physicians who are qualified by training and experience in the safe use and handling of radionuclides and whose experience and training have been approved by the appropriate government agency authorized to license the use of radionuclides.

HOW SUPPLIED: Sodium iodide I 123 for oral administration is supplied in aqueous solution in glass vials of 1mCi and in capsules of 100 µCi.



REGISTERED NUCLEAR MEDICINE TECHNOLOGIST

MEDICAL CENTER HOSPITALS offers a great multitude of settings for your talent as a NUCLEAR MEDICINE TECHNOLOGIST. Excellent salary and benefits with features such as:

**PAID ANNUAL LEAVE
PAID LIFE INSURANCE
TUITION REIMBURSEMENT**

As key professionals in a sophisticated medical environment, Nuclear Medicine Technologists are an integral part of our spirited health care team. We offer the opportunity for personal and professional growth in the largest non-profit civilian hospital in Virginia. Work with one of the most progressive Radiology Departments in the country...Contact Personnel -

PAULETT MASON

**MEDICAL CENTER HOSPITAL
NORFOLK GENERAL DIVISION**

600 Gresham Drive
Norfolk, Virginia 23507

(804) 628-3831

An Equal Opportunity Employer

Technical Service Specialist

As the world leader in the field of nuclear medicine, New England Nuclear is continually expanding to meet the steady increase in customer demand. We currently have open a new position for a Technical Service Specialist in our Radiopharmaceutical Development Department.

The primary responsibility for this position will be to respond to a wide variety of technical product inquiries. These inquiries are related to In-Vivo or In-Vitro use of Radiopharmaceutical products. Additional responsibilities include assisting NEN customers in areas of nuclear instrumentation, licensing, and radiopharmaceutical quality control. This individual will assume a key responsibility in keeping the NEN management team aware of customer and product related information pertinent to product acceptability.

The ideal candidate for this position will have a BS degree or CNMT Certification in Nuclear Medicine Technology plus 3 to 5 years as a Nuclear Medicine Technologist. Experience in Nuclear Cardiology is desirable. Strong communication skills are essential.

Interested candidates should submit a resume with salary history to Susan Kinney, New England Nuclear, 601 Treble Cove Road, Billerica, MA 01862.



New England Nuclear

...A BETTER WORLD

an equal opportunity employer

NUCLEAR CARDIOLOGY TECHNOLOGIST

Expanding nuclear medicine dept. with new MDS computer and portable and stationary cameras seeks person with excellent background and experience in nuclear cardiology/computers to begin new service. Planned computer research project. Very progressive 238-bed general acute non-profit community hospital with 90% occupancy. Culturally rich sea coast community. Requires nuclear medicine registration and demonstrated leadership skills. Salary range \$1,500 to \$1,900 per month. Contact:

**Warren Rowe, Director of Personnel
SAN PEDRO PENINSULA HOSPITAL
1300 West 7th Street
San Pedro, California 90732
(213) 832-3311 Ext. 3131**

NUCLEAR MEDICINE TECHNOLOGIST

Registered or registry eligible technologist for full time position in modern 410 bed acute care hospital. St. Mary's is located in a city of 100,000 midway between St. Louis and Chicago. Contact:

**Personnel Department
St. Mary's Hospital
1800 E. Lake Shore Drive
Decatur, IL
(217) 429-2966**

JUST PUBLISHED!

Single Photon Emission Computed Tomography and Other Selected Computer Topics


The 21 papers in this volume include an overview of single photon emission computed tomography, and numerous papers that describe and evaluate specific systems and techniques. Angled collimators, seven-pin-hole collimators, and Fourier tomography are included in the methods and systems described. Performance specifications of various instruments are discussed and compared. The volume also provides a review of the various algorithms and computer systems available. Papers on the use of computers in the analysis of cardiac studies, data acquisition and storage methods, Anger camera evaluation, and whole body counting round out the book.

244 pages. Illustrated. Proceedings of a Society of Nuclear Medicine Computer Council Symposium, January 24, 1980, Miami Beach, FL. Ronald R. Price, Program Chairman.

Order from: Book Order Dept., Society of Nuclear Medicine, 475 Park Ave. South, New York, NY 10016. List price—\$27.00. SNM members in good standing may purchase the book for \$18.00. \$2.50 per book for postage and handling must be included with all orders. Purchase order or prepayment in U.S. funds drawn on a U.S. bank is required.

ISBN 0-932004-06-7

EXCEPTIONAL



Left: Erect stress test in progress. Note that pedal unit has been lowered to eliminate contact with patient's knees when camera head is positioned.

Above: Supine stress test in progress. Far side of table is unobstructed to easily accommodate a gamma camera. Adjustment wheels control vertical and horizontal travel of pedal unit for patients of all sizes.

Cardiac Stress Table & Ergometer System

**Permits all patient positions, from supine through upright,
plus standard imaging procedures**

VERSATILE • PRACTICAL • COST EFFECTIVE

Here is the most versatile stress imaging table available. It permits routine stress tests of patients in *any* position, from supine to upright. Whatever your nuclear cardiology requirements are, this exceptional system provides them quickly and easily...full gamma camera clearance, complete mobility, vibration-free rigidity, positive (but comfortable) patient restraint, easy-to-make positioning adjustments, unobstructed access to the patient and controls, choice of Collins or Quinton ergometers, and much more. The unit also doubles as a standard imaging table.

Send for full details. Ask for Bulletin 2890-B

**VICTOREEN
NUCLEAR ASSOCIATES**




100 Voice Road
Carle Place, N.Y. 11514
(516) 741-6360



Patient positioned for standard imaging. Legs and feet are supported comfortably.





Power. To push around.

Mobile Power You Can Control.

Take it anywhere—and interface it with any gamma camera.

Newest of our Spectrum One nuclear medicine computers and cameras, the MCS-560 makes analysis of cardiology and nuclear medicine procedures easy.

Easy To Use.

The MCS-560 is the most powerful mobile nuclear medicine computer anywhere. Sophisticated analysis is easy with: conversational programs pushbutton protocols unique MEDI-BASIC programming language built-in ECG Isolator/Detector unique tomographic reconstruction programs.

Backed by our own dedicated nuclear products service team, we're building our one-source reputation for excellence.

TECHNICARE


Technicare Corporation
29100 Aurora Road
Solon, Ohio 44139
(216) 248-1800

A **Johnson & Johnson** Company

Formerly Ohio-Nuclear, Inc.

Minitec[®]

(Technetium Tc 99m)
Generator



Hot where it matters Cold where it counts



Designed for minimum exposure

- Unique construction (no exposed tubing) and thick shielding (1 3/8" lead) provide high shielding-to-activity ratio.
- Another 1 1/2" of lead shielding provided by Maxi-Shield™.
- Built-in handle reduces hand exposure during carrying.
- A comparison study of radiation exposure from the three leading generator systems is available from your Squibb Representative.



Convenient

- Available in potencies of 220, 440, 880, 1330, 1770 and 2220 mCi.
- Tuesday-calibrated generators delivered Wednesday a.m.
- Saturday-calibrated generators delivered Monday a.m.

Easy to operate

- After a few simple preparations, Minitec elutes automatically and quickly. Subsequent elutions are even simpler.
- Small-volume, high-concentration eluates give maximum flexibility for varying applications.



On-time delivery

- Express ground transportation and special air system assure on-time deliveries.

Dependable service

- Specially trained Technical Associates and Squibb Technical Customer Service provide prompt personal attention when needed.

See next page for brief summary.

Minitec[®] (Technetium Tc 99m) Generator

Medotopes[®]



SQUIBB[®]

Minitec[®]

(Technetium Tc 99m)
Generator



Hot where it matters
Cold where it counts

MINITEC[®]
Technetium Tc 99m
GENERATOR

DESCRIPTION: Minitec (Technetium Tc 99m) Generator consists of a specially designed lead-shielded alumina column containing adsorbed fission-produced Mo 99. Tc 99m, the short-lived daughter of Mo 99, is obtained as sterile sodium pertechnetate Tc 99m by periodic elutions of the generator with an isotonic saline solution.

INDICATIONS AND USAGE: Sodium pertechnetate Tc 99m is indicated in ADULTS as an agent for brain imaging including cerebral radionuclide angiography, thyroid imaging, salivary gland imaging, placenta localization, and blood pool imaging including radionuclide angiography. (For use of sodium pertechnetate Tc 99m as a diagnostic radiopharmaceutical in CHILDREN, consult package insert.)

CONTRAINDICATIONS: None known.

WARNINGS: Radiopharmaceuticals should not be administered to patients who are pregnant or to nursing mothers unless the expected benefit to be gained outweighs the potential hazards.

Since sodium pertechnetate Tc 99m is excreted in human milk during lactation, formula-feedings should be substituted for breast-feedings.

Ideally, examinations using radiopharmaceuticals, especially those elective in nature, of a woman of childbearing capability should be performed during the first few (approximately 10) days following the onset of menses.

Radiation risks associated with the use of sodium pertechnetate Tc 99m are greater in children than in adults and, in general, the younger the child the greater the risk owing to greater absorbed radiation doses and longer life expectancy. These greater risks should be taken firmly into account in all benefit-risk assessments involving children.

IMPORTANT: Since material obtained from the generator may be intended for intravenous administration, aseptic technique must be strictly observed in all handling. Only the eluent provided should be used to elute the generator. Do not administer material eluted from the generator if there is any evidence of foreign matter.

PRECAUTIONS: In the use of any radioactive material, care should be taken to insure minimum radiation exposure to the patient and occupational workers consistent with proper patient management. At the time of administration the solution should be crystal clear.

ADVERSE REACTIONS: No adverse reactions specifically attributable to the use of sodium pertechnetate Tc 99m have been reported.

For full prescribing information, consult package insert.

HOW SUPPLIED: Minitec (Technetium Tc 99m) Generator is available in potencies of 220, 440, 880, 1330, 1770, or 2220 millicuries Mo 99 at calibration time. The generator is supplied with vials of sterile, nonpyrogenic eluent; a sterile needle adapter assembly and evacuated sterile collecting vials. Other accessories including lead shields, reference standard solutions, and a whole vial assay kit are available on request for use with the Minitec (Technetium Tc 99m) Generator.

SQUIBB[®]

MIRD PAMPHLETS AVAILABLE (Medical Internal Radiation Dose)

PAMPHLETS

- 1 (Revised) A revised schema for calculating the absorbed dose from biologically distributed radionuclides. (\$5.25)
- 5 (Revised) Estimates of specific absorbed fractions for photon sources uniformly distributed in various organs of a heterogeneous phantom. (\$7.75)
- 10 Radionuclide decay schemes and nuclear parameters for use in radiation-dose estimation. (\$8.00)
- 11 'S' absorbed dose per unit cumulated activity for selected radionuclides and organs. (\$11.00)
- 12 Kinetic models for absorbed dose calculations. (\$5.25)

SUPPLEMENTS

- 1 Includes 3 pamphlets: "Schema for absorbed dose calculations for biologically distributed radionuclides"; "Energy deposition in water by photons from point isotropic sources"; and "Absorbed fractions for photon dosimetry." (\$1.50)
- 3 Includes the original pamphlet #5: "Estimates of absorbed fractions for monoenergetic photon sources uniformly distributed in various organs of a heterogeneous phantom." (\$1.50)
- 5 Includes 2 pamphlets: "Distribution of absorbed dose around point sources of electrons and beta particles in water and other media"; and "Absorbed fractions for small volumes containing photon-emitting radioactivity." (\$1.50)
- 6 Includes pamphlet 9: "Radiation dose to humans from ⁷⁵Se-L-Selenomethionine." (\$3.00)

SPECIAL OFFER

All available MIRD pamphlets and supplements for only \$25.00 plus \$4.00 for shipping and handling.

Attractive binders for the pamphlets and supplement #1 are available at \$4.50 each.

MIRD Pamphlets and Supplements may be ordered from: Book Order Dept., Society of Nuclear Medicine, 475 Park Avenue South, New York, NY 10016. All orders must be prepaid or accompanied by a purchase order. Checks must be in U.S. funds only, please.

Mail to: Book Order Dept., Society of Nuclear Medicine, 475 Park Avenue South, New York, NY 10016. Make checks payable to: Society of Nuclear Medicine, Inc., U.S. funds only, please.

PAMPHLETS	SUPPLEMENTS	SPECIAL OFFER
___ 1 (\$5.25)	___ 1 (\$1.50)	___ \$25.00 plus
___ 5 (\$7.75)	___ 3 (\$1.50)	___ \$4.00 for shipping
___ 10 (\$8.00)	___ 5 (\$1.50)	___ & handling. (Does
___ 11 (\$11.00)	___ 6 (\$3.00)	___ not include binder.)
___ 12 (\$5.25)		

BINDERS ___ \$4.50 each

SHIPPING AND HANDLING CHARGES

1 item	\$1.00	10-19 items	\$6.00
2 items	2.00	20-29 items	8.00
3 items	3.00	30-39 items	10.00
4-9 items	4.00		

TOTAL \$ _____

SHIPPING & HANDLING CHARGES \$ _____

FOREIGN ORDERS ADD \$7.50 \$ _____

TOTAL ENCLOSED \$ _____

SEND TO:

NAME _____

ADDRESS _____

ZIP _____

JNM 10/80



Introducing a better image.

**Toshiba's GCA-50A
Low Energy Mobile Gammacamera
with 1.8mm resolution.**

TOSHIBA
MEDICAL SYSTEMS

The 37 photomultiplier tubes and superior electronics in Toshiba's all new GCA-50A provide an exceptional intrinsic detector resolution of 1.8mm to yield a more functional diagnostic image.

In addition, the large 12.25" crystal diameter and high count rate of more than 200 kcps (40% window width) reduce scanning times and optimize the effective use of your unit.

Reliability is also exceptional. The GCA-50A design is based on the internationally proven performance of Toshiba's GCA-402. But now you have the flexibility of a totally mobile unit along with complete options like a cardiac synchronizer, 4- or 16-frame multi-image formatting, and computer capability.

So if you want to speed scanning time while improving your images and the overall effectiveness of the nuclear medicine group, contact your Toshiba representative now. Ask for full information on the GCA-50A Mobile Gammacamera.



Toshiba Medical Systems
1154 Dominguez
Carson, CA 90745

A division of Toshiba International Corporation

Help your cardiologist study heart kinetics non-invasively with Brattle-gated scintiphotos.



RAO, DIASTOLE



RAO, SYSTOLE



LAO, DIASTOLE




LAO, SYSTOLE

The RAO view shows akinesis of the lower antero-lateral wall and apex; and contraction of the inferior wall and high up the antero-lateral wall. The LAO view shows good contrac-

tion posteriorly and akinesis of the septal aspect of the chamber. Patient was injected IV with 20mCi of ^{99m}Tc -labelled Human Serum Albumin. The agent was prepared using the New

England Nuclear Electrolysis Kit for labelling HSA. Write or call for a portfolio of Brattle-gated lung, liver and heart studies.



No knobs, no meters, no errors

The spartan panel above tells the second-best part of our story. If you want to photograph peak systole, press the SYSTOLE button. If, say, you want systole only at full expiration, press the EXPIRATION button as well. If only breathing is relevant, don't press the heart button.

The Brattle is connected to the patient and to your gamma (or x-ray or ultrasonic) camera. Whenever the patient is in the selected phase, both the scope and the scaler on your gamma camera are gated ON, and film is exposed. Otherwise, they are OFF.

Brattles lock onto patients — and stay locked on

It doesn't matter if the patient's heart rate and breathing depth change while he's under the collimator be-

cause we stay right with him. Brattles contain an ECG to track heart, a plethysmograph to track respiration, and a tiny computer to deduce systole and diastole times from the heart signal. And because it's all built in, your operator need not be a physiologist.

We don't cover our tracks — we print them

The panel lights flash whenever the patient reaches the selected phases; and pushing the RECORDER-ON button gets you an ECG tracing marked with breathing and camera-on times. You can verify function before, during and after exposure.

A single pair of axillary electrodes captures both heart and breath
It's easy. And we supply disposable, pre-filled electrodes.

Some Brattles have been in clinical use for over three years — in community and major hospitals

More than half of our instruments are in community hospitals and the list is growing rapidly. Upon request, we'll supply names of happy users in your area.

What's the next step? Get in touch

Ask your NEN man about Brattles and HSA Kits. He can show you a portfolio of clinical pictures and arrange to have one of our people give you a demo. Or write or call us direct. We'll send you brochures on this and other models, and will give you your own set of clinical pictures and a bibliography on gated scintigraphy. If you wish, we'll even make you a Brattle owner. (This is the best part of our story.)

Brattle Instrument Corporation

243 Vassar Street • Cambridge, Massachusetts 02139 • 617-661-0300

The CentRIA® System Tablet

**not solid phase...
not liquid phase...**

Separation so rapid...
so specific...so complete...
it's revolutionary!

If you need cost-effective,
simple RIA, with no reagent
"lock-in," and flexibility to
meet your laboratory needs,
call: **800-431-1944** toll free.

In New York State call
collect: **914-769-9452**



270 Marble Avenue
Pleasantville, N.Y. 10570