

THE ASSAY OF IODINE-123

F. R. Hudson, H. I. Glass, and S. L. Waters
Hammersmith Hospital, London, England

Compounds labeled with ^{123}I are becoming more common in scintillation imaging due to the convenient combination of the 13.3-hr half-life and the 159-keV gamma radiation. Here attention is drawn to the particularly large summing effects observed when ^{123}I is counted in a well counter. These could lead to errors in the assay of doses if the analyzer windows used are too narrow.

The properties of ^{123}I make it very suitable for use in scintillation camera studies and other nuclear medicine investigations. Its 13.3-hr half-life ensures considerably reduced patient exposure to radiation compared with ^{131}I , and the principal gamma-ray at 159 keV can be detected with high efficiency. Only limited use has been made of ^{123}I until recently, due to its restricted availability. However, an increasing interest is now being shown in its use (1-4) and it seems desirable to comment on the strong sum peaking that may arise when samples containing ^{123}I are counted in a well counter.

MATERIALS AND METHODS

The ^{123}I used in this laboratory is prepared in the Medical Research Council cyclotron by the method of Silvester et al (5). The principal contaminants found in the preparation are ^{124}I (half-life 4.2 days), ^{125}I (60 days), and ^{126}I (13 days). The typical concentrations of these isotopes relative to ^{123}I are 5%, 0.9%, and 0.5%, respectively, at approximately 5 hr after the end of bombardment. These levels rise significantly if there is any substantial delay in the use of the preparation.

Assay of doses for injection is frequently carried out with a NaI well counter and the intercomparison of doses and blood samples in the course of an investigation may involve the counting of samples at different geometries to avoid excessive deadtime

losses. For a 0.5- μCi dose of ^{123}I dispensed in 2 ml of saline into a polythene vial for counting in a standard well counter (Fig. 1), very different spectra were observed when the vial is supported at, for example, the level of the top face of the well (Fig. 2) and when the vial is fully inserted into the well (Fig. 3).

The spectrum of Fig. 2 shows the peaks of the 159-keV ^{123}I gamma-ray, the 28-keV Te x-ray, and 74-keV Pb x-ray. Figure 3 shows additional peaks at 61 keV and 187 keV. (The Pb x-ray is no longer visible because the sample is almost totally surrounded by the NaI crystal.) The additional peaks have been identified by Hupf, Eldridge, and Beaver (6) as arising from x-ray/x-ray coincidences at 61 keV and from x-ray/gamma-ray coincidences at

Received May 21, 1975; revision accepted Sept. 25, 1975.
For reprints contact: S. L. Waters, Medical Physics Dept., M.R.C. Cyclotron Unit, Hammersmith Hospital, London W12 0HS, U.K.




FIG. 1. Sodium iodide well counter with sample vial fully inserted.




FIG. 2. Energy spectrum for ^{123}I sample with vial base level with top of crystal face.




FIG. 3. Energy spectrum for ^{123}I sample with vial fully inserted.

187 keV. The x-ray/x-ray coincidences arise from summing the x-ray that follows the ^{123}I electron capture and the x-ray that follows the internal conversion of the 159-keV gamma-ray of ^{123}Te . The summing effects are unusually strong in ^{123}I and a calibration has been carried out to show the magnitude of the errors that might arise if sources are counted under different geometric conditions with narrow analyzer window settings.

A sample of ^{123}I was dispensed into a standard polythene counting vial in 2 ml of saline. As the sample was withdrawn from the well, counts were taken in a 145–175-keV window and in a 145–205-keV window. The ratio of the two counts is plotted against the distance of the center of the sample from the base of the well in Fig. 4.

RESULTS AND DISCUSSION

Examination of Fig. 4 shows that errors up to 50% may occur if samples are counted under dif-

ferent geometric conditions with narrow analyzer window settings. This could happen if activities are assayed in different volumes; absorption of the 28-keV x-ray within a sample or sample container could also affect the strength of the summing effects. Remember that the entrance of the well counter is different for samples placed in the well and for samples supported above the well. This difference may be significant for low-energy counting.

The way in which the detector-sample geometry affects this formation of sum peaks is illustrated in Fig. 5. This shows why the intensity of the 187-keV sum peak falls steadily when a sample is withdrawn from the well and why it is difficult to define the sum peak unambiguously if the sample is counted at the face of a standard $7.5 \times 7.5\text{-cm}$ NaI crystal.

Examination of the ^{123}I decay scheme in the *Nuclear Data Tables* (7) shows that ^{123}I decays to ^{123}Te by electron capture. The 28-keV Te x-ray is thus emitted in coincidence with the 159-keV gamma-ray.




FIG. 4. Variation of K with D, where K is ratio of counts in 145–175-keV window to counts in 145–205-keV window and D is height of sample's center over well base.




FIG. 5. Geometry dependence of summing effects illustrated by decay event (left) deep in counter well which produces summing and by decay event (right) near counter surface which does not produce sum peak.

The summing of these photons is particularly strong as the gamma-ray and the x-ray occur in 82.9% and in 86.5% of decays, respectively. The sum peak at 187 keV is clearly separated from the primary gamma-ray peak. This separation occurs only in nuclides of mass number greater than approximately 100. The sum peak for ^{123}I is also easy to observe since it lies in a clean region of the spectrum.

Summing of this magnitude is not common among nuclides in regular use in nuclear medicine as it is unusual for the x-ray intensity to be so high. However, ^{203}Pb , proposed by Stark et al (8) as a potentially useful agent for cisternography, shows similar effects. Clearly, therefore, when a comparison is made between a standard and an unknown containing a summing nuclide, they should be counted in the same geometry whatever window system might have been chosen. If the count rate of a specimen produces deadtime problems, then the specimen should be diluted and not just placed at a less efficient geometry farther away from the crystal.

ACKNOWLEDGMENTS

Thanks are due to M. Thakur and D. Williams (M.R.C.

Cyclotron Unit) for the preparation of the ^{125}I and to S. Bowring and J. Nosil for helpful comments.

REFERENCES

1. SHORT MD, GLASS HI, CHISHOLM GD, et al: Gamma-camera renography using ^{125}I -hippuran. *Br J Radiol* 46: 289-294, 1973
2. ATKINS HL, KLOPPER JF, LAMBRECHT RM, et al: A comparison of technetium 99m and iodine 123 for thyroid imaging. *Am J Roentgenol Radium Ther Nucl Med* 117: 195-201, 1973
3. MCKEIGHAN RE, MUEHLEHNER G, MOYER RA: Gamma camera collimator considerations for imaging ^{125}I . *J Nucl Med* 15: 328-331, 1974
4. NISHIYAMA H, SODD VJ, BERKE RA, et al: Evaluation of clinical value of ^{125}I and ^{123}I in thyroid disease. *J Nucl Med* 15: 261-265, 1974
5. SILVESTER DJ, SUGDEN J, WATSON IA: Preparation of ^{125}I by α -particle bombardment of natural antimony. *Radiochem Radioanal Lett* 2: 17-20, 1969
6. HUPF HB, ELDRIDGE JS, BEAVER JE: Production of iodine-123 for medical applications. *Int J Appl Radiat Isot* 19: 345-351, 1968
7. MARTIN MJ, BLICHERT-TOFT PH: *Nuclear Data Tables*, A8, 1-2 Oct 1970, New York, Academic, p 88
8. STARK VJ, HARPER PV, LATHROP KA, et al: Evaluation of ^{203}Pb as an imaging agent for cisternography. *J Nucl Med* 13: 468-469, 1972

23rd ANNUAL MEETING

THE SOCIETY OF NUCLEAR MEDICINE

June 8-11, 1976

Dallas Convention Center

Dallas, Texas

Members and nonmembers of the Society of Nuclear Medicine are invited to attend the 23rd Annual Meeting. The convenience and comfort of attendees will be emphasized as much as possible, with good food service, snack bars, and lounge areas readily available in the Convention Center. This year's program will include lectures by distinguished scientists in specialty areas as well as selected research papers of merit. Due to the popularity of the teaching sessions, the format has been expanded so that teaching is not confined to the early morning hour, but is distributed throughout the daily program.

A wide variety of nuclear medicine topics in clinical research, clinical practice, and basic science will be covered, including the following: Bone/Joint, Cardiovascular, Computer/Data Analysis, Computerized Axial Tomography, Dosimetry, Endocrine/Metabolism, Gastroenterology, Hematology, Instrumentation, In Vitro Assays, Neurology, Oncology, Pediatrics, Pulmonary, Radiopharmaceuticals, and Renal/Electrolytes. Scientific sessions for the Technologist Section will be held concurrently.

Registration forms are available from the Society of Nuclear Medicine, 475 Park Avenue South, New York, N.Y. 10016.