

IN A FOG??

using aerosols to determine the patency of the pulmonary airway system? Use a gas (that's what the airway system is for), and Xenon (127 or 133) are gases which are safe, economical and easy to administer with the XENAMATIC™ 3000.

- Shielded for Xe 127 and Xe 133 (radiation profile available on request).
- World's only system that allows you to study patients on Ventilators.
- Largest and most efficient Xenon trap with a built-in monitor alarm system.
- Built-in O₂ monitor with digital display and control.
- A rebreathing system that saves Xenon.
- Low breathing resistance so you can study sick patients.
- Semi-automatic operation.
- Remote Control Capability.

Get out of the FOG-making business, and call today for more information on putting gases where gases belong, with the XENAMATIC.

Also available, Model 2000.

For more information, please call or write,

Circle Reader Service No. 32




DIVERSIFIED DIAGNOSTIC PRODUCTS, INC.


11603 Windfern
Houston, TX 77064
713-955-5323

Announcing A Big *STEP*[™] Forward In Nuclear Cardiac Imaging.

Conventional thallium image shows inferior diaphragmatic attenuation artifact.



STEP image eliminates artifact, clearly showing normal perfusion.



Picker proudly introduces STEP[™]—Simultaneous Transmission Emission Protocol. The first and only non-uniform attenuation correction device for 360° cardiac SPECT.

This new technology, combined with your PRISM[™] triple-head

system, eliminates artifacts due to bone, breast, diaphragmatic, and other soft-tissue attenuation that may be mistaken for coronary

artery disease. And since STEP does not add to acquisition time, patient throughput remains high.

So take a big STEP into the future of nuclear

cardiac imaging. Call Picker today at 1-800-323-0550 for more information on STEP.

Picker International, 595 Miner Road, Cleveland, Ohio 44143.

 **PICKER**
MORE THAN IMAGES. INSIGHT.

DOES YOUR WIPE TEST COUNTER MEASURE UP?




The CAPRAC® Wipe Test Counter from Capintec measures up from the outset because it gets right down to basics - it saves time and improves your results.

The NaI drilled well crystal lets you perform your wipe tests using short counting times while detecting very low levels of activity (seconds for 1 nCi). All this with the accuracy to meet NRC and state regulations.

Another very important measure is "what's on the wipe". That's no prob-

CAPRAC SN:000123
May 16, 1991 15:53

W-01



	875.	± 42.	cpm
CF: 2.50	2.19	± .10	kdpm
I 131	2.42	± .14	kdpm
TC99m	222.	± 24.	dpm
GA67	1375.	± 61.	dpm
IN111	1440.	± 64.	dpm
TL201	439.	± 32.	dpm

lem with CAPRAC because of its energy discrimination and isotope identification through gamma spectroscopy, specific isotopes can be identified at the touch of a button.


The final measure? Permanence! CAPRAC excels here, too, with its optional printer that allows permanent printed records of all data and gamma energy histograms.

CAPRAC has many more features you should know about. Please contact us for complete details.



CAPINTEC, INC.

6 ARROW ROAD, RAMSEY, N.J. USA 07446
TOLL FREE: (800) 631-3826 OR
(201) 825-9500 FAX: (201) 825-1336






*Helix*TM
The Revolution Never Stops

Available with
All Digital
Detectors

Elscint

Slip-Ring Nuclear Tomography

Continuous high speed orbiting - possible only with **Slip-Ring** technology - revolutionizes SPECT imaging. For the first time, nuclear tomographic slices evolve in real-time, while you watch. Increased throughput and optimized image quality are immediate benefits. And Helix expands horizons to innovations like Dynamic Tomography and Whole-Body SPECT, opening the door to future applications.



Real-Time Automatic Body Contouring

Top clinical performance is matched by outstanding ergonomic design. Totally automated **real-time** body-contouring eliminates pre-scan set-up. And effortless, simultaneous dual-collimator exchange makes Helix the technologist's best friend.

All-Digital Camera Design

Whole-Body and SPECT: The dual-head Helix excels in both. **Real-time** body contouring maximizes Whole-Body **and** tomographic resolution, consistently, in every scan. Ultra-flared fan-beam collimation yields better than 6.0 mm resolution, quadrupling brain SPECT efficiency. And innovative **all-digital** camera design ensures unsurpassed imaging in every application.



HelixTM

A Continuous Revolution in Nuclear Imaging




Elscint

The Intelligent Image

Elscint/U.S.A.: (201) 342-2020, 1-800-ELSCINT.

Elscint/Belgium: (2) 720.92.46 □ Elscint/Brazil: (11) 869-4644 □ Elscint/Canada: (416) 474-1229 □ Elscint/Central & Eastern Europe, Austria: (1) 9855-681
□ Elscint/France: (1) 48-57-08-18 □ Elscint/Germany: (61) 22-7070 □ Elscint/Hong Kong: (5) 292231 □ Elscint/Israel: (04) 310310 □ Elscint/Italy: (2) 39320603
□ Elscint/Mexico: (5) 254-5939 □ Elscint/Spain: (3) 209.21.99 □ Elscint/U.K.: (923) 239511.



How to recognize a candidate for Cardiolite®

The shape of your patients may help you recognize the potential for soft-tissue attenuation, especially in fleshy figures.

For female and large-chested or obese male patients, Cardiolite comes through with higher photon energy (140 keV) to provide images with greater anatomical detail. Clear images can enhance interpretive confidence—which may reduce false-positives and equivocal cases.

Cardiolite also offers the unique advantage of direct measurement of both myocardial perfusion and ventricular function from one study.

So the next time you're faced with imaging female and large-chested or obese male patients, use Cardiolite and reduce soft-tissue attenuation.

Cardiolite®

Kit for the preparation of Technetium Tc99m Sestamibi

*To reduce soft-tissue attenuation
Cardiolite comes through*



Stress testing should be performed only under the supervision of a qualified physician in a laboratory equipped with appropriate resuscitation and support apparatus. There have been infrequent reports of signs and symptoms consistent with seizure and severe hypersensitivity after administration of Tc99m Sestamibi.

Please see brief summary of prescribing information on adjacent page.

© 1994, DuPont Pharma

Brief Summary

Cardiolite®

Kit for the preparation of Technetium Tc99m Sestamibi

FOR DIAGNOSTIC USE

DESCRIPTION: Each 5ml vial contains a sterile, non-pyrogenic, lyophilized mixture of: Tetrakis (2-methoxy isobutyl isonitrile) Copper (I) tetrafluoroborate - 1.0mg Sodium Citrate Dihydrate - 2.6mg L-Cysteine Hydrochloride Monohydrate - 1.0mg Mannitol - 20mg Stannous Chloride, Dihydrate, minimum (SnCl₂•2H₂O) - 0.025mg Stannous Chloride, Dihydrate, (SnCl₂•2H₂O) - 0.075mg Tin Chloride (Stannous and Stannic) Dihydrate, maximum (as SnCl₂•2H₂O) - 0.086mg

Prior to lyophilization the pH is 5.3-5.9. The contents of the vial are lyophilized and stored under nitrogen.

This drug is administered by intravenous injection for diagnostic use after reconstitution with sterile, non-pyrogenic, oxidant-free Sodium Perchnetate Tc99m Injection. The pH of the reconstituted product is 5.5 (5.0-6.0). No bacteriostatic preservative is present.

The precise structure of the technetium complex is Tc99m[MIBI]₆⁺ where MIBI is 2-methoxy isobutyl isonitrile.

INDICATIONS AND USAGE: CARDIOLITE, Kit for the preparation of Technetium Tc99m Sestamibi is a myocardial perfusion agent that is useful in the evaluation of ischemic heart disease. CARDIOLITE, Kit for the preparation of Technetium Tc99m Sestamibi is useful in distinguishing normal from abnormal myocardium and in the localization of the abnormality, in patients with suspected myocardial infarction, ischemic heart disease or coronary artery disease. Evaluation of ischemic heart disease or coronary artery disease is accomplished using rest and stress techniques.

CARDIOLITE, Kit for the preparation of Technetium Tc99m Sestamibi, is also useful in the evaluation of myocardial function using the first pass technique.

Rest-exercise imaging with Tc99m Sestamibi in conjunction with other diagnostic information may be used to evaluate ischemic heart disease and its localization.

In clinical trials, using a template consisting of the anterior wall, inferior-posterior wall and isolated apex, localization in the anterior or inferior-posterior wall in patients with suspected angina pectoris or coronary artery disease was shown. Disease localization isolated to the apex has not been established. Tc99m Sestamibi has not been studied or evaluated in other cardiac diseases.

It is usually not possible to differentiate recent from old myocardial infarction or to differentiate recent myocardial infarction from ischemia.

CONTRAINDICATIONS: None known.

WARNINGS: In studying patients in whom cardiac disease is known or suspected, care should be taken to assure continuous monitoring and treatment in accordance with safe, accepted clinical procedure. Infrequently, death has occurred 4 to 24 hours after Tc99m Sestamibi use and is usually associated with exercise stress testing (See Precautions).

PRECAUTIONS:

GENERAL

The contents of the vial are intended only for use in the preparation of Technetium Tc99m Sestamibi and are not to be administered directly to the patient without first undergoing the preparative procedure.

Radioactive drugs must be handled with care and appropriate safety measures should be used to minimize radiation exposure to clinical personnel. Also, care should be taken to minimize radiation exposure to the patients consistent with proper patient management.

Contents of the kit before preparation are not radioactive. However, after the Sodium Perchnetate Tc99m Injection is added, adequate shielding of the final preparation must be maintained.

The components of the kit are sterile and non-pyrogenic. It is essential to follow directions carefully and to adhere to strict aseptic procedures during preparation.

Technetium Tc99m labeling reactions involved depend on maintaining the stannous ion in the reduced state. Hence, Sodium Perchnetate Tc99m Injection containing oxidants should not be used.

Technetium Tc99m Sestamibi should not be used more than six hours after preparation.

Radiopharmaceuticals should be used only by physicians who are qualified by training and experience in the safe use and handling of radionuclides and whose experience and training have been approved by the appropriate government agency authorized to license the use of radionuclides.

Stress testing should be performed only under the supervision of a qualified physician and in a laboratory equipped with appropriate resuscitation and support apparatus.

The most frequent exercise stress test endpoints, which resulted in termination of the test during controlled Tc99m Sestamibi studies (two-thirds were cardiac patients) were:

Fatigue	35%
Dyspnea	17%
Chest Pain	16%
ST-depression	7%
Arrhythmia	1%

Carcinogenesis, Mutagenesis, Impairment of Fertility

In comparison with most other diagnostic technetium labeled radiopharmaceuticals, the radiation dose to the ovaries (1.5rads/30mCi at rest, 1.2 rads/30mCi at exercise) is high. Minimal exposure (ALARA) is necessary in women of childbearing capability. (See Dosimetry subsection in DOSAGE AND ADMINISTRATION section.)

The active intermediate, [Cu(MIBI)]₂BF₄, was evaluated for genotoxic potential in a battery of five tests. No genotoxic activity was observed in the Ames, CHO/HPRT and sister chromatid exchange tests (all *in vitro*). At cytotoxic concentrations (≥ 20μg/ml), an increase in cells with chromosome aberrations was observed in the *in vitro* human lymphocyte assay. [Cu(MIBI)]₂BF₄ did not show genotoxic effects in the *in vitro* mouse micronucleus test at a dose which caused systemic and bone marrow toxicity (9mg/kg, > 600 × maximal human dose).

Pregnancy Category C

Animal reproduction and teratogenicity studies have not been conducted with Technetium Tc99m Sestamibi. It is also not known whether Technetium Tc99m Sestamibi can cause fetal harm when administered to a pregnant woman or can affect reproductive capacity. There have been no studies in pregnant women. Technetium Tc99m Sestamibi should be given to a pregnant woman only if clearly needed.

Nursing Mothers
Technetium Tc99m Perchnetate is excreted in human milk during lactation. It is not known whether Technetium Tc99m Sestamibi is excreted in human milk. Therefore, formula feedings should be substituted for breast feedings.

Pediatric Use

Safety and effectiveness in children below the age of 18 have not been established.

ADVERSE REACTIONS: During clinical trials, approximately 8% of patients experienced a transient metallic or bitter taste immediately after the injection of Technetium Tc99m Sestamibi. A few cases of transient headache, flushing and non-itching rash have also been attributed to administration of the agent. Cases of angina, chest pain, and death have occurred (See WARNINGS and PRECAUTIONS). The following adverse reactions have been rarely reported: signs and symptoms consistent with seizure occurring shortly after administration of the agent; transient arthritis in the wrist joint; and severe hypersensitivity, which was characterized by dyspnea, hypotension, bradycardia, asthma and vomiting within two hours after a second injection of Technetium Tc99m Sestamibi.

DOSAGE AND ADMINISTRATION: The suggested dose range for I.V. administration in a single dose to be employed in the average patient (70kg) is:
370-1110MBq (10-30mCi)

The dose administered should be the lowest required to provide an adequate study consistent with ALARA principles (see also PRECAUTIONS).

When used in the diagnosis of myocardial infarction, imaging should be completed within four hours after administration.

The patient dose should be measured by a suitable radioactivity calibration system immediately prior to patient administration. Radiochemical purity should be checked prior to patient administration.

Parenteral drug products should be inspected visually for particulate matter and discoloration prior to administration whenever solution and container permit.

Store at 15-25°C before and after reconstitution.

RADIATION DOSIMETRY: The radiation doses to organs and tissues of an average patient (70kg) per 1110MBq (30mCi) of Technetium Tc99m Sestamibi injected intravenously are shown in Table 4.

Table 4. Radiation Absorbed Doses from Tc99m Sestamibi

Organ	Estimated Radiation Absorbed Dose			
	REST		STRESS	
	2.0 hour void	4.8 hour void	2.0 hour void	4.8 hour void
	rads/30mCi	mGy/1110MBq	rads/30mCi	mGy/1110MBq
Breasts	0.2	2.0	0.2	1.9
Gallbladder Wall	2.0	20.0	2.0	20.0
Small Intestine	3.0	30.0	3.0	30.0
Upper Large Intestine Wall	5.4	55.5	5.4	55.5
Lower Large Intestine Wall	3.9	40.0	4.2	41.1
Stomach Wall	0.6	6.1	0.6	5.8
Heart Wall	0.5	5.1	0.5	4.9
Kidneys	2.0	20.0	2.0	20.0
Liver	0.6	5.8	0.6	5.7
Lungs	0.3	2.8	0.3	2.7
Bone Surfaces	0.7	6.8	0.7	6.4
Thyroid	0.7	7.0	0.7	6.8
Ovaries	1.5	15.5	1.6	15.5
Testes	0.3	3.4	0.4	3.9
Red Marrow	0.5	5.1	0.5	5.0
Urinary Bladder Wall	2.0	20.0	4.2	41.1
Total Body	0.5	4.8	0.5	4.8

Organ	STRESS			
	2.0 hour void	4.8 hour void	2.0 hour void	4.8 hour void
	rads/30mCi	mGy/1110MBq	rads/30mCi	mGy/1110MBq
Breasts	0.2	2.0	0.2	1.8
Gallbladder Wall	2.8	28.9	2.8	27.8
Small Intestine	2.4	24.4	2.4	24.4
Upper Large Intestine Wall	4.5	44.4	4.5	44.4
Lower Large Intestine Wall	3.3	32.2	3.3	32.2
Stomach Wall	0.5	5.3	0.5	5.2
Heart Wall	0.5	5.6	0.5	5.3
Kidneys	1.7	16.7	1.7	16.7
Liver	0.4	4.2	0.4	4.1
Lungs	0.3	2.6	0.2	2.4
Bone Surfaces	0.6	6.2	0.6	6.0
Thyroid	0.3	2.7	0.2	2.4
Ovaries	1.2	12.2	1.3	13.3
Testes	0.3	3.1	0.3	3.4
Red Marrow	0.5	4.6	0.5	4.4
Urinary Bladder Wall	1.5	15.5	3.0	30.0
Total Body	0.4	4.2	0.4	4.2

Radiopharmaceutical Internal Dose Information Center, July 1990, Oak Ridge Associated Universities, P.O. Box 117, Oak Ridge, TN 37831, (615) 576-3449.

HOW SUPPLIED: Du Pont Radiopharmaceutical's CARDIOLITE®, Kit for the Preparation of Technetium Tc99m Sestamibi, is supplied as a 5ml vial in kits of two (2), five (5) and thirty (30) vials, sterile and non-pyrogenic.

Prior to lyophilization the pH is between 5.3-5.9. The contents of the vials are lyophilized and stored under nitrogen. Store at 15-25°C before and after reconstitution. Technetium Tc99m Sestamibi contains no preservatives. Included in each two (2) vial kit are one (1) package insert, six (6) vial shield labels and six (6) radiation warning labels. Included in each five (5) vial kit are one (1) package insert, six (6) vial shield labels and six (6) radiation warning labels. Included in each thirty (30) vial kit are one (1) package insert, thirty (30) vial shield labels and thirty (30) radiation warning labels.

The U.S. Nuclear Regulatory Commission has approved this reagent kit for distribution to persons licensed to use byproduct material pursuant to section 35.11 and section 35.200 of Title 10 CFR Part 35, to persons who hold an equivalent license issued by an Agreement State, and, outside the United States, to persons authorized by the appropriate authority.



Marketed by
Du Pont Radiopharmaceutical Division
The Du Pont Merck Pharmaceutical Co.
331 Treble Cove Road
Billerica, Massachusetts, USA 01862
For ordering Tel. Toll Free: 800-225-1572
All other business: 800-362-2668
(For International, call 617-350-9332)

SPECT BRAIN IMAGING CLINICAL FELLOWSHIP

Department of Radiology
Section of Nuclear Medicine



BENEFIT

This program is designed for nuclear medicine physicians, radiologists, technologists and referring physicians. It is intended to educate participants about the clinical utility of SPECT brain imaging with agents such as Ceretec® and Neurolite®.

Objectives include:

- Development of interpretation skills for brain images.
- Appreciation of clinical applications of SPECT brain imaging.
- Knowledge of image acquisition and reconstruction.
- Appreciation of factors that influence image quality.
- Knowledge of quality control techniques for SPECT.

SPONSORSHIP:

This program is sponsored by the Medical College of Wisconsin.

TUITION:

The tuition fee of \$650 includes the course syllabus, handouts, breaks, breakfasts, lunches, and other amenities involved in making this a pleasant learning experience. Maximum enrollments have been established. Cancellations prior to the course will be refunded, less a \$30 administrative fee.

CREDIT:

The Medical College of Wisconsin is accredited by the Accreditation Council for Continuing Medical Education to sponsor continuing medical education for physicians.

Accordingly, the Medical College of Wisconsin designates this continuing medical education activity as meeting the criteria for 13.00 hours in Category I toward the Physician's Recognition Award of the American Medical Association.

Nuclear Medicine Technologists who attend the SPECT Brain imaging Clinical Fellowship are eligible for 1.0 VOICE credit.

Register me for the following dates: (Please indicate a second choice)

- September 12-13, 1994 November 14-15, 1994

I will need reservations for _____ Sunday and Monday
night / _____ only on Monday night,
I will need a _____ single / _____ double room.

A check in the amount of \$650 should accompany this registration form and be made payable to the Medical College of Wisconsin. Telephone registrations must be confirmed by check within 10 days.

Name _____

Address _____

City/State/Zip _____


Office Phone _____

- work address home address

Registrations and payment should be sent to:


LisaAnn Trembath
SPECT Brain Imaging Fellowship Coordinator
Nuclear Medicine Division
Medical College of Wisconsin
8700 W. Wisconsin Avenue
Milwaukee, WI 53226 • (414) 777-3756

DON'T USE A FOIL COLLIMATOR...



FOR BETTER IMAGE QUALITY

USE A MICRO CAST COLLIMATOR...



FROM VAN MULLEKOM NUCLEAR FIELDS

NUCLEAR FIELDS (USA) CORP.

Phone: (708)299-8450 Fax: (708)299-8452

Circle Reader Service No. 132

WHAT DOES A SMALL CLINIC IN THE SOUTHEAST HAVE IN COMMON WITH A MAJOR MEDICAL CENTER IN THE MIDWEST



They both subscribe to
Nuclear Medicine Update.

Today, over 90,000 copies of *Nuclear Medicine Update* are sent annually to both the nuclear medicine community's referring physicians and Medicare Carrier Medical Directors. Published jointly by the Central Chapter of The Society of Nuclear Medicine and the American College of Nuclear Physicians, *Nuclear Medicine Update* offers your referring staff timely features about nuclear medicine procedures. Join the growing family of *Nuclear Medicine Update* subscribers, now!

For subscription information, please contact:

University Nuclear Medicine, Ltd.
3651 Red Bud Court, Downers Grove, IL 60515-1352
1-800-382-4541 • FAX 1-708-971-8103

*XPert*TM

Make an *expert* decision



Elscint

***XPert* Aided Diagnosis**

A built-in knowledge-base and advanced algorithms empower *XPert* to analyze data intelligently, infusing its processing modules with NM technical expertise. Propelled by superscalar micro-computing power, *XPert* helps you expand the frontiers of nuclear imaging.


XPert's Toolbox includes interactive graphical tools for high precision lesion delineation. And smart image structure interpretation totally automates SPECT reconstruction, homing in on target tissues, without operator intervention.

An *XPert* Display of Power

A multi-processor array of Intel Pentium and RISC number crunchers energizes *XPert* with 122 Mflops / 150 MIPS for 30 msec/slice SPECT reconstruction speed. A 30 Mpixel/sec graphic engine with 1280x1024 display reveals lesions with remarkable sharpness. Advanced clinical macro-programming expands *XPert's* diagnostic power with the vast repertoire of CLIP programs, developed by thousands of Apex users over the last decade.

***XPert* Link, Universal Connectivity**

XPert is a great communicator. Reaching out beyond its total link with Elscint imagers, it networks equally well with others. And *XPert* provides a transparent digital connection to fine-resolution laser multi-imagers and color printers. Superior PACS capabilities yield optimal equipment-use, boosting departmental productivity and cutting equipment costs.



*XPert*TM

Knowledge is power... the power of the expert.



Elscint

The Intelligent Image



Elscint/U.S.A.: (201) 342-2020, 1-800-ELSCINT.

Elscint/Belgium: (2) 720.92.46 □ Elscint/Brazil: (11) 869-4644 □ Elscint/Canada: (416) 474-1229 □ Elscint/Central & Eastern Europe, Austria: (1) 9855-681

Elscint/France: (1) 48-57-08-18 □ Elscint/Germany: (61) 22-7070 □ Elscint/Hong Kong: (5) 292231 □ Elscint/Israel: (04) 310310 □ Elscint/Italy: (2) 39320603

Elscint/Mexico: (5) 254-5939 □ Elscint/Spain: (3) 209.21.99 □ Elscint/U.K.: (923) 239511.

CELEBRATE NUCLEAR MEDICINE WEEK OCTOBER 2-8 1994



Nuclear Medicine Week—
October 2 through 8—is the prime time to
demonstrate pride in your profession — and to
make the profession's presence known both
among the public and other health care
professionals.

Under the sponsorship of the Society of
Nuclear Medicine and SNM's Technologist
Section, Nuclear Medicine Week offers an
excellent opportunity to educate, to stimulate,
and to promote the successes of nuclear
medicine. This week also gives you a specific
time to spotlight your facility to referring
physicians, potential patients, and to anyone
else in your community who could benefit
from nuclear medicine.

To help enhance the visibility of nuclear
medicine facilities, the Technologist Section,
has designed a striking new poster to mark this
year's event. Buttons, stickers, and guidelines

are also available to assist you with your celebration. We guarantee that this year's
sensational design — will draw attention and spur positive comment.

*P.S. Don't Forget Syncor's PR Star Contest — details in "Guidelines." Packet. Look for
details in upcoming journals. Be a Public Relations Star and win prizes for yourself
and your institution.*

CELEBRATE NUCLEAR MEDICINE WEEK

The following materials are available for promoting Nuclear Medicine Week in your area. One poster, sticker, and a button, all in full color, have been designed for this year.

POSTERS — \$5.00 each, 4-9 posters are \$4.50 each, 10 or more \$4.00 each.

I would like _____ posters × \$ _____ \$ _____

BUTTONS — \$1.00 each

I would like to order ____ buttons \$ _____

STICKERS — \$.25 each (same design as the button)

I would like to receive ____ stickers.
(Minimum order is 10 stickers) _____
\$ _____

Total \$ _____

I would like to order a free set of “Guidelines for Promoting Nuclear Medicine Week.”

Payment must be enclosed with your order. Payments must be made in U.S. dollars drawn on U.S. banks. No foreign funds will be accepted. Make checks payable to

The Society of Nuclear Medicine

Orders will be sent out by 1st class mail or UPS. Orders received after Sept. 1, 1994 will be assessed a 15% surcharge, payable before shipment, to ensure timely delivery.

Name _____

Hospital/Company _____

Address _____

City _____ State _____ Zip _____

Telephone _____

Please return this form to:
Nuclear Medicine Week
The Society of Nuclear Medicine
136 Madison Avenue
New York, NY 10016-6760

Policy—The *Journal of Nuclear Medicine* accepts classified advertisements from medical institutions, groups, suppliers, and qualified specialists in nuclear medicine. Acceptance is limited to Positions Open, Positions Wanted, and Equipment. We reserve the right to decline, withdraw, or modify advertisements.

Rates for Classified Listings—\$22.00 per line or fraction of line (approx. 50 characters per line, including spaces). Please allow 28 characters for the first line which will appear in capital letters. Special rates for SNM members on Positions Wanted. \$10.00 per line. *Note: Box numbers are available for the cost of the 2 lines required.*

Rates for Display Ads—Agency commissions are offered on display ads only

Full page	\$1500	Quarter page	\$700
Half Page	950	Eighth page	550

Publisher-Set Charges—Page \$150; half page \$100; quarter page \$75; eighth page \$50.

Terms—Payment must accompany order. Make checks payable, in U.S. dollars on U.S. banks only, to: The Society of Nuclear Medicine.

Deadline—First of the month preceding the publication date (January 1 for February issue). Please submit classified listings typed double spaced. No telephone orders are accepted.

Send Copy to—Classified Advertising Department, The Society of Nuclear Medicine, 136 Madison Avenue, New York, NY 10016-6760, (212) 889-0717 FAX: (212) 545-0221

Positions Available

Physician

Position in Mid-Atlantic medium size growing community hospital, combined Nuclear Medicine and Pathology practice (AP/CP). Send CV to: Box 801, The Society of Nuclear Medicine, 136 Madison Avenue, New York, NY 10016.

LOCUMS: 100% NM practice 6-10 wks/yr. Excellent Dept. Send CV to Dr. Cheng, 3118 Colyar Dr., Chattanooga, TN 37404, (615) 495-8736.

NORTHERN CALIFORNIA - The Kaiser Permanente Medical Center in Santa Clara, CA is seeking a BCBE Nuclear Medicine Physician for a career opportunity with the nation's leading HMO. All aspects of nuclear medicine services are provided to a 250,000 member prepaid patient population. Clinical and administrative experience required. Internal medicine background preferred. Academic opportunities are available. Competitive salary, generous benefits and a comprehensive retirement program. For more information, send CV with cover letter to: Diane Butler, The Permanente Medical Group, Inc. Physician Recruitment, Dept. 68, 1814 Franklin, 4th floor, Oakland, CA 94612. EOE.

Radiologist

RADIOLOGIST/NUCLEAR MEDICINE- 5 person NY/NJ group seeking radiologist with special competency in Nuclear Medicine. Interest in mammography desired but not essential. Young, progressive group located in 400 bed hospital with nearby imaging center. Send CV to: James Heimann, M.D., 5 Franklin Ave., Belleville, NJ 07109; (201) 450-2038, (201) 751-2011.

Technologist

TECHNOLOGIST TRAINING PROGRAM
Positron Emission Tomography Technologist Training Program. A busy positron diagnostic center is offering a six-month training program to certified (NMTCB or AART) nuclear medicine technologists. Included in the training program is a 24-hour lecture series covering topics such as physics, chemistry, computer applications and patient care as applied to PET imaging, plus cyclotron operation and radiopharmaceutical preparation. Much of the program will be in the form of closely supervised on-the-job training. The student will be instructed in the operation of the positron camera, associated computers and the cyclotron. Tuition is charged. Interested parties should contact the Department of Nuclear Medicine, William Beaumont Hospital, Royal Oak, MI 48073-6769, (810) 551-4126.

Computers in Nuclear Medicine: A Practical Approach

Kai Lee, PhD

Computers have become an indispensable tool in nuclear medicine. This is the book for those who wish to acquire a basic understanding of how computers work and the processing techniques used to obtain diagnostic information from radionuclide images. The text gives a thorough description of the hardware components of a nuclear medicine computer system and explains the principles behind many common image processing techniques. The following topics are discussed in detail:

- Functions and components of a computer system
- Mass storage devices
- Input and output devices
- Computer software
- Nuclear medicine image acquisition methods
- Methods of qualitative image analysis
- Quantitative image analysis
- Nuclear cardiology
- Quantitative data analysis
- Single-photon emission computed tomography
- Selecting a computer for nuclear medicine

The book is illustrated throughout to help the reader conceptualize the topics as they are discussed.

290 pp, 6 x 9, softcover

\$30 member (+ \$2.50) Total \$32.50


\$45 nonmember (+ \$2.50) Total \$47.50

In Canada, add \$5.00; elsewhere, add \$20.00
for shipping and handling

To order, send payment to:

Bookmasters
State Rt. 42, RD11
Mansfield, OH 44903
1 (800) 247-6553

Confidence Builder . . .



Build a solid foundation as you prepare for the NMTCB Examination.

Increase the effectiveness of your study time.

SNM's Review of Nuclear Medicine Technology is the best single study aid you can own as you prepare for certification exams. Current, authoritative, thorough—the *Review* is a valuable addition to the libraries of students and specialists alike. Practical appendices cover—

- Test-taking techniques*
- Sample questions and answers*
- Pertinent NRC regulations.*

**To order, call toll-free, Bookmasters, Inc.,
1-800-247-6553**

(Outside the U.S. 419-281-1802).

The Finest Line of Cardiac Gates Available



For over fourteen years, Advanced Medical Research, now known as AccuSync Inc., has been serving the cardiac health care industry with the finest line of cardiac gates available in today's market.

Our dedication to service and commitment to provide you with a reliable product have built the reputation of our gates.


With a complete line of models available, you are able to choose the gate which best corresponds to your specific requirements.

The AccuSync 5L, our top model (featured at left) includes CRT monitor (visual) and Strip Chart Recorder (hard copy).

Model Specifications:

- Auto/Manual trigger control
- No delay
- ECG output
- Audio indicator
- Trigger pulse LED
- Isolation amplifier for patient safety
- Compatible with all computers

AccuSync models 5L, 6L and 1L are CSA and ETL (UL544) approved



Model	Strip Chart	CRT Monitor	HR/R-R Int	Trigger
5L	•	•	•	•
6L		•	•	•
1L	•		•	•
3L			•	•
4M				•

Accessory and optional products available:
The AccuAmp 5, the 5 lead system available for AccuSync 5L, 6L, and 1L, transmits information through fiber optic link. Patient cables, lead wires, and BNC cables available for AccuSync models.