## **Correct Use of Dose Calibrator Values**

**FO THE EDITOR:** In a recent article by Salako and DeNardo (1) the nuthors discussed the relative accuracy of two recommended dose calibrator settings for assaying solutions of  $^{90}$ Y using a commercially available lose calibrator. The first setting, "775 × 100," is recommended by Capintec, Inc. (Pittsburgh, PA) (2) for measuring  $^{90}$ Y in a standard National Institute of Standards and Technology (NIST)-style 5-ml glass impoule. This setting was known to result in inaccurate assays, which rompted the second calibration factor, "48 × 10," to be experimentally etermined by Coursey et al. (3) in the same sample geometry with direct etermination of the solution activity using several detection methods. In the present article (1) the authors attempt to arrive at a multiplicative factor that would somehow salvage the erroneous "775" setting, recommending a calibration factor of "775 × 70."

While it is laudable that these workers recognize the importance of using the correct dose calibrator setting, we feel that there are several points that need to be clarified:

- 1. Dose calibrator settings depend strongly on the geometry of the sample, especially for low-energy photon emitters. The calibrator settings recommended by most dose calibrator manufacturers are valid only for the standard NIST geometry; that is, a 5-ml flame-sealed, thin-walled glass ampoule. Radiopharmaceuticals are usually shipped from the manufacturer in vials or single-dose plastic syringes, both of which have different photon absorption characteristics than the geometry for which the calibration factor was determined. The result is that a discrepancy of 10% or more can be observed if the incorrect calibration factor is applied for a particular geometry. It is imperative that the correct calibration factor be applied for the particular container used in the radioassay.
- 2. The dose calibrator is sensitive not only to the type of container used but also to the filling volume of solution in the container. This effect is clearly seen in the glass vial data presented in Table 3 from Salako and DeNardo's article (1). As the vial is filled, the apparent activity appears to decrease. This is most likely due to photon absorption by the source liquid, leading to a lower response in the dose calibrator.
- 3. The multiplicative factor applied to calibration factors that fall outside the normal operational range of the instrument potentiometer (such as "48 × 10," "775 × 100," etc.) are usually chosen so as to provide a convenient way to calculate the activity from the dose calibrator display. One can easily find any combination of dial settings and multiplicative factors that can give the correct activity, and they would all be equally appropriate. However, most people find it easier to multiply or divide multiples of 10. Therefore, despite the fact that Salako and DeNardo (1) have apparently determined a multiplicative factor for the "775" setting that will provide the correct activity, we still recommend the use of the "48 × 10" setting for <sup>90</sup>Y on the basis of ease of use.
- 4. The experimental design of the study was somewhat flawed in that one of the aims as outlined in the beginning of the article was to investigate the variability of shipments of <sup>90</sup>Y using a dose calibrator. The correct approach to the problem should have been to first determine the calibration factor using a standardized source of <sup>90</sup>Y. After establishing the correct setting, the sources could then have been assayed. As for the liquid scintillation measurements, it is not known why the author's liquid scintillation activity measurements differed from the manufacturer's value by about 12%, but relying on the manufacturer's stated activity only introduced an additional uncertainty. A better approach would have been to trace the activities

of the manufacturer's shipments against a standardized sample of  $^{90}\mbox{Y}_{.}$ 

In summary, it should be kept in mind that dose calibrators are extremely useful and reliable instruments as long as the correct calibration factor is applied for the radionuclide of interest in the geometry for which the calibration factor was determined. When attempting to determine a new calibration factor, it is imperative that a standardized solution of the radionuclide be used for all measurements.

## REFERENCES

- Salako QA, DeNardo SJ. Radioassay of yttrium-90 using the radionuclide dose calibrator. J Nucl Med 1997;38:723-726.
- Capintec, Inc. Radioisotope calibrator owner's manual, revised ed. Pittsburgh, PA: Capintec, Inc.; 1986.
- Coursey BM, Calhoun JM, Cessna JT. Radioassays of yttrium-90 used in nuclear medicine. Nucl Med Biol 1993;20:693-699.

Brian E. Zimmerman Bert M. Coursey Jeffrey T. Cessna

Ionizing Radiation Division, Physics Laboratory National Institute of Standards and Technology Gaithersburg, Maryland

REPLY: We acknowledge the interest of Drs. Zimmerman, Coursey and Cessna in our article on radioassay of <sup>90</sup>Y by radionuclide calibrator (1). Our work was actually a follow-up to that of Coursey et al. (2) in accordance with their experimental plan and recommendations. In their work, Coursey et al. (2) standardized a solution of <sup>90</sup>Y supply by a tritium-tracing method (3), then used standard solutions prepared from this by serial dilution to characterize radionuclide calibrators, liquid and solid scintillators, and Cerenkov counters. In consonant with their results, Coursey et al. (2) advised prospective <sup>90</sup>Y users that "alternatively, solutions may be standardized by liquid-scintillation counting using the method" they described, "and those working solutions may be used to calibrate the radionuclide calibrators." In our work (1), we prepared counting solutions (10 replicates) from a commercial sample of 90YCl, solution after the serial dilution technique of Coursey et al. (2), then established the radioactivity contents (MBq) in these test solutions by liquid scintillation counting (LSC) against tritium standards. Standard solutions were then prepared from the commercial 90YCl3 solution and used to calibrate our radionuclide calibrator. Our results corroborated those of Coursey et al. (2) in terms of the right radionuclide calibrator dial setting for <sup>90</sup>Y calibration. Consequently, we regard the comments of Zimmerman et al., in their letter, as pertinent to the import of our message in the publication under reference 1. The following are our respective reactions to their comments:

1. We support the statement on the importance of using the right calibration factor for radioassay of <sup>90</sup>Y, however, it is important to stress that calibration factors and container correction factors are not the same. Coursey et al. (2) re-established the appropriateness of calibrator dial setting 48 (with multiplicative factor 10) for measuring a <sup>90</sup>Y source contained in a thin-walled glass ampoule. We have also confirmed this same dial setting 48 × 10 as a correct calibration factor for calibrating <sup>90</sup>Y sources contained in 1-ml glass vials. The essential point might be to recommend a particular container (e.g., the thin-walled glass since the reference <sup>90</sup>Y source is contained in this) as a reference container to which other containers are related with appropriate container correction factors, but this we believe many <sup>90</sup>Y suppliers are doing already. Therefore, the statement of Zimmerman et al. that "it is imperative that the correct calibration factor be applied for the particular container used in the radioassay"

- should not be construed to mean that different calibrator dial settings should be used for radioassay of different <sup>90</sup>Y sources.
- We agree with the explanation for why dose calibrator readings decrease for <sup>90</sup>Y source in glass vials after a certain level of filling or dilution.
- 3. We have actually studied <sup>90</sup>Y calibration at dial 775 (in addition to dial 48) not only because it was originally recommended by Capintec, Inc. (Pittsburgh, PA) (4) and used by many laboratories, but also because our major <sup>90</sup>Y supplier then, always calibrated their shipments at this dial with a multiplicative factor of 60. In addition to corroborating the accuracy of measurement at the recommended dial 48 with a multiplicative factor of 10 (2), we have also found that the more accurate multiplicative factor for dial 775 is 70, not 60.
- 4. Zimmerman and his colleagues were probably misled into believing that our experimental design was flawed by one of our statement objectives, which was to investigate variabilities in 90 calibration among suppliers. But our experiment was approached as suggested by Zimmerman et al. in their letter. Further, before setting up our experiment, we contacted Dr. Coursey by phone and he clarified the issues involved in the standardization of 90 solutions by the use of 90 reference material or LSC. The 90 reference material was not available from the National Institute of Science and Technology (NIST) at the time, and we followed the "alternative standardization" advice of Coursey et al. (2). The results of our experiment firmly established for us the use of calibrator dial 48 (with multiplicative factor of 10) to calibrate all 90 sources. Before this study, we have often observed disagreements between our calibration measurements

at check-ins  $(48 \times 10)$  and stated calibrations of respective suppliers. On inquiry, it appeared the suppliers were using different methods of  $^{90}$ Y calibration in their laboratories. After the experiment, we retroactively compared the check-in values in our records with suppliers' calibration methods and values, then advised each supplier appropriately.

In conclusion, we agree on the use of calibrator dial 48 (with a multiplicative factor of 10) as being convenient and now popular for calibrating <sup>90</sup>Y radiation. The results of our study further supported the expectation of achieving higher accuracy when calibration factors are investigated with standardized solutions from <sup>90</sup>Y reference material than from liquid scintillation counting.

## REFERENCES

- Salako QA, DeNardo SJ. Radioassay of yttrium-90 using radionuclide dose calibrator. J Nucl Med 1997;38:723-726.
- Coursey BM, Calhoun JM, Cessna JT. Radioassay of yttrium-90 used in nuclear medicine. Nucl Med Biol 1993;20:693-699.
- Coursey BM, Garcia-Torano E, Golas DB, et al. The standardization and decay scheme of rhenium-186. Appl Radiat Isot 1991;42:865-870.
- Capintec, Inc. Radioisotope calibrator owner's manual, revised ed. Pittsburgh, PA: Capintec, Inc.; 1986.

Qansy A. Salako

Division of Medical Imaging and Radiopharmaceutical Drug Products

U.S. Food and Drug Administration\*

Rockville, Maryland

This letter was written by Dr. Salako in his private capacity. No official support or endorsement by the Food and Drug Administration is intended or should be inferred.