



## This sponge

By eliminating the disadvantages of earlier methods, the Triosorb Sponge has achieved a real breakthrough in thyroid testing. **It is an in vitro test unmatched in accuracy, speed and convenience.**

**Accuracy:** Because factors such as red blood cells and exogenous iodine have been eliminated from consideration in the Triosorb Test, it is unmatched in accuracy.

**Speed:** With only 3 washes and no need for double pipettings, shakers, or incubators, the Triosorb Test can be

# revolutionized thyroid testing!

more rapidly performed than any other T-3 test.

**Convenience:** Triosorb is in a disposable kit ready for immediate use at room temperature, making it the simplest and most convenient thyroid function test to perform.

McAdams\* reported that "The resin sponge (Triosorb) technique is superior to the erythrocyte method for performing the  $I^{131}$  T3 test in terms of simplicity, convenience and elimination of errors characteristic of the erythrocyte procedure."

**Triosorb is available to all doctors, hospitals and clinical laboratories—AEC licensing is not required. Because Triosorb will enable far more screenings to be performed, this procedure may soon become as standard as today's blood counts and urinalyses.**




\*McAdams, G. B. and Reinfrank, R. F., Jnl. Nuclear Med., 5:112, Feb., 1964.

**TRIOSORB®**  
**T-3 DIAGNOSTIC KIT**  
**ABBOTT LABORATORIES NORTH CHICAGO, ILL.**

FOR BIOLOGICAL RESEARCH

...CLINICAL DIAGNOSIS

...SPECIAL STUDIES



## ARMAC® provides rapid, accurate measurement of radioactivity

The Armac Scintillation Detector is ideally suited for *in vivo* research using gamma-emitting tracers. Its 4-pi counting geometry and high counting efficiency require the administration of only small doses of radioactivity, yet assure significant measurements in short counting periods.

- In biological research Armac can be used to count radioactivity in small laboratory animals in the study of retention and excretion of various gamma-emitting isotopes. Because only small doses of radioactivity are required, sacrificing is unnecessary and the animal acts as its own control for repeated experimentation.
- In clinical research Armac is useful in nutritional

studies and other applications where changes in the amount of circulating radioactivity in the blood can be directly related to other body functions. These changes can be measured accurately and rapidly using the blood-rich portion of the patient's forearm positioned in the Armac's counting chamber.

- In special studies Armac can be used to quickly and accurately determine the presence of low levels of radioactivity in bulk samples of meat, milk, water, blood, wastes, soil, or tissue with little or no sample preparation. The photograph at right, above, shows an Armac Detector counting radioactivity in a live fish in flowing stream water. (Courtesy Oak Ridge National Laboratory.)

A complete counting system includes the Model 446 Armac Scintillation Detector and any of several Packard Spectrometers, control units and ratemeters. There are more than ten Armac Detector Systems available to meet your specific counting requirements. Choose visual or recorded data output; digital or analog data presentation. Any system can be readily modified or expanded to meet changing research requirements. Special systems incorporating automatic sample changers are also available.

For detailed information contact your Packard Sales Engineer, or write for Bulletin 1024.

**Packard**

**PACKARD INSTRUMENT COMPANY, INC.**

2200 WARRENVILLE ROAD • DOWNERS GROVE, ILLINOIS 60515

TELEPHONE: 312/969-6000



The Volk Radiochemical Company has been producing radiopharmaceuticals and radiochemicals for human and research use for the past ten years. This is our only business. Continued technological improvements and innovations have been constant goals as exemplified by our pioneering production of Iodine-125 as a useful medical isotope and our introduction of the "Silver Saddle" which removes free iodide from iodinated organic compounds such as Hippuran.

*Twenty separate individual isotopic dose products are available in standard and special potencies.*

**CONVENIENT**—to the world's busiest airport (Chicago's O'Hare) insures delivery tomorrow via Air Express to any one of 21,000 U. S. cities.

Collect call ordering  
to our  
Skokie, Illinois laboratory  
(area code 312 673-3760)



or Burbank, California  
laboratory  
(area code 213 849-6023)  
assures you that your orders  
are handled as specified.


*Write for the 1965 Volk Radiomedicines catalog.*



**RADIOCHEMICAL COMPANY**

CHICAGO—8260 Elmwood Ave., Skokie, Ill.  
LOS ANGELES—803 N. Lake St., Burbank, Calif.  
WASHINGTON—P.O. Box 335, Silver Spring, Md.  
NEW YORK—P.O. Box 345, New York, N. Y.

## T3 THYROID TEST



MODEL GWS - 78

### **TRANSISTORIZED GAMMA WELL SYSTEM**

**\$1495**

A NEW, LOW COST, Transistorized Gamma Well Scintillation System that you can afford. Includes a precision scaler, 2" Na I (T1) well crystal, scintillation probe with 2" Photomultiplier Tube, and manual shielded well scintillation sample changer. This entire system is available for immediate delivery. Also available are various timers for preset time operation (Model TM-1 mechanical timer \$125; transistorized electronic timer Model TM-12 \$500).

**SCALER - 300-3000 adjustable high voltage  
1 usec resolving time / 999,999 total count / Aural speaker  
5 uA input sensitivity / Current mode operation**




Manufactured exclusively for  
NUCLEAR SUPPLIES by the world's  
first company to successfully  
build transistorized nuclear equipment,  
Kobe Industries Corporation of Japan.

Complete line of instruments,  
modules, automatic systems and  
accessories listed in new catalog.  
Write to P.O. Box 312, Encino,  
California, or Telephone,  
Area Code 213; 787-1722.



# Large Crystal Radioisotope Scanners


OHIO-NUCLEAR  
Model 54



## FEATURES:

- Soundly engineered—research proven in use.
- 200 inches (500 cm) per minute scanning speeds.
- Choice of above or beneath table scanning.
- Dual Scanner capability—now or later.
- Gammagraphic Photoscan & Dot-tapper.
- Prompt service by capable factory engineers.
- 5-inch crystal scanner price—\$16,750.
- 8-inch crystal scanner price—\$20,700.

*For Full Information WRITE OR CALL US • WE BELIEVE YOU WILL BE GLAD YOU DID.*




**OHIO-NUCLEAR, INC.**


1725 FALL AVENUE

CLEVELAND, OHIO

216 - 621-8477



**organs**  
currently demonstrable by  
Pickier MagnaScanner techniques



 **MEDOTOPES®**  
**SQUIBB RADIOPHARMACEUTICALS**

Medotopes reflect the latest developments in nuclear medicine. All provide the utmost in safety and convenience. All have unique packaging safeguards so that direct contact is never required. Exclusive lead shield enclosures are fitted with bottle caps that unscrew automatically. Saf-Tag® vials and bottles are carefully encased and double protected by transparent, shatterproof plastic coatings, and shipping cartons have convenient "pull-tab" openers. And, each preparation is custom-handled, each delivery custom-routed by Squibb Traffic Service. Access to three major airports expedites shipment.

Squibb Radiopharmaceuticals are available to the AEC-licensed physician. For full information, write to Professional Service Dept., Squibb, 745 Fifth Avenue, New York 22, N. Y.

**SQUIBB**  
Squibb Quality—the Priceless Ingredient  
OWNED BY SQUIBB 


# A Valuable Addition To Your Professional Library

*Journal of*

**NUCLEAR**



*Official Publication  
Society of Nuclear Medicine*



## AN IMPORTANT NEW JOURNAL

*featuring*

Original articles in clinical medicine, basic and clinical medical research, physics and chemistry dealing with the use of isotopes in humans, and articles on related subjects. The latter includes dosimetry, instrumentation, protection, techniques, biologic effects contributing to the use or effects of isotopes in clinical medicine or the clinical effects of ionizing radiation.

George E. Thoma, M.D., St. Louis—*Editor*

G. O. Broun, Jr., M.D., St. Louis, Titus C. Evans, Ph.D., Iowa City,

Neil Wald, M.D., Pittsburgh, Eugene L. Saenger, M.D., Cincinnati—*Associate Editors*

**The Journal of NUCLEAR MEDICINE**  
333 North Michigan Avenue, Chicago, Illinois 60601  
\$20.00 per year, U.S.    \$21.00 Foreign

Name.....  
Address.....  
City..... State ..... Zip Code.....

*Please remit by check or money order.*

**Published  
Monthly**



## Announcement to Authors Preliminary Notes

Space will be reserved in each issue of THE JOURNAL OF NUCLEAR MEDICINE for the publication of one preliminary note concerning new original work that is an important contribution in Nuclear Medicine.

Selection of the preliminary note shall be on a competitive basis for each issue. One will be selected after careful screening and review by the Editors. Those not selected will be returned immediately to the authors without criticism. Authors may resubmit a rejected or revised preliminary note for consideration for publication in a later issue. The subject material of all rejected manuscripts will be considered confidential.

The text of the manuscript should not exceed 1200 words. Either two illustrations, two tables, or one illustration and one table will be permitted. An additional 400 words of text may be submitted if no tables or illustrations are required. Only the minimum number of references should be cited.

Manuscripts should be mailed to the Editor, Dr. George E. Thoma, St. Louis University Medical Center, 1402 South Grand Blvd., St. Louis, Missouri 63104. They must be received before the first day of the month preceding the publication month of the next issue, *e.g.*, preliminary notes to be considered for the November, 1965 issue must be in the hands of the Editor before October 1, 1965.


---

## Advertising Index

<b>Abbott Laboratories</b> North Chicago, Illinois Inside Front Cover, i	<b>Packard Instruments</b> Downers Grove, Illinois ..... iii
<b>Nuclear-Chicago</b> Des Plains, Illinois .... Back Cover	<b>Picker X-Ray Corp.</b> White Plains, New York ..... x
<b>Nuclear Consultants</b> St. Louis, Missouri Inside Back Cover	<b>Squibb, E. R. &amp; Sons</b> New York, New York ..... xi
<b>Nuclear Supplies</b> Encino, California ..... viii	<b>Tracerlab</b> Waltham, Massachusetts ..... xx
<b>Ohio Nuclear Corp.</b> Cleveland, Ohio ..... ix	<b>Volk Radiochemical</b> Chicago, Illinois ..... v

**undergraduate  
simplicity**

**post-doctoral  
performance**



## **Tracerlab's new multi-mission scaler**

The new model 132 M manual scaler is a versatile, convenient instrument that even a near-novice can use with professional aplomb — and get results!

What's more, this instrument utilizes high-reliability printed circuit boards and transistor circuitry throughout and features performance characteristics you'd hardly expect at its most reasonable price:

High speed (500 kc/sec) with 800-nanosecond resolving time, 3.5-millivolt sensitivity. Voltage-and-charge-sensitive inputs with 0.5% linearity discriminator and 3000:1 input sensitivity range. Other features include 0-3000 vdc power supply with better than 0.003% regulation, 0-60 minute electromechanical timer with 0.1 second accuracy, infinite preset times and visual readout. All this adds up to the very latest in scaler techniques — the 132 M — from the first family in nuclear instrumentation.



**T R A C E R L A B**

A Division of Laboratory For Electronics, Inc.  
WALTHAM, MASSACHUSETTS 02154

• Film Badge Service • Health Physics • Bioassays • Sources • Nuclear Instrumentation • Radiochemicals  
Radioactive Waste Disposal • Radiation Monitoring Instrumentation • Isotope Applications

# SINGLE SOURCE

For the New Short  
Half-Life or Low Energy  
Isotopes

## Chlormerodrin-Hg<sup>197</sup>

## Chlormerodrin-Hg<sup>203</sup>

 Proven brain and renal scanning agents. Hg<sup>197</sup> available at lower cost through NCC's "Controlled Decay" program.

## 1-Mercuri-2-Hydroxypropane-Hg<sup>197</sup>


 For improved spleen scanning. Much simpler and more convenient than heat treating tagged red cells. MHP-Hg<sup>203</sup> also available for experimental use.

## Iodine-I<sup>125</sup>

 I-<sup>125</sup> is a longer-life isotope than I-<sup>131</sup>, and offers significantly lower radiation dosage to the patient. Available as tagged albumin and other iodinated compounds.

## Technetium-99m

*For Scanning Brain...Thyroid...Liver...Lung*



### TechneKow® generator

Supplies ample material as sodium pertechnetate for daily use. Economical for as few as seven or eight brain scans per week. New #1810 Shielded Dispenser available — simple semi-automated method of obtaining sterile solution.

### "Instant Technetium"

 Sodium pertechnetate in sterile, pyrogen-free solution permits immediate use without the necessity for "milking," sterilization, or calibration.

### Technetium-Sulfide Colloid

 In sterile, pyrogen-free form, reported as providing excellent liver scans. Also usable for lung scanning following the technique of Quinn. 0517



**NUCLEAR CONSULTANTS  
CORPORATION**

Box 6172, Lambert Field - St. Louis, Mo. 63145  
314 PErshing 9-8927

CHICAGO • CLEVELAND • HOUSTON • LOS ANGELES • NEW YORK • SAN FRANCISCO • WASHINGTON, D.C.

# Every one of these prize-winning engineers is a "yes" man.

Here is a group of unabashed "yes" men. Nuclear-Chicago mechanical engineers—and prize-winners, too, for their design and engineering of instruments for clinical diagnosis.

With their colleagues in electronic engineering, they've gained a reputation for affirmative responses to the clinician's questions.

Such as: Is there an instrument that can handle large patient loads, faster than a scanner? Make dynamic studies of isotope distribution? Answer: Yes, our Pho/Gamma™ Scintillation Camera.

Or: Is there a unit that gives sharp, clear photo and dot scans? Without a lot of guesswork or dial-twisting? Answer: Yes, our Pho/Dot® Isotope Scanner.

Other "yes" men—scientists, engineers, technicians — are working for Nuclear-Chicago, too. They're saying "yes" and providing the right instrument at the right time to fulfill a wide variety of needs. In science, medicine, education, and industry.

And about those prizes. There are two, both won in the American Iron and Steel Institute's 1964-65 *Design in Steel* competition. One is the Citation of Excellence in Engineering of Commercial Equipment—won by the Pho/Dot engineering team. The other is the Citation of Excellence for Design of Commercial Equipment—awarded to the Pho/Gamma team.

Not bad for a bunch of "yes" men.

NUC: D-5-229




## NUCLEAR-CHICAGO

A DIVISION OF NUCLEAR-CHICAGO CORPORATION  
333 Howard, Des Plaines, Ill. 60018 U.S.A.

In Europe: Donker Curtiusstraat 7  
Amsterdam W, The Netherlands

*Left to right: William Guth,  
Joseph Kus, Ray Meeder,  
Donald Hausser, and  
Peter Vischer. Also in the picture:  
Model 6401 Pho/Gamma  
Scintillation Camera—electronics  
console at left and detector  
(with back-shield removed)  
at right.*



*Scientists and engineers interested in challenging career opportunities are invited to contact our personnel director.*