



This sponge

By eliminating the disadvantages of earlier methods, the Triosorb Sponge has achieved a real breakthrough in thyroid testing. **It is an in vitro test unmatched in accuracy, speed and convenience.**

Accuracy: Because factors such as red blood cells and exogenous iodine have been eliminated from consideration in the Triosorb Test, it is unmatched in accuracy.

Speed: With only 3 washes and no need for double pipettings, shakers, or incubators, the Triosorb Test can be

revolutionized thyroid testing!

more rapidly performed than any other T-3 test.

Convenience: Triosorb is in a disposable kit ready for immediate use at room temperature, making it the simplest and most convenient thyroid function test to perform.

McAdams* reported that "The resin sponge (Triosorb) technique is superior to the erythrocyte method for performing the I^{131} T3 test in terms of simplicity, convenience and elimination of errors characteristic of the erythrocyte procedure."

Triosorb is available to all doctors, hospitals and clinical laboratories—AEC licensing is not required. Because Triosorb will enable far more screenings to be performed, this procedure may soon become as standard as today's blood counts and urinalyses.




501207

*McAdams, G. B. and Reinfrank, R. F., *Jrnl. Nuclear Med.*, 5:112, Feb., 1964.

TRIOSORB®
T-3 DIAGNOSTIC KIT
ABBOTT LABORATORIES NORTH CHICAGO, ILL.

MODEL 12T-015



**TRANSISTORIZED
POWER SUPPLY**

Continuously adjustable 400-3000 volts
 Low voltage: Line regulation better than $\pm 0.05\%$
 High voltage: Line regulation $\pm 0.003\%$


MODEL 01T-002



**TRANSISTORIZED
LINEAR AMPLIFIER**

Input sensitivity: $0.05 \mu\text{A}$
 Clipping: Double DL or RC
 Current pulse can be integrated to provide voltage pulse output
 Rise time better than 50 nano sec.


MODEL 03T-002



**TRANSISTORIZED
PULSE HEIGHT ANALYZERS**

Three models available:
 Integral Discriminator
 Manual
 Automatic Scan


MODEL 11T-001



**TRANSISTORIZED
LINEAR RATE METER**

6 ranges—4 time constants
 Alarm circuit
 10mv recorder output
 Range suppress
 Input impedance: 50 ohms
 Input sensitivity: $10\text{-}500 \mu\text{A}$

MODEL PA-20



**TRANSISTORIZED
PREAMPLIFIER**

Gain: Adjustable, 70 to 400
 Rise time: 200 nano seconds
 Noise: Less than $0.05 \mu\text{A}$
 Power requirement: $-9\text{V}, 3 \text{mA}$

KOBE KOGYO

 CORPORATION

Manufactured exclusively for NUCLEAR SUPPLIES by the world's first company to successfully build transistorized nuclear equipment, Kobe Kogyo Corporation of Japan.

Complete line of instruments, modules, automatic systems and accessories listed in new catalog. Write to P.O. Box 312, Encino, California, or Telephone, Area Code 213, 787-1722.



News...

from



- 1. NCC acquires Volk Laboratories.** We have acquired the Chicago (Skokie) radiopharmaceutical laboratories of Volk Radiochemical Company — now we can offer direct shipment to Wisconsin, Minnesota, Iowa, Michigan, and other midwestern states. Plus several new products.
- 2. NCC opens New York Laboratory.** We have opened a new laboratory in New York City, offering fast direct delivery to the entire East Coast.

These expanded facilities mean additional laboratories and shipping points — more flexibility so we can meet your needs exactly. You'll get the product you want, delivered when you want it.


To place your order, call the nearest office, collect if you wish.



FOR BIOLOGICAL RESEARCH

... CLINICAL DIAGNOSIS

... SPECIAL STUDIES



ARMAC® provides rapid, accurate measurement of radioactivity

The Armac Scintillation Detector is ideally suited for *in vivo* research using gamma-emitting tracers. Its 4-pi counting geometry and high counting efficiency require the administration of only small doses of radioactivity, yet assure significant measurements in short counting periods.

- In biological research Armac can be used to count radioactivity in small laboratory animals in the study of retention and excretion of various gamma-emitting isotopes. Because only small doses of radioactivity are required, sacrificing is unnecessary and the animal acts as its own control for repeated experimentation.
- In clinical research Armac is useful in nutritional

studies and other applications where changes in the amount of circulating radioactivity in the blood can be directly related to other body functions. These changes can be measured accurately and rapidly using the blood-rich portion of the patient's forearm positioned in the Armac's counting chamber.

- In special studies Armac can be used to quickly and accurately determine the presence of low levels of radioactivity in bulk samples of meat, milk, water, blood, wastes, soil, or tissue with little or no sample preparation. The photograph at right, above, shows an Armac Detector counting radioactivity in a live fish in flowing stream water. (Courtesy Oak Ridge National Laboratory.)


A complete counting system includes the Model 446 Armac Scintillation Detector and any of several Packard Spectrometers, control units and ratemeters. There are more than ten Armac Detector Systems available to meet your specific counting requirements. Choose visual or recorded data output; digital or analog data presentation. Any system can be readily modified or expanded to meet changing research requirements. Special systems incorporating automatic sample changers are also available.

For detailed information contact your Packard Sales Engineer, or write for Bulletin 1024.



PACKARD INSTRUMENT COMPANY, INC.
2200 WARRENVILLE ROAD • DOWNERS GROVE, ILLINOIS 60515
TELEPHONE: 312/969-6000

ALL RIGHT, 3M, SO TELL ME ALL ABOUT
RADIATION SOURCES
THAT OUTLAST RADIUM AND Co-60—
OFFER GREATER SAFETY AND ECONOMY!



Mail to Nuclear Products, 3M Company,
Dept. BRB-55, St. Paul, Minn. 55119.

TELL ME MORE about 3M Brand Medical Radiation
Sources.

Name _____

Address _____

City, State, Zip _____

The coupon will bring you full information on cells and needles made with 3M Brand Cs-137 Radiating Microspheres. They offer isodose contours similar to radium and Co-60, with these advantages . . .

The 3M devices outlast their Co-60 counterparts by 6 times, retain 85% of initial activity for over 7 years and cost less in the long run. They cost only a fifth as much as radium devices. Outlast them, too, because there's no radon emission to impair service life.

And the 3M sources are *safer!* No internal build-up of toxic gases that can cause rupture during sterilization or


therapy. *Safer* because the actual isotope is confined in 3M Radiating Microspheres that render it chemically and physiologically inert. *Safer*, too, because the 3M cesium sources offer a softer, easier-to-confine gamma radiation than radium or Co-60, thus need less shielding.

Available in a variety of sizes and levels of radioactivity so you can use these devices with established dosage tables and applicators. Send us the coupon for full information.

Nuclear Products 

Large Crystal Radioisotope Scanners


OHIO-NUCLEAR
Model 54



FEATURES:

- Soundly engineered—research proven in use.
- 200 inches (500 cm) per minute scanning speeds.
- Choice of above or beneath table scanning.
- Dual Scanner capability—now or later.
- Gammagraphic Photoscan & Dot-tapper.
- Prompt service by capable factory engineers.
- 5-inch crystal scanner price—\$16,750.
- 8-inch crystal scanner price—\$20,700.

For Full Information WRITE OR CALL US • WE BELIEVE YOU WILL BE GLAD YOU DID.




OHIO-NUCLEAR, INC.

1725 FALL AVENUE

CLEVELAND, OHIO

216 - 621-8477



NUC-D-4-264

RENOGRAPHY PLUS

THE SYSTEM. Renaltron™ IV—integrated and compact—specifically designed for renography. Detects, analyzes, displays and records data more accurately and conveniently than you may have thought possible.

PLUS USAGE. Easily applied to cardiac output, hepatic function, and circulation studies.

FEATURING. Analyzer output for your tape recorder. Lets you preserve raw data for further analysis.

NEW EASE. Mobile stand, independently adjustable detector arms, rotating console on swing-away supporting arm—all this makes for facility of probe adjustment, patient observation, and electronics operation.

LOOKING AHEAD. Solid state circuitry assures reliable long-term performance. Construction provides for future mounting of third detector arm and probe.

MODEL 4130. Two 2-inch collimated scintillation detectors, mobile stand with two adjustable arms, movable shelf for control console, H. V. supply, dual ratemeter, dual single-channel analyzers, and dual-pen servo chart recorder. Write us for details.




NUCLEAR-CHICAGO

A DIVISION OF NUCLEAR-CHICAGO CORPORATION

333 Howard, Des Plaines, Ill. 60018 U.S.A.

In Europe: Donker Curtiusstraat 7
Amsterdam W, The Netherlands

Scientists and engineers interested in challenging career opportunities are invited to contact our personnel director.



organs
currently demonstrable by
Picker MagnaScanner techniques



 **MEDOTOPES®**
SQUIBB RADIOPHARMACEUTICALS

Medotopes reflect the latest developments in nuclear medicine. All provide the utmost in safety and convenience. All have unique packaging safeguards so that direct contact is never required. Exclusive lead shield enclosures are fitted with bottle caps that unscrew automatically. Saf-Tag® vials and bottles are carefully encased and double protected by transparent, shatterproof plastic coatings, and shipping cartons have convenient "pull-tab" openers. And, each preparation is custom-handled, each delivery custom-routed by Squibb Traffic Service. Access to three major airports expedites shipment.

Squibb Radiopharmaceuticals are available to the AEC-licensed physician. For full information, write to Professional Service Dept., Squibb, 745 Fifth Avenue, New York 22, N. Y.

SQUIBB 
Squibb Quality—the Priceless Ingredient
SQUIBB DIVISION **Olin**

HARSHAW NaI (TI) CRYSTALS

pace Health Physics developments



IN WHOLE BODY COUNTING

- Crystals to 11½" Diameter
- Measured background spectrum
- Resolution and Stability Performance Guaranteed

IN SCANNING

- Best energy resolution
- Small Detectors
- Multi-Crystal Mosaics for Dynamic Scanning
- Large Diameter Crystals for Pinhole Camera Techniques

IN GAMMA DOSE RATE DETECTORS

- New Development in Scintillator-Photoconductor Module for Dose Rate Measurement
- High Sensitivity
- Fast Response Time

IN THERMOLUMINESCENT DOSIMETRY

- Activated LiF Offering
- High Sensitivity
- Reproducibility
- Linear Energy Response

IN CHARGED PARTICLE DETECTION

- Surface Barrier Diodes offering minimum size for Probe Devices
- Li Drift Detectors for Beta Detection

Detailed information can be obtained from:



Crystal-Solid State Division
THE HARSHAW CHEMICAL COMPANY
1845 East 97th Street • Cleveland, Ohio 44106 • Telephone 216 721-8368
Utrecht, Netherlands—Contact Harshaw-Van Der Hoorn N. V.
Frankfurt, W. Germany—Contact Harshaw Chemie GmbH



NEW! A STANDARD ION CHAMBER For the Laboratory



Model 11-140

For standardization purposes and accurate measurements of Gamma and Beta sources

Also used for:

- Calibration of "Generator-Milking" Radioisotopes:
TECHNETIUM-99^m,
Iodine-132, Gallium-68,
Strontium-87^m
and others.
- Re-Checking shipments from suppliers, in original containers:

For complete details —
Send for 12-page brochure, NE-3



ATOMIC ACCESSORIES, INC.
811 West Merrick Road
Valley Stream, N.Y.



Subsidiary of Baird-Atomic, Inc..

CLINICAL CONSULTANT

**Radiation Protection Division
Department of National Health and Welfare**

Ottawa, Ontario

Salary up to \$15,100

A medical doctor is required immediately to serve as an adviser in all matters relating to the regulatory control of the use of radioisotopes in medicine.

The doctor must: be fully registered, or qualified to register, to practice medicine in a province of Canada; have five years of acceptable experience in medical practice or a combination of post-graduate training in this field, three of which shall have been in medical activities similar or closely related to the work to be performed.

**For further information and application forms,
please write IMMEDIATELY to the
CIVIL SERVICE COMMISSION OF CANADA, OTTAWA 4,
referring to Circular 65-597 (HW-HRS-9).**

ADVERTISING INDEX

Abbott Laboratories North Chicago, Illinois Inside Front Cover, i	Nuclear Consultants Corporation St. Louis, Missouri Inside Back Cover, v
Atomic Accessories, Inc Valley Stream, New York xii	Nuclear Supplies Encino, California iii
Civil Service Commission Ottawa, Canada xiii	Ohio-Nuclear Cleveland, Ohio viii
Harshaw Chemical Company Cleveland, Ohio xii	Packard Instrument Chicago, Illinois vi
3M Nuclear Products St. Paul, Minnesota vii	Picker-X-Ray Corporation White Plains, New York x
Nuclear-Chicago Des Plaines, Illinois Back Cover, ix	Squibb, E. R. and Sons New York, New York xi
	Tracerlab Waltham, Massachusetts xviii


A Valuable Addition To Your Professional Library

Journal of

NUCLEAR



*Official Publication
Society of Nuclear Medicine*



AN IMPORTANT NEW JOURNAL

featuring

Original articles in clinical medicine, basic and clinical medical research, physics and chemistry dealing with the use of isotopes in humans, and articles on related subjects. The latter includes dosimetry, instrumentation, protection, techniques, biologic effects contributing to the use or effects of isotopes in clinical medicine or the clinical effects of ionizing radiation.

George E. Thoma, M.D., St. Louis—*Editor*
G. O. Broun, Jr., M.D., St. Louis, Titus C. Evans, Ph.D., Iowa City,
Neil Wald, M.D., Pittsburgh, Eugene L. Saenger, M.D., Cincinnati—*Associate Editors*


The Journal of NUCLEAR MEDICINE		
333 North Michigan Avenue, Chicago, Illinois 60601		
\$20.00 per year, U.S. \$21.00 Foreign		
Name		
Address		
City.....	State.....	Zip Code.....
<i>Please remit by check or money order.</i>		

**Published
Monthly**



**undergraduate
simplicity**

**post-doctoral
performance**



Tracerlab's new multi-mission scaler

The new model 132 M manual scaler is a versatile, convenient instrument that even a near-novice can use with professional aplomb — and get results!

What's more, this instrument utilizes high-reliability printed circuit boards and transistor circuitry throughout and features performance characteristics you'd hardly expect at its most reasonable price:

High speed (500 kc/sec) with 800-nanosecond resolving time, 3.5-millivolt sensitivity. Voltage-and-charge-sensitive inputs with 0.5%-linearity discriminator and 3000:1 input sensitivity range.

Other features include 0-3000 vdc power supply with better than 0.003% regulation, 0-60 minute electromechanical timer with 0.1 second accuracy, infinite preset times and visual readout. All

this adds up to the very latest in scaler techniques — the 132 M — from the first family in nuclear instrumentation.




TRACERLAB

A Division of Laboratory For Electronics, Inc.
WALTHAM, MASSACHUSETTS 02154

Richmond, California • Houston, Texas • Malines, Belgium • Sales Offices in Principal Cities
Film Badge Service • Health Physics • Bioassays • Sources • Nuclear Instrumentation • Radiochemicals
Radioactive Waste Disposal • Radiation Monitoring Instrumentation • Isotope Applications

SINGLE SOURCE



**For the New Short
Half-Life or Low Energy
Isotopes**

More and more users are looking to Nuclear Consultants Corporation for the latest developments in radiopharmaceuticals. Whatever your need, whenever it arises, you can depend on NCC for demonstrated leadership and service.

Chlormerodrin-Hg¹⁹⁷

Chlormerodrin-Hg²⁰³

 Proven brain and renal scanning agents. Hg¹⁹⁷ available at lower cost through NCC's "Controlled Decay" program.

1-Mercuri-2-Hydroxypropane-Hg¹⁹⁷

 For improved spleen scanning. Much simpler and more convenient than heat treating tagged red cells. MHP-Hg²⁰³ also available for experimental use.

Iodine-I¹²⁵

 I-¹²⁵ is a longer-life isotope than I-¹³¹, and offers significantly lower radiation dosage to the patient. Available as tagged albumin and other iodinated compounds.

Technetium-99m

For Scanning Brain...Thyroid...Liver...Lung



TechneKow[®] generator

Supplies ample material as sodium pertechnetate for daily use. Economical for as few as seven or eight brain scans per week. New #1810 Shielded Dispenser available — simple semi-automated method of obtaining sterile solution.

"Instant Technetium"



Sodium pertechnetate in sterile, pyrogen-free solution permits immediate use without the necessity for "milking," sterilization, or calibration.

Technetium-Sulfide Colloid



In sterile, pyrogen-free form, reported as providing excellent liver scans. Also usable for lung scanning following the technique of Quinn.


6517



**NUCLEAR CONSULTANTS
CORPORATION**

Box 6172, Lambert Field • St. Louis, Mo. 63145
314 PErshing 9-8927

CHICAGO • CLEVELAND • HOUSTON • LOS ANGELES • NEW YORK • SAN FRANCISCO • WASHINGTON, D.C.



NUC:D-5-226

DOWNTIME: ZERO

CIRCUIT SLIP-UPS ARE OUT. The ultimate in reliability: operating circuits in Nuclear-Chicago's Pho/Dot® isotope scanner. Premium, ultra-stable, solid-state components. Heat-isolated critical circuits. Regulated power supplies.

RUGGED, TROUBLE-FREE CONSTRUCTION IS IN. In Pho/Dot, extra strength for rigidity and long-term operation. Scanning motions controlled by low-friction drives for reproducible results. Heavy-duty casters for mobility. Leveling jacks to simplify horizontal positioning.

SECOND-BEST IS OUT. Control of Pho/Dot quality every step of the way. Test after test of each subassembly. Exhaustive tests of the complete instrument. Precision and reliability built in—to stay.

DEPENDABLE SERVICE EVERYWHERE IS IN. The Pho/Dot guarantee: one-year, all-inclusive coverage of parts, labor, shipping, preventive maintenance, and travel time.

PHO/DOT FACTS. Many of your colleagues know. Your Nuclear-Chicago sales engineer can tell you even more about Pho/Dot's advanced design and foolproof operation. Or write us.



NUCLEAR-CHICAGO

A DIVISION OF NUCLEAR-CHICAGO CORPORATION
313 Howard, Des Plaines, Ill. 60018 U.S.A.

In Europe: Donker Curtiusstraat 7
Amsterdam W, The Netherlands

Scientists and engineers interested in challenging career opportunities are invited to contact our personnel director.