¹¹C-Dihydrotetrabenazine β-Cell Imaging

he prevalence of diabetes is dramatically high, and the acquired type II form is particularly expected to further increase considerably (from 6.6% of the world population [285 million people] in 2010 to 7.8% [439 million people] in 2030) (1). Diabetes is caused by a disrupted secretion of insulin and a subsequent increase in plasma glucose levels. Insulin is produced by the β -cells of the 1-2 million islets of Langerhans in a pancreas. These islets represent only about 2%-3% of the pancreas cell mass, which consists mostly of exocrine tissue. Autopsy studies have shown that the relative β -cell volume is reduced by more than 50% in diabetic subjects,

See page 1439

and animal models in which more than 50% of the β -cell volume had been destroyed have shown impaired glucose homeostasis (2).

Thus, an important tool for diabetes research would be noninvasive in vivo quantification of β -cell mass, particularly longitudinal quantification for monitoring new antidiabetic therapies. Research on imaging of β -cell mass has focused on various techniques (PET, SPECT, MRI, optical imaging) using, among others, peptides, antibodies, manganese, annexin derivatives, sulfonylureas, and 6-fluoro-L-dopa (*3*,*4*). The challenge of imaging β -cell mass is the small dimensions (~0.1mm diameter) of the islets of Langerhans, combined with the relatively small fraction of β -cells (~2%) in the pancreas, requiring that tracers or contrast agents have a high selectivity for β -cells versus exocrine tissue.

A study that mapped gene expression showed that the vesicular monoamine transporter 2 (VMAT2) is expressed by β -cells but not by exocrine pancreatic cells (5). ¹¹C-dihydrotetrabenazine, a PET tracer that had been developed to visualize VMAT2 in the brain (6), was applied in PET studies of diabetes and showed less pancreatic uptake in diabetic rats than in control rats (5) and a longitudinal decrease in the pancreatic standardized uptake value ($\leq 65\%$) of spontaneous diabetic rats (7). Encouraged by these results, researchers have synthesized and are evaluating several ¹⁸F derivatives of dihydrotetrabenazine (8.9).

The preclinical results were partially confirmed in clinical studies that showed lower ¹¹C-dihydrotetrabenazine binding potentials in type 1 diabetic patients than in control subjects. Patients with longstanding type 1 diabetes are, however, expected to have a near-total loss of β-cell mass, whereas the observed ¹¹C-dihydrotetrabenazine binding potential was still 86% of the value in control subjects (10). It was hypothesized that the discrepancy between preclinical and clinical studies was caused by modeling artifacts introduced by using the kidney as a reference tissue, by the use of VMAT2 expression in other cell types of the pancreas, and by species differences in tracer kinetics or radiometabolite accumulation.

Using ex vivo autoradiography and in vitro tissue-binding competition studies, Fagerholm et al. now show in this issue of *The Journal of Nuclear Medicine* that non-VMAT2 binding of ¹¹C-dihydrotetrabenazine to exocrine pancreas exceeds VMAT2 binding to β -cells (*11*). Interestingly, the non-VMAT2 binding (and thus the total pancreas binding) of ¹¹C-dihydrotetrabenazine is also decreased in streptozotocin diabetic rats. In vitro autoradiography of the human pancreas confirmed that dihydrotetrabenazine binding represents mainly nonspecific binding.

Similar results have been published recently for 9-¹⁸F-fluoroethyl-(+)-dihydrotetrabenazine by Eriksson et al. (12), who observed specific binding of the tracer to isolated islet tissue but found that specific islet binding did not sufficiently exceed nonspecific binding in the abundant exocrine pancreatic tissue to allow β -cell mass quantification.

In the evaluation of future PET tracers of β -cell mass (VMAT2), it will be important to include a thorough highresolution ex vivo and in vitro evaluation of β -cell specificity focusing not only on high specific binding to β -cells but also on low binding of the tracer or its metabolites to the abundant surrounding exocrine pancreatic tissue.

Guy M. Bormans

Laboratory for Radiopharmaceutical Sciences K.U. Leuven Leuven, Belgium

REFERENCES

- International Diabetes Federation Web site. Diabetes and Impaired Glucose Tolerance: Global Burden—Prevalence and Projections, 2010 and 2030. Available at: http://www.diabetesatlas.org/ content/diabetes-and-impaired-glucose-tolerance. Accessed June 30, 2010.
- Matveyenko AV, Butler PC. Relationship between beta-cell mass and diabetes onset. *Diabetes Obes Metab.* November 2008;(10 suppl 4):23–31.
- Moore A. Advances in beta-cell imaging. Eur J Radiol. 2009;70:254–257.
- Lin M, Lubag A, McGuire MJ, et al. Advances in molecular imaging of pancreatic beta cells. *Front Biosci.* 2008;13:4558–4575.
- Simpson NR, Souza F, Witkowski P, et al. Visualizing pancreatic beta-cell mass with [¹¹C] DTBZ. Nucl Med Biol. 2006;33:855–864.
- Koeppe RA, Frey KA, Vander Borght TM, et al. Kinetic evaluation of [¹¹C]dihydrotetrabenazine by dynamic PET: measurement of vesicular monoamine transporter. J Cereb Blood Flow Metab. 1996;16: 1288–1299.

Received May 25, 2010; revision accepted Jun. 3, 2010.

For correspondence or reprints contact: Guy M. Bormans, Laboratory for Radiopharmaceutical Sciences, K.U. Leuven, Herestraat 49, Leuven B-3000, Belgium.

E-mail: guy.bormans@pharm.kuleuven.be COPYRIGHT © 2010 by the Society of Nuclear

Medicine, Inc. DOI: 10.2967/jnumed.110.077602

- Souza F, Simpson N, Raffo A, et al. Longitudinal noninvasive PET-based beta cell mass estimates in a spontaneous diabetes rat model. *J Clin Invest.* 2006;116:1506–1513.
- Kung HF, Lieberman BP, Zhuang ZP, et al. In vivo imaging of vesicular monoamine transporter 2 in pancreas using an ¹⁸F epoxide derivative of tetrabenazine. *Nucl Med Biol.* 2008;35:825–837.
- Kung MP, Hou C, Lieberman BP, et al. In vivo imaging of beta-cell mass in rats using ¹⁸F-FP-(+)-DTBZ: a potential PET ligand for studying diabetes mellitus. J Nucl Med. 2008;49:1171–1176.
- Goland R, Freeby M, Parsey R, et al. ¹¹Cdihydrotetrabenazine PET of the pancreas in subjects with long-standing type 1 diabetes and in healthy controls. *J Nucl Med.* 2009;50:382–389.
- Fagerholm V, Mikkola KK, Ishizu T, et al. Assessment of islet specificity of dihydrotetrabenazine radiotracer binding in rat and human pancreas. *J Nucl Med.* 2010;51: 1439–1446.
- Eriksson O, Jahan M, Johnström P, et al. In vivo and in vitro characterization of [¹⁸F]-FE-(+)-DTBZ as a tracer for beta-cell mass. *Nucl Med Biol.* 2010;37:357–363.