

SUPPLEMENTAL TABLE 1: Overview of PET/MR imaging systems described in literature

Short name	Company/university	Year	Usage*	MR field strength (T)	Scintillator†	Detector technology‡	Crystal rings	Blocks/ring	Crystal block	Crystal size (mm ³)	Axial FOV (cm)	Transaxial FOV (cm)	Spatial resolution (mm)	Sensitivity (%)	Timing resolution (ns)	Energy resolution (%)	Reference
1 "McPET I"	University of California Los Angeles, CA (US)	1997	P-R-Sim	0.2	LSO	MC-PMT/fibers	1	48	16§	2x2x10	1.0	<3.8	2.1	N.A.	20	41%	(1)
2 "McPET II"	University of California Los Angeles, CA (US)	1997	P-R-Sim	0.2 and 9.4	LSO	MC-PMT/fibers	1	72	24§	2x2x5	0.2	<5.4	2.1	N.A.	26	45%	(2)
3 Kings	Kings College London (UK)	2005	P-R-Sim	3	LSO	PSPMT/fibers	1	8	1x4x4	2x3x5	0.3	5.6	1.4-3.2	0.95	20.9	70%	(3)
4 Split Magnet	University of Cambridge Cambridge (UK)	2006	P-R-Sim	1	LSO	PSPMT/fibers	48	24	12x12	1.52x1.52x10	7.2	<14.7	1.6	3.6	3.6	27%	(4)
5 West Virginia	West Virginia University Morgantown, WV (US)	2007	P-R-Sim	3	LSO	PSPMT/fibers	20	2	20x20	2.5x2.5x15	5.0	8.0	2.65	N.A.	N.A.	85%	(5)
6 Kobe City I	Kobe City College of Technology, Kobe (JP)	2009	P-R-Sim	0.15	MLS	PSPMT/fibers	3	32	2x2x2	2.5x3.5x3.5	0.5	N.A.	2.8-3.0	0.25	6	29%/50%	(6)
7 Fieldcycled	Western Ontario London, ON (CA)	2009	P-R-Seq	0.30	BGO	PSPMT/fibers	8	2	8x8	6.2x5.6x30	-5.0	N.A.	N.A.	N.A.	N.A.	N.A.	(7)
8 Kobe City II	Kobe City College of Technology, Kobe (JP)	2010	P-R-Sim	0.30	LGSO	PSPMT/fibers	11	16	11x9x2	1.9x2.2x6/7	2.1	8.0	2.4-2.9	1.5	5	40%	(8)
9 Nagoya	Nagoya University Nagoya (JP) Mediso Ltd. Budapest (HU)	2012	P-R-Sim	0.30	LGSO	PSPMT/fibers	13	16	11x13x2	0.9x1.3x5/6	2.1	5.6	1.2	1.5	3.1	24%	(9)
10 "nanoScan"		2013	P-Co-Seq	1	LYSO	PSPMT	81	12	39x81	1.12x1.12x13	9.4	4.5/9.4/12.0	1.5-2.1	8.4	N.A.	N.A.	(10)
11 UCD	University of California Davis, CA (US)	2006	P-R-Sim	7	LSO	PSAPD/fibers	8	16	8x8	1.43x1.43x6	1.2	3.5	1.1-1.7	0.35	2.5-7.5	25±5.5%	(11)
12 Tübingen I	University of Tübingen Tübingen (DE)	2007	P-R-Sim	7	LSO	APD	12	10	12x12	1.6x1.6x4.5	1.9	4.0	1.9	0.23	8	14.6%-29.9%	(12)
13 "ratCap"	Brookhaven National Laboratory, Upton, NY (US)	2011	P-R-Sim	9.4	LSO	APD	8	12	4x8	2.2x2.2x5	1.8	3.8	1.2	0.3	13.9	13%	(13)
14 Tübingen II	University of Tübingen Tübingen (DE)	2013	P-R-Sim	7	LSO	APD	45	16	15x15	1.5x1.5x10	7.2	7.2	>1.5	N.A.	N.A.	N.A.	(14)
15 Sogang I	Sogang University Seoul (KR)	2011	P-R-Sim	3	LYSO	SiPM	4	16	4x4	3x3x10	1.3	<7.0	2.8-4.1	0.65	1.5	17.6%	(15)
16 Seoul	Seoul National University, Seoul (KR)	2012	P-R-Sim	3	LGSO	SiPM	20	12	20x18	1.5x1.5x7	3.2	13.6	1.0-1.5	0.12	1.23	14%	(16)
17 RWTH	RWTH Aachen Aachen (DE)	2012	P-R-Sim	3	LYSO	SiPM/digital SiPM fibers/	22	10	22x22	1.3x1.3x10	3.0	16.0	2.2	0.6	N.A.	12.9%	(17)
18 Eulji	Eulji University Gyeonggi (KR)	2012	P-R-Sim	3	LYSO	SiPM	6	12	6x6	2.47x2.74x20	N.A.	N.A.	N.A.	N.A.	2	18%	(18)
19 Kobe City III	Kobe City College of Technology, Kobe (JP)	2012	P-R-Sim	0.15	LGSO	SiPM phoswich	11	16	11x9x2	1.1x1.2x5/6	13.2	8.0	1.6	0.6	3.1	27%	(19)
20 Sogang II	Sogang University Seoul (KR)	2013	P-R-Sim	3	LYSO	SiPM	4	72	4x4	3x3x20	12.9	25.0	3.1-6.6	0.33	4.2	18.1%	(20)
22 "Ingenuity"	Koninklijke Philips NV Eindhoven (NL)	2011	C-Co-Seq	3	LYSO	PMT	44	28	23x44	4x4x22	18	60	4.6-5.3	0.7	0.52	11.6%	(21)
21 "brainPET"	Siemens AG München (DE)	2008	C-R-Sim	3 and 9.4	LSO	APD	72	32	12x12	2.5x2.5x20	19.3	32	1.8-6.4	7.2	4.9	17.1%	(22,23)
23 "mMR"	Siemens AG München (DE)	2012	C-Co-Sim	3	LSO	APD	64	56	8x8	4x4x20	25.8	59.4	4.3-6.6	1.5	2.93	14.5%	(24)
24 GE¶	GE Healthcare Waukesha, WI (US)	2014	C-Co-Sim	3	LYSO	SiPM	45	112	4x9	3.95x5.3x25	25	60	3.9	2.4	0.40	10.5%	(25)

*Code specifying the field (P: preclinical, C: clinical), the status (R: research, Co: commercial) and operation (Sim: simultaneous, Seq: sequential).

†Abbreviations: BGO: bismuth germanate; LGSO: lutetium gadolinium oxyorthosilicate; LSO: lutetium oxyorthosilicate; LYSO: lutetium yttrium oxyorthosilicate; MLS: mixed lutetium silicates.

‡Abbreviations: APD: avalanche photodiode; MC: multichannel; PMT: photomultiplier tube; PS: position sensitive; SiPM: silicon photomultiplier.

§Describes the number of crystals per MC-PMT.

||A variant of the ratCap has recently been introduced commercially by MR Solutions Ltd.

¶This scanner has been presented at the RSNA and SNM conferences but is not yet commercially available.

REFERENCES

1. Shao Y, Cherry SR, Farahani K, et al. Simultaneous PET and MR imaging. *Phys Med Biol.* 1997;42:1965–1970.
2. Shao Y, Cherry SR, Farahani K, et al. Development of a PET detector system compatible with MRI/NMR systems. *IEEE Trans Nucl Sci.* 1997;44:1167–1171.
3. Mackewn JE, Halsted P, Charles-Edwards G, et al. Performance evaluation of an MRI-compatible pre-clinical PET system using long optical fibers. *IEEE Trans Nucl Sci.* 2010;57:1052–1062.
4. Hawkes RC, Fryer TD, Siegel S, Ansorge RE, Carpenter TA. Preliminary evaluation of a combined microPET-MR system. *Technol Cancer Res Treat.* 2010;9:53–60.
5. Raylman RR, Majewski S, Velan SS, et al. Simultaneous acquisition of magnetic resonance spectroscopy (MRS) data and positron emission tomography (PET) images with a prototype MR-compatible, small animal PET imager. *J Magn Reson.* 2007;186:305–310.
6. Yamamoto S, Hatazawa J, Imaizumi M, et al. A multi-slice dual layer MR-compatible animal PET system. *IEEE Trans Nucl Sci.* 2009;56:2706–2713.
7. Bindseil GA, Gilbert KM, Scholl TJ, Handler WB, Chronik BA. First image from a combined positron emission tomography and field-cycled MRI system. *Magn Reson Med.* 2011;66:301–305.
8. Yamamoto S, Imaizumi M, Kanai Y, et al. Design and performance from an integrated PET/MRI system for small animals. *Ann Nucl Med.* 2010;24:89–98.
9. Yamamoto S, Watabe H, Kanai Y, et al. Development of a flexible optical fiber based high resolution integrated PET/MRI system. *Med Phys.* 2012;39:6660–6671.
10. Nagy K, Tóth M, Major P, et al. Performance evaluation of the small-animal nanoScan PET/MRI system. *J Nucl Med.* 2013;54:1825–1832.
11. Wu Y, Catana C, Farrell R, et al. PET performance evaluation of an MR-compatible PET insert. *IEEE Trans Nucl Sci.* 2009;56:574–580.
12. Judenhofer MS, Wehr HF, Newport DF, et al. Simultaneous PET-MRI: a new approach for functional and morphological imaging. *Nat Med.* 2008;14:459–465.
13. Maramraju SH, Smith SD, Junnarkar SS, et al. Small animal simultaneous PET/MRI: initial experiences in a 9.4 T microMRI. *Phys Med Biol.* 2011;56:2459–2480.
14. Wehr HF, Hossain M, Lankes K, et al. Simultaneous PET-MRI reveals brain function in activated and resting state on metabolic, hemodynamic and multiple temporal scales. *Nat Med.* 2013;19:1184–1189.
15. Kang J, Choi Y, Hong KJ, et al. A small animal PET based on GAPDs and charge signal transmission approach for hybrid PET-MR imaging. *J Instrum.* 2011;6:P08012–P08012.
16. Yoon HS, Ko GB, Kwon SI, et al. Initial results of simultaneous PET/MRI experiments with an MRI-compatible silicon photomultiplier PET scanner. *J Nucl Med.* 2012;53:608–614.
17. Wehner J, Weissler B, Dueppenbecker P, et al. PET/MRI insert using digital SiPMs: investigation of MR-compatibility. *Nucl Instrum Meth A.* 2014;734:116–121.
18. Hong SJ, Kang HG, Ko GB, Song IC, Rhee J-T, Lee JS. SiPM-PET with a short optical fiber bundle for simultaneous PET-MR imaging. *Phys Med Biol.* 2012;57:3869–3883.
19. Yamamoto S, Watabe T, Watabe H, et al. Simultaneous imaging using Si-PM-based PET and MRI for development of an integrated PET/MRI system. *Phys Med Biol.* 2012;57:N1–N13.
20. Hong KJ, Choi Y, Jung JH, et al. A prototype MR insertable brain PET using tileable GAPD arrays. *Med Phys.* 2013;40:042503.
21. Zaidi H, Ojha N, Morich M, et al. Design and performance evaluation of a whole-body Ingenuity TF PET-MRI system. *Phys Med Biol.* 2011;56:3091–3106.
22. Kolb A, Wehr HF, Hofmann M, et al. Technical performance evaluation of a human brain PET/MRI system. *Eur Radiol.* 2012;22:1776–1788.
23. Shah NJ, Oros-Peusquens A-M, Arrubla J, et al. Advances in multimodal neuroimaging: hybrid MR-PET and MR-PET-EEG at 3 T and 9.4 T. *J Magn Reson.* 2013;229:101–115.
24. Delso G, Fürst S, Jakoby B, et al. Performance measurements of the Siemens mMR integrated whole-body PET/MR scanner. *J Nucl Med.* 2011;52:1914–1922.
25. Levin C, Glover G, Deller T, McDaniel D, Peterson W, Maramraju SH. Prototype time-of-flight PET ring integrated with a 3T MRI system for simultaneous whole-body PET/MR imaging [abstract]. *J Nucl Med.* 2013;54(suppl 2):45P.